thermoMACS™ Separation Unit

User manual

Order no. 130-091-136

130-091-137/02

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

Miltenyi Biotec Inc.
12740 Earhart Avenue
Auburn CA 95602, USA
Phone 800 FOR MACS,
+1 530 888 8871
Fax +1 530 888 8925
macs@miltenyibiotec.com

www.miltenyibiotec.com

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use.

Table of Contents

▲ Warnings and precautions.....	2
▲ Warnhinweise und Vorsichtsmaßnahmen.....	3
▲ Advertisements et précautions.....	4
▲ Avvertimenti e precauzioni.....	5
▲ Advertencias y precauciones.....	6
▲ 1. Introduction.....	7
1.1 Description of the thermoMACS™ Separator.....	7
1.2 Applications.....	7
1.3 The MACS® Technology.....	8
▲ 2. Technical service.....	9
▲ 3. The thermoMACS™ Separation Unit.....	10
3.1 Additional materials required.....	10
3.2 Equipment supplied with the thermoMACS™ Separation Unit.....	10
3.3 Key components of the thermoMACS™ Separation Unit.....	11
▲ 4. Installation.....	12
▲ 5. Maintenance.....	12
▲ 6. Instructions for use.....	13
▲ 7. Status of the thermoMACS™ Separator.....	13
▲ 8. Isolation and processing of molecules with the thermoMACS™ Separator.....	14
8.1 Tips & hints.....	14
8.2 Special protocols.....	15
▲ 9. thermoMACS™ accessories.....	16
▲ 10. Related products.....	16
▲ 11. Technical data and specifications.....	17

Warnings and precautions

- ▲ The thermoMACS™ Separator is for research applications only. Read and follow all operating instructions carefully before proceeding with installation and operation of the system.
- ▲ The thermoMACS Separator is equipped with an extremely powerful magnet. Its magnetic field can damage computer screens, watches, floppy discs, and other objects that react to magnetic fields. Never allow the thermoMACS Separator to be closer than 20 cm to any magnetically sensitive object.
- ▲ Electric devices pose the risk of an electric shock. To reduce the risk of an electric shock, take special care while handling fluids. Caution must be taken that fluids do not enter the interior of the device. Clean up spills immediately. Unplug the power cord before cleaning the thermoMACS Separator.
- ▲ Use the thermoMACS Separator only with the enclosed power supply. Use of a non-specified power supply may damage the device.

Limited warranty

- If damage is due to the misuse of the thermoMACS Separator, all warranties shall be void. This includes damage to the device itself as well as any consequential damages.
- The appropriate use of the device is specified in the thermoMACS user manual.
- Warranties are void if the damage is due to the use of separation reagents which are not produced and distributed by Miltenyi Biotec. This includes damage to the device itself as well as any consequential damages.
- Under all circumstances, liability is limited to the value of the thermoMACS Separator. Any liability beyond the value is excluded.
- The thermoMACS Separation Unit is warranted for 12 months from the date of delivery.

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1. Introduction

1.1 Description of the thermoMACS™ Separator

The thermoMACS™ Separator is a heatable magnetic separator specifically developed for molecular biology applications. The thermoMACS Separator, in combination with μ Columns and μ MACS™ MicroBeads, is used for the isolation of highly pure biomolecules such as mRNA, sequence-specific DNA or proteins. In addition to applications at room temperature, the thermoMACS Separator allows molecular biological reactions to be carried out at increased temperatures of 37 °C or 42 °C such as enzymatic reactions. The elution volume from a μ Column may be as low as 50–150 μ L and therefore allows a highly concentrated target fraction to be collected, which is convenient for further downstream applications. The thermoMACS Separator can hold up to four μ Columns to process four samples simultaneously.

1.2 Applications

Solid phase enzymatic reactions on μ Columns are an easy and elegant method for purification of biomolecules. While the molecule is immobilized in the column, it is incubated with an enzymatic solution. To ensure optimal reaction conditions, the column can be heated to 37 °C and 42 °C, respectively, which covers the reaction conditions for most of the enzymes used in molecular biology applications. After the reaction, buffer and reaction components can be easily and efficiently washed out, thereby allowing serial reactions to be carried out in a single μ Column. Expensive and time-consuming purifications between single enzyme reactions is not necessary. Consequently, there is no loss of biomolecules due to the transfer from one tube to the next. This can significantly improve sensitivity of downstream assays like (real-time) PCR or microarray analysis after mRNA isolation and in-column cDNA synthesis or labeling, respectively.

Furthermore, enzymes may be immobilized in the column. This allows the enzymatic modification of substrates and subsequent elution of these while the enzyme remains on the column.

Performing an enzymatic reaction using biomolecules in a column matrix – such as the MACS™ Columns – is a new application and is filed for a patent.

7

1. Introduction

1.3 The MACS® Technology

For magnetic purification, biomolecules are labeled with target specific μ MACS™ MicroBeads. Subsequently they are separated in a μ Column that is placed in a strong magnetic field. μ MACS MicroBeads are optimized for the isolation of biomolecules, such as RNA or proteins. They are conjugated to specific molecules e.g., proteins like Protein A, Protein G, epitope tag specific antibodies, Streptavidin or single-stranded oligonucleotides, and can be used for magnetic labeling of virtually any target molecule. The colloidal behavior of the MicroBeads facilitates and accelerates the binding of μ MACS MicroBeads to the target molecules. The MicroBead labeled sample is then applied to the column which is placed in the magnet. While the labeled molecules are retained in the magnetic matrix, the non-target molecules are effectively rinsed out of the column. As the columns flow is by gravity, tedious centrifugation steps are eliminated from the procedure.

For elution of the target molecule an appropriate elution buffer which disrupts the binding between the target molecule and its immobilized ligand, is applied onto the column. While the target molecule is collected in the eluate, the probe remains within the column. A second option is the release of the target molecule–probe complex by removing the column from the magnetic field. In this case, the complete complex can be eluted using a buffer of choice.

8

2. Technical service

Germany/Austria/Switzerland

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach, Germany
Phone: +49 2204 8306-0
Fax: +49 2204 85197
E-mail: macs@miltenyibiotec.de

USA/Canada

Miltenyi Biotec Inc.
12740 Earhart Avenue
Auburn, CA 95602, USA
Phone: +1 530 888 8871
Fax: +1 530 745 2899
Freephone: 800 FOR MACS
Freephone: +1 800 367 6227
E-mail: macs@miltenyibiotec.com

Australia

Miltenyi Biotec Australia Pty. Ltd.
Phone: +61 02 8877 7400
Fax: +61 02 9889 5044
E-mail: macs@miltenyibiotec.com.au

Benelux

Miltenyi Biotec B.V.
Phone: 0800 94016
Fax: 0800 99626
E-mail: macs@miltenyibiotec.nl

China

Miltenyi Biotec Shanghai Office
Phone: +86 21 62351005
Fax: +86 21 62350953
E-mail: macs@miltenyibiotec.com.cn

France

Miltenyi Biotec
Phone: +33 1 56 98 16 16
Fax: +33 1 56 98 16 17
E-mail: macs@miltenyibiotec.fr

Italy

Miltenyi Biotec S.r.l.
Phone: +39 51 64 60 411
Fax: +39 51 64 60 499
E-mail: macs@miltenyibiotec.it

Singapore

Miltenyi Biotec Asia Pacific Pte Ltd.
Phone: +65 6238 8183
Fax: +65 6238 0302
E-mail: macs@miltenyibiotec.com.sg

Spain

Miltenyi Biotec S.L.
Phone: +34 91 512 1290
Fax: +39 91 512 1291
E-mail: macs@miltenyibiotec.es

UK

Miltenyi Biotec Ltd.
Phone: +44 1483 799 800
Fax: +44 1483 799 811
E-mail: macs@miltenyibiotec.co.uk

For countries not listed, please contact
Miltenyi Biotec GmbH, Germany.

9

3. The thermoMACS™ Separation Unit

3.1 Equipment supplied with the thermoMACS™ Separation Unit

thermoMACS™ Separator

μ MACS™ Sealing Solution, 1 mL (additional Sealing Solution can be purchased separately # 130-091-160)

Power supply

Main power supply cable

Low voltage power supply cable with temperature sensor

User manual

3.2 Additional materials required

MACS® MultiStand (# 130-042-303)

The MACS™ MultiStand is necessary for setting up μ MACS™, MiniMACS™, QuadroMACS™, OctoMACS™, and thermoMACS™ Separation Units. MACS Separation Units are magnetically attached to the metal stand. The MultiStand can hold one thermoMACS Separator.

μ Columns (# 130-042-701)

μ Columns are specifically designed for small scale molecular biology applications in combination with μ MACS MicroBeads. The magnetically labeled molecules are efficiently retained in the column matrix while unlabeled material is washed away. Elution of the target molecules can be performed in a convenient, small volume, e.g., 50–100 μ L. μ Columns can be used for a large variety of applications such as RNA isolation, (co-)immunopurification, tagged protein isolation, etc.

Reagents

A complete list of MACS™ molecular products and protocols is available at www.miltenyibiotec.com

Rack

Please use an appropriate rack for the collection tubes.

10


3. The thermoMACS™ Separation Unit

3.3 Key components of the thermoMACS™ Separation Unit

The key components of the thermoMACS™ Separation Unit are shown in figure 1.

1. The thermoMACS Separator generates a high gradient magnetic field to bind magnetically labeled biomolecules in the μ Column. The magnetic field is identical to that in the μ MACS™ Separator.
2. Heating elements are present in each column slot, heating up the reaction chamber in the μ Column to pre-defined temperatures of 37 °C or 42 °C.
3. A microprocessor-based electronic board and thermocouples in the heating elements control the temperature of the reaction chamber in the μ Column. This maintains a constant temperature throughout the whole separation process, which ensures stable and highly reproducible experimental results.
4. The temperature sensor integrated in the power supply cable of the thermoMACS Separator measures the room temperature. This allows precise amendment of the temperature in the column relative to variations in the local environment.

The thermoMACS™ Separation Unit (Figure 1)

11

4. Installation

1. Place the MACS™ MultiStand on a clean and stable laboratory bench.
 2. Fix the thermoMACS™ Separator to the MultiStand by holding the rear surface against the stand. It will remain attached to the MultiStand by internal magnets on the rear surface of the thermoMACS Separator.
 3. Take the low voltage power supply cable and insert the output into the plug on the left side of the thermoMACS Separator. The power supply cable is correctly orientated, when the red dot on the plug meets the red dot on the connector.
- ▲ **Note:** Use the thermoMACS Separator only with the enclosed power supply. Using a non-specified power supply may damage the device.
4. Connect one end of the main power supply cable with the power supply and plug it into an appropriate electrical outlet.

The thermoMACS Separator is now ready for separation. The column temperature is turned to heat off (yellow LED) as the default setting.

▲ **Note:** The thermoMACS Separator and the temperature sensor should not be used in direct proximity to any heat source or strong illuminants.

5. Maintenance

Clean up spills immediately. Unplug the power cable before cleaning the thermoMACS Separator. Use a detergent containing wash solution (e.g. 1% SDS solution) and/or alcohol (e.g. 70% ethanol solution) for disinfection and cleaning.


▲ **Note:** The thermoMACS Separation Unit should not be submerged under water, or cleaned in a dish washer, ultrasonic cleaner, autoclave, or similar appliances.

12

6. Instructions for use

The keypad (fig. 2) is used to select the temperature of the reaction chamber in the μ Column. Gently touching the control button on the keypad changes the temperature as indicated by the three light emitting diodes (LED): the red LED for 42 °C, the green LED for 37 °C, and the yellow LED for heat off, i.e. for working at room temperature. When the target temperature is changed via the control button to 37 °C or 42 °C, the corresponding temperature LED blinks until the target temperature is reached. Then the LED is permanently illuminated.

The keypad (Figure 2)

7. Status of the thermoMACS™ Separator

LED signals	Status of the thermoMACS™ Separator
No LED is on	The thermoMACS Separator may not be connected to the power supply.
The yellow LED is on	The thermoMACS Separator is in standby mode and the heating elements are deactivated.
The yellow LED blinks fast	The thermoMACS Separator has a hardware error. Please contact Technical Service.
The green LED blinks slowly	The 37 °C heating mode is on but the target temperature has not been reached.
The green LED is on	The target temperature of 37 °C has been reached.
The red LED blinks slowly	The 42 °C heating mode is on but the target temperature has not been reached.
The red LED is on	The temperature of 42 °C has been reached.

13

8. Isolation and processing of molecules

8.1 Tips & hints

▲ General considerations

Before starting an application on the thermoMACS™ Separator, consider whether the assay requires immediate heating, or heating after biomolecule separation:

Immediate heating:

Pre-heat the column in the thermoMACS Separator for 5–10 min before applying the buffer for equilibration. Use of pre-heated buffer is recommended. Apply the reaction solution directly after flow-through of the equilibration buffer. Immediately add 1 μ L of μ MACS™ Sealing Solution onto the μ Column.

Heating after biomolecule separation:

First, carry out the separation according to the standard μ MACS protocol at room temperature (thermoMACS Separator „off“). Before starting the heated reaction, set the thermoMACS Separator to the required temperature. If not otherwise specified in the μ MACS reagent user manual, wait for 5–10 minutes, then add the reaction solution. Use of pre-heated buffer for the reaction is recommended. Immediately add 1 μ L of μ MACS Sealing Solution.

▲ Warm-up time

The thermoMACS Separator requires 5–10 minutes to reach 37 °C or 42 °C. Permanent illumination of the LED indicates that the target temperature is reached. If the μ Column is inserted in the thermoMACS Separator before it is switched on, the column matrix will reach 37 °C or 42 °C after 5–10 minutes. If the μ Column is inserted in a pre-heated thermoMACS Separator, allow an additional five minute adjustment period to ensure that the matrix in the column is heated to the correct temperature.

▲ Buffer temperature

The use of cold buffer on a pre-warmed column is not recommended due to the possibility of air bubble formation in the column matrix.

▲ μ MACS Sealing Solution

The μ MACS Sealing Solution is a unique solution that creates a closed reaction chamber in the μ Column, thereby preventing evaporation of the

14

8. Isolation and processing of molecules

reaction mixture at reaction temperatures above room temperature. After the reaction solution is added on the heated μ Column, 1 μ L of the μ MACS Sealing Solution should be applied on top of the column matrix. Make sure to pipette the μ MACS Sealing Solution directly onto the iron sphere matrix in the μ Column so that it can spread and completely cover the surface of the matrix. If the μ MACS Sealing Solution is applied only on the side wall of the column reservoir it will not run down to the matrix due to its small volume. The volume of the applied μ MACS Sealing Solution should not exceed 2 μ L.

Generally, the μ MACS Sealing Solution only needs to be applied once. It stays on top of the column matrix and does not enter the matrix, because its density is lower than the buffer densities. When buffer is applied on top of a sealed column, the μ MACS Sealing Solution lets the buffer pass into the matrix and reseals the surface. However, the function of the μ MACS Sealing Solution can be impaired when working with detergents or other amphipathic substances (see paragraph „If amphipathic substances...“ below).

The μ MACS Sealing Solution does not interfere with the immobilized biomolecules in the column reaction chamber nor does it contaminate the elution fraction as long as working with detergents or other amphipathic substances is avoided (see paragraph „If amphipathic substances...“ below).

If amphipathic substances including fatty acids, phospholipids, and detergents, especially ionic detergents like SDS, are used after sealing the column with the μ MACS Sealing Solution, the effectiveness of the sealing may be compromised. For additional incubation steps we recommend applying a further 1 μ L of the μ MACS Sealing Solution.

If the current application is performed at room temperature (thermoMACS Separator switched to „heat off“), no μ MACS Sealing Solution is necessary.

8.2 Applications & protocols

For specific applications please refer to www.miltenyibiotec.com.

15

9. thermoMACS™ accessories

MACS™ Product	Quantity	Order no.
MACS MultiStand	1 piece	130-042-303
μ Column	20 columns	130-042-701
μ MACS Sealing Solution	1 mL	130-091-160

10. Related products

Please refer to www.miltenyibiotec.com for a complete list of MACSmolecular products.

MACS Product	Application	No. of tests per kit	Order no.
mRNA isolation and cDNA synthesis			
μ MACS mRNA Isolation Kit - Small Scale	mRNA isolation from cells, tissue and blood	20 10	130-075-201 130-090-276
μ MACS mRNA Isolation Kit - Large Scale	mRNA isolation from cells, tissue and blood	8 4	130-075-101 130-090-277
μ MACS mRNA Isolation Kit - For Total RNA	mRNA isolation from total RNA	8	130-075-102
μ MACS One-step cDNA Kit	mRNA isolation and in-column cDNA synthesis	20	130-091-902
μ MACS One-step cDNA Labeling Kit	mRNA isolation and in-column cDNA labeling	20	130-092-443
μ MACS One-step T7 Template Kit	mRNA isolation and in-column amplification	20	130-092-866
Protein and molecular complex isolation			
μ MACS Protein A MicroBeads	Immunopurification of proteins	20-40	130-071-001
μ MACS Protein G MicroBeads	Immunopurification of proteins	20-40	130-071-101
μ MACS Epitope Tag Protein Isolation Kits	Isolation of His-, HA-, c-myc-, GST-, and GFP-tagged proteins	40	
μ MACS Streptavidin Kit	Isolation of biotinylated molecules and their targets	20	130-074-101
μ MACS FactorFinder Kit	Isolation of transcription factors	20	130-092-317

16

11. Technical data and specifications

- Size (thermoMACS Separator): 191×40×33 mm
- Weight (thermoMACS Separator): 500 g
- Input voltage of the thermoMACS Separator: 12 V
Power consumption: 52 W
- Input voltage of the power supply (enclosed with the thermoMACS Separator): 100-240 VAC, ~ 50-60 Hz
Power consumption: 52 W
- ▲ Note:** Use the thermoMACS Separator only with the enclosed power supply. Use of a non-specified power supply may damage the device.
- The thermoMACS Separator is labeled as a protection class III device.
- Conditions of operation: 10 °C to 30 °C with 0% to 85% humidity.
- μ Column reaction temperatures: 37 °C, 42 °C, room temperature (heat off).
- The thermoMACS Separator can process up to four samples in parallel.
- Standards: The thermoMACS Separation Unit meets the intent of the EMC Directive 89/336/EEC and the Low Voltage Directive (LVD) 73/23/EEC. Compliance was demonstrated by conformance to the following specifications which have been listed in the Official Journal of the European Communities.

thermoMACS™ Separator: Power Supply:

EMC: EN 55022, Class B; ICE 61000-4-2 ICE 61000-4-3	EMC: EN 55022, Class B EN 61000-3-2 EN 61000-3-3 EN 61000-4-2 EN 61000-4-3 EN 61000-4-4 EN 61000-4-6 EN 61000-4-11
Safety: EN 61010/A2	Safety: UL 1950 CSA 22.2 EN 60950

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the Technical Specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS is a registered trademark of Miltenyi Biotec GmbH.

thermoMACS and μ MACS are trademarks of Miltenyi Biotec GmbH.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

© 2007 Miltenyi Biotec GmbH. Printed in Germany.

17

18