

## Contents

1. Description
  - 1.1 Background information
  - 1.2 Applications
  - 1.3 Reagent and instrument requirements
2. Protocol
  - 2.1 Reagent preparation
  - 2.2 Automated dissociation using the gentleMACS™ Dissociator
    - 2.2.1 Automated dissociation of up to 7 days old EBs
    - 2.2.2 Automated dissociation of 8–20 days old EBs
  - 2.3 Manual dissociation
    - 2.3.1 Manual dissociation of up to 7 days old EBs
    - 2.3.2 Manual dissociation of 8–20 days old EBs

## 1. Description

<b>Components</b>	6 vials, containing 2.5 mL Enzyme P 2×50 mL Buffer X 1.5 mL Buffer Y 1 vial Enzyme A (lyophilized powder) 1 mL Buffer A
<b>Size</b>	For 50 digestions.  The specified number of digestions is valid when embryoid bodies at several stages generated from a starting cell number of $1 \times 10^6/10$ mL (mouse) and $4 \times 10^6/10$ mL (human) are used.
<b>Storage</b>	Upon arrival immediately store Enzyme P in aliquots at $-20^\circ\text{C}$ . Store all other components at $2-8^\circ\text{C}$ upon arrival. Reconstitute Enzyme A before the date indicated on the box label. For information about reconstitution and storage after reconstitution of the lyophilized component refer to chapter 2.1.

## 1.1 Background information

Embryoid body (EB) formation is a crucial step in many ES or iPS cell differentiation protocols. For subsequent cell analysis or isolation of specific cell populations of the EB, single-cell suspensions are a prerequisite, for example, to achieve the highest possible purity and recovery during cell separations with MACS® Technology. The gentleMACS™ Dissociator provides optimized programs to attain single-cell suspensions from *in vitro* generated EBs.

EBs are cellular aggregates generated from differentiating pluripotent stem cells grown in suspension in defined cell numbers and media. Cell aggregation is initiated by suspension culture in hanging drops or as mass suspension culture in non-tissue culture plates or spinner flasks. During the three dimensional differentiation process, *in vivo* embryonic development is recapitulated to a limited extent and differentiated cell types of all three germ layers are generated. Embryoid bodies are a well studied model system for investigating early stages of development and deriving distinct cell types from pluripotent stem cells.

## 1.2 Applications

- Dissociation of EBs into single-cell suspensions for subsequent cell separations using MACS Technology.
- Cultivation of EB-derived cells.
- Phenotyping or enumeration of individual EB-derived cell populations by flow cytometry.
- Dissociation of pluripotent stem cell-derived neurospheres.

▲ **Note:** The Embryoid Body Dissociation Kit is not suitable for the dissociation of tissue-derived neurospheres. In this case use the Neurosphere Dissociation Kit (P) (# 130-095-943) or the Neurosphere Dissociation Kit (T) (# 130-095-944) depending on the antigen epitope of interest.

## 1.3 Reagent and instrument requirements

- Buffer: Dulbecco's phosphate-buffered saline (DPBS) without  $\text{Ca}^{2+}$  and  $\text{Mg}^{2+}$
- EB culture medium
- Pre-Separation Filters, 70  $\mu\text{m}$  (# 130-095-823)
- (Optional) gentleMACS Dissociator (# 130-093-235)
- (Optional) gentleMACS M Tubes (# 130-093-236, # 130-096-335)
- (Optional) ART® 1000 REACH™ pipet tips (Molecular BioProducts, Inc.) for removal of dissociated material from the closed M Tubes.

## 2. Protocol

▲ The dissociation protocols have been optimized for embryoid bodies generated from a starting cell number of  $1 \times 10^6/10$  mL (mouse) and  $4 \times 10^6/10$  mL (human).

▲ For cell culture experiments subsequent to embryoid body dissociation, all steps should be performed under sterile conditions.

▲ Partially dissociated EBs may adhere unspecifically to plastic surfaces such as the tube walls. Thus, avoid unnecessary rotation of tubes during dissociation and keep the EBs covered with enzyme solution at all times.

### 2.1 Reagent preparation

1. Prepare Enzyme A by reconstitution of the lyophilized powder in its vial with 1 mL of Buffer A. Do not vortex. Prepare aliquots of appropriate volume and avoid repeated freeze-thaw-cycles. Store aliquots at  $-20^\circ\text{C}$ . This solution is stable for 6 months after reconstitution. For cell culture experiments subsequent to tissue dissociation, Enzyme A should be sterile filtered.

2. Prepare 1950  $\mu\text{L}$  enzyme mix 1 by adding 50  $\mu\text{L}$  of Enzyme P to 1900  $\mu\text{L}$  of Buffer X and vortex. Pre-heat the mixture to  $37^\circ\text{C}$  for 10–15 minutes before use.

▲ **Note:** For enumeration or magnetic cell separation of CD309<sup>+</sup> (KDR<sup>+</sup>, Flk-1<sup>+</sup>, VEGFR2) cells add 10  $\mu\text{L}$  of Enzyme P to 1940  $\mu\text{L}$  of Buffer X.

3. Prepare 30  $\mu\text{L}$  enzyme mix 2 by adding 20  $\mu\text{L}$  of Buffer Y to 10  $\mu\text{L}$  of Enzyme A.
4. Add 30  $\mu\text{L}$  enzyme mix 2 to 1950  $\mu\text{L}$  enzyme mix 1.

Enzyme mix 1		Enzyme mix 2	
Enzyme P 50 $\mu\text{L}$	Buffer X 1900 $\mu\text{L}$	Buffer Y 20 $\mu\text{L}$	Enzyme A 10 $\mu\text{L}$

### 2.2 Automated dissociation using the gentleMACS™ Dissociator

▲ For details on the use of the gentleMACS™ Dissociator, refer to the gentleMACS Dissociator user manual. The protocol is not compatible with the gentleMACS Octo Dissociator or the gentleMACS Octo Dissociator with Heaters.

#### 2.2.1 Automated dissociation of up to 7 days old EBs

1. Transfer up to 20 mL of EB-containing cell culture medium into the gentleMACS M Tube. Centrifuge at  $300 \times g$  for 2 minutes. Discard supernatant.
2. Add 10 mL of DPBS to the gentleMACS M Tube. Resuspend EBs carefully. Centrifuge at  $300 \times g$  for 2 minutes. Discard supernatant.
3. Add 1980  $\mu\text{L}$  of the combined pre-heated enzyme mix 1 and 2 to the gentleMACS M Tube and tightly close the lid.
4. Incubate sample for 10 minutes at  $37^\circ\text{C}$  without agitation and ensure EBs are covered by the enzyme solution.
5. Invert the M Tube in one movement and attach it upside down onto the sleeve of the gentleMACS Dissociator. Keep the tube upside down during the following incubation and dissociation steps 6–9.
6. Run the gentleMACS Program **EB\_01**.
7. Once the program is finished incubate sample for 5 minutes at  $37^\circ\text{C}$  without agitation.

8. Attach M Tube upside down onto the sleeve of the gentleMACS Dissociator.
9. Run the gentleMACS Program **EB\_02**.
10. Once the program is finished add 8 mL of cell culture medium or DPBS to the M Tube.
11. Resuspend sample and apply the cell suspension to a Pre-Separation Filter,  $70 \mu\text{m}$ , placed on a suitable (13–15 mL) tube.  
▲ **Note:** Dissociated EBs can be removed from the closed M Tube by pipetting through the septum-sealed opening in the center of the cap of the M Tube. Use ART 1000 REACH 1000  $\mu\text{L}$  pipette tips.
12. Wash Pre-Separation Filter,  $70 \mu\text{m}$ , with 3 mL of EB culture medium or DPBS.
13. Discard filter and centrifuge cell suspension at  $300 \times g$  for 5 minutes. Aspirate supernatant completely.
14. Resuspend cells with appropriate buffer or medium to the required volume for further applications.

#### 2.2.2 Automated dissociation of 8–20 days old EBs

1. Transfer up to 20 mL of EB-containing cell culture medium into the gentleMACS M Tube. Centrifuge at  $300 \times g$  for 2 minutes. Discard supernatant.
2. Add 10 mL of DPBS to the gentleMACS M Tube. Resuspend EBs carefully. Centrifuge at  $300 \times g$  for 2 minutes. Discard supernatant.
3. Add 1980  $\mu\text{L}$  of the combined pre-heated enzyme mix 1 and 2 to the gentleMACS M Tube and tightly close the lid.
4. Incubate sample for 10 minutes at  $37^\circ\text{C}$  without agitation and ensure EBs are covered by the enzyme solution.
5. Invert the M Tube in one movement and attach it upside down onto the sleeve of the gentleMACS Dissociator. Keep the tube upside down during the following incubation and dissociation steps 6–12.
6. Run the gentleMACS Program **EB\_01**.
7. Once the program is finished incubate sample for 10 minutes at  $37^\circ\text{C}$  without agitation.
8. Attach M Tube upside down onto the sleeve of the gentleMACS Dissociator.
9. Run the gentleMACS Program **EB\_02**.
10. Once the program is finished incubate sample for 5 minutes at  $37^\circ\text{C}$  without agitation.
11. Attach M Tube upside down onto the sleeve of the gentleMACS Dissociator.
12. Run the gentleMACS Program **EB\_02**.
13. Once the program is finished add 8 mL of cell culture medium or DPBS to the M Tube.
14. Resuspend sample and apply the cell suspension to a Pre-Separation Filter,  $70 \mu\text{m}$ , placed on a suitable (13–15 mL) tube.  
▲ **Note:** Dissociated EBs can be removed from the closed M Tube by pipetting through the septum-sealed opening in the center of the cap of the M Tube. Use ART 1000 REACH 1000  $\mu\text{L}$  pipette tips.

15. Wash Pre-Separation Filter, 70  $\mu$ m, with 3 mL of cultivation medium or DPBS.
16. Discard filter and centrifuge cell suspension at 300 $\times$ g for 5 minutes. Aspirate supernatant completely.
17. Resuspend cells with appropriate buffer or medium to the required volume for further applications.
7. Pipette up and down for 1 minute using a 1000  $\mu$ L pipette.
8. Incubate for 5 minutes at 37 °C without agitation.
9. Pipette up and down for 1 minute using a 1000  $\mu$ L pipette.
10. Add 8 mL of cell culture medium or DPBS.
11. Pass the cell suspension through a Pre-Separation Filter, 70  $\mu$ m, placed on a suitable (13–15 mL) tube.
12. Wash the Pre-Separation Filter with 3 mL of EB culture medium or DPBS.
13. Discard filter and centrifuge cell suspension at 300 $\times$ g for 5 minutes. Aspirate supernatant completely.
14. Resuspend in an appropriate buffer or medium required for your downstream application.

### 2.3 Manual dissociation

▲ The protocol below is compatible with a wide range of EB-based differentiation protocols. However, depending on cell type of interest and the protocol used for EB generation, pipetting and incubation times may have to be optimized.

▲ Dissociation by pipetting may give variable results, depending on the user. For optimal standardization, the use of the gentleMACS Dissociator is recommended (refer to section 2.2)

#### 2.3.1 Manual dissociation of up to 7 days old EBs

1. Transfer up to 20 mL of EB-containing cell culture medium into a 50 mL tube. Centrifuge at 300 $\times$ g for 2 minutes. Discard supernatant.
2. Add 10 mL of DPBS and resuspend EBs carefully. Centrifuge at 300 $\times$ g for 2 minutes. Discard supernatant.
3. Add 1980  $\mu$ L of the combined pre-heated enzyme mix 1 and 2.
4. Incubate for 10 minutes at 37 °C without agitation and ensure EBs are covered by the enzyme solution.
5. Start mechanical dissociation by carefully pipetting up and down using a 1000  $\mu$ L pipette for 1 minute.
6. Incubate for 5 minutes at 37 °C without agitation.
7. Pipette up and down for 1 minute using a 1000  $\mu$ L pipette.
8. Add 8 mL of cell culture medium or DPBS.
9. Pass the cell suspension through a Pre-Separation Filter, 70  $\mu$ m, placed on a suitable (13–15 mL) tube.
10. Wash the Pre-Separation Filter with 3 mL of EB culture medium or DPBS.
11. Discard filter and centrifuge cell suspension at 300 $\times$ g for 5 minutes. Aspirate supernatant completely.
12. Resuspend in an appropriate buffer or medium required for your downstream application.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

#### 2.3.2 Manual dissociation of 8–20 days old EBs

1. Transfer up to 20 mL of EB-containing cell culture medium into a 50 mL tube. Centrifuge at 300 $\times$ g for 2 minutes. Discard supernatant.
2. Add 10 mL of DPBS and resuspend EBs carefully. Centrifuge at 300 $\times$ g for 2 minutes. Discard supernatant.
3. Add 1980  $\mu$ L of the combined pre-heated enzyme mix 1 and 2.
4. Incubate sample for 10 minutes at 37 °C without agitation and ensure EBs are covered by the enzyme solution.
5. Start mechanical dissociation by carefully pipetting up and down using a 1000  $\mu$ L pipette for 1 minute.
6. Incubate sample for 10 minutes at 37 °C without agitation.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

gentleMACS, MACS, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.