

## Contents

1. Description
  - 1.1 Principle of the Umbilical Cord Dissociation Kit
  - 1.2 Background information
  - 1.3 Applications
  - 1.4 Reagent and instrument requirements
2. Protocols
  - 2.1 Reagent preparation
  - 2.2 Umbilical cord dissociation protocol
    - 2.2.1 Dissociation of frozen tissue
    - 2.2.2 Dissociation of fresh tissue

## 1. Description

This product is for research use only.

<b>Components</b>	6 vials, containing <b>1 vial of Enzyme D (lyophilized powder)</b> <b>1 vial of Enzyme B (lyophilized powder)</b> <b>1 vial of Enzyme A (lyophilized powder)</b> <b>1 vial of 2.5 mL Enzyme P</b> <b>2 vials of 30 mL Buffer L</b>
<b>Size</b>	For 25 digestions. The specified number of digestions is valid when digesting umbilical cord pieces that do not exceed 0.5 g in weight following the protocols in chapter 2.2.
<b>Storage</b>	Upon arrival immediately store Enzyme P in aliquots at -20 °C. Store all other components at 2-8 °C upon arrival. Reconstitute all components before the date indicated on the box label. For information about reconstitution and storage after reconstitution of the lyophilized components refer to chapter 2.1.

### 1.1 Principle of the Umbilical Cord Dissociation Kit

Human umbilical cord can be dissociated to single-cell suspensions by combining mechanical dissociation with enzymatic degradation of the extracellular adhesion proteins which maintain the structural integrity of tissues.

In a first step, the umbilical cord is soaked with enzymes which degrade the extracellular matrix. In a second step, single-cells are freed from the extracellular matrix by using the gentleMACS™ Dissociators. After dissociation remaining particles are removed by filtration using MACS® SmartStrainer.

Cells should be processed immediately for downstream applications, such as cell culture, cell separation, or cell analysis.

### 1.2 Background information

The Umbilical Cord Dissociation Kit, human enables the gentle and efficient generation of single-cell suspensions from human umbilical cord. The kit is not certified as GMP product and has been particularly developed to enable a quality control of frozen umbilical cord samples for umbilical cord tissue banks by detecting mesenchymal stem cells (MSCs) directly after dissociation by flow cytometry without long-standing cultivation of the cells. Additionally, the cells can be also used in downstream applications such as cell culture or molecular studies.

### 1.3 Applications

- Dissociation of fresh or frozen human umbilical cord pieces for subsequent cell separations using MACS Technology.
- Cultivation of umbilical cord cell populations, e.g., mesenchymal stem cells or endothelial cells.
- Phenotyping or enumeration of umbilical cord cell populations by flow cytometry or fluorescence microscopy.

### 1.4 Reagent and instrument requirements

- PBS: phosphate-buffered saline pH 7.2
- DMEM (# 130-091-437)
- MACS SmartStrainers (100 µm) (# 130-098-463)
- gentleMACS Dissociator (# 130-093-235), gentleMACS Octo Dissociator (# 130-095-937), or gentleMACS Octo Dissociator with Heaters (# 130-096-427)
- gentleMACS C Tubes (# 130-093-237, # 130-096-334)
- (Optional) 70 % ethanol
- (Optional) MACS Tissue Storage Solution (# 130-100-008)
- (Optional) ART® 1000 REACH™ pipet tips (Molecular BioProducts, Inc.) for removal of dissociated material from the closed C Tubes.
- (Optional for expansion of MSCs) StemMACS™ MSC Expansion Media Kit XF, human (# 130-104-182) or StemMACS MSC Expansion Media, human (# 130-091-680)
- (Optional for differentiation of MSCs) StemMACS AdipoDiff Media, human (# 130-091-677), StemMACS ChondroDiff Media, human (# 130-091-679), or StemMACS OsteoDiff Media, human (# 130-091-678)

## 2. Protocols

▲ For details on the use of the gentleMACS Dissociators, refer to the gentleMACS Dissociator user manuals.

▲ For cell culture experiments subsequent to tissue dissociation, all steps should be performed under sterile conditions.

### 2.1 Reagent preparation

1. Prepare aliquots of appropriate volume of Enzyme P to avoid repeated freeze-thaw-cycles. Store aliquots at  $-20^{\circ}\text{C}$ . This solution is stable for 6 months.
2. Prepare Enzyme D by reconstitution of the lyophilized powder in the vial with 3 mL Buffer L supplied with the kit. Do not try to resuspend by pipetting or vortexing. Invert vial after closing and wait 5–10 minutes to dissolve the pellet. Prepare aliquots of appropriate volume to avoid repeated freeze-thaw-cycles. Store aliquots at  $-20^{\circ}\text{C}$ . This solution is stable for 6 months after reconstitution. For cell culture experiments subsequent to tissue dissociation, Enzyme D should be sterile filtered prior to aliquoting.
3. Prepare Enzyme A by reconstitution of the lyophilized powder in the vial with 1 mL Buffer L supplied with the kit. Do not vortex. Prepare aliquots of appropriate volume to avoid repeated freeze-thaw-cycles. Store aliquots at  $-20^{\circ}\text{C}$ . This solution is stable for 6 months after reconstitution.
4. Prepare Enzyme B by reconstitution of the lyophilized powder in the vial with 1 mL Buffer L supplied with the kit. Do not vortex. Prepare aliquots of appropriate volume to avoid repeated freeze-thaw-cycles. Store aliquots at  $-20^{\circ}\text{C}$ . This solution is stable for 6 months after reconstitution. For cell culture experiments subsequent to tissue dissociation, Enzyme B should be sterile filtered prior to aliquoting.

### 2.2 Umbilical cord dissociation protocol

#### 2.2.1 Dissociation of frozen tissue

1. Cryopreservation:  
Wash umbilical cord in an appropriate buffer. Cut umbilical cord into small pieces (<5 mm diameter) to allow quick penetration of cryopreservation buffer. Add a cryoconservation medium containing 10% DMSO and transfer sample in a cryopreservation tube or bag. Freeze sample slowly at  $-70^{\circ}\text{C}$ , e.g., by using an appropriate styrofoam box. After 24 hours transfer sample into a liquid nitrogen tank.
2. Thawing of samples:  
Switch on a water bath at  $37^{\circ}\text{C}$ . Thaw the sample (about 0.5 g of tissue) completely using the water bath. Transfer the sample into a tube containing 45 mL of DMEM. Invert the tube for 1 minute and transfer the sample onto a MACS SmartStrainer (100  $\mu\text{m}$ ).
3. Add 2.2 mL of Buffer L into the gentleMACS C Tube.
4. Transfer the retained sample from the strainer into the C Tube.
5. Add 100  $\mu\text{L}$  of Enzyme D, 62.5  $\mu\text{L}$  of Enzyme P, 4  $\mu\text{L}$  of Enzyme B, and 10  $\mu\text{L}$  of Enzyme A into the C Tube.  
▲ **Note:** Enzyme A, B, and D can be premixed before addition into the C Tube. Do not premix any of the enzymes with Enzyme P.

6. Tightly close C Tube beyond the first resistance.  
If using the heating function of the gentleMACS Octo Dissociator with Heaters, attach the C Tube upside down onto the sleeves of the gentleMACS Octo Dissociator with Heaters. It has to be ensured that the sample material is still located in the area of the rotor/stator before starting the gentleMACS program. Run the gentleMACS Program **37C\_h\_UCDK\_1** and continue with step 10.
7. Incubate sample in a water bath at  $37^{\circ}\text{C}$  for 3 hours.
8. Attach C Tube upside down onto the sleeve of the gentleMACS Dissociator.  
▲ **Note:** It has to be ensured that the sample material is still located in the area of the rotor/stator before starting the gentleMACS program.
9. Run the gentleMACS program **h\_cord\_01**.
10. After termination of the program, detach C Tube from the gentleMACS Dissociator.
11. Perform a short centrifugation step to collect the sample material at the tube bottom.
12. Apply sample to a MACS SmartStrainer (100  $\mu\text{m}$ ).  
▲ **Note:** Dissociated tissue can be removed from the closed C Tube by pipetting through the septum-sealed opening in the center of the cap of the C Tube. Use ART 1000 REACH 1000  $\mu\text{L}$  pipette tips.
13. Wash the filter with 5 mL of DMEM.  
▲ **Note:** For maximum cell recovery, rinse the used C Tube with wash buffer before transfer to the filter.
14. Discard the filter, invert tube after closing and centrifuge sample at  $300\times g$  for 15 minutes at room temperature.
15. Aspirate supernatant completely and resuspend cells with medium or an appropriate buffer to the required volume for further applications.
16. Process cells immediately for further applications.  
▲ **Note:** MSCs, endothelial cells, and leukocytes can be distinguished after dissociation by the surface markers CD90, CD73, CD34, and CD45 (MSCs are CD90<sup>+</sup>CD73<sup>+</sup>CD34<sup>+</sup>CD45<sup>-</sup>; endothelial cells are CD34<sup>+</sup>CD73<sup>+</sup>CD45<sup>-</sup>CD90<sup>-</sup>). If using the MSC Phenotyping Kit, be aware that CD105 is degraded during the dissociation step and that endothelial cells cannot be distinguished from leukocytes by this kit.  
▲ **Note:** For expansion of MSCs we recommend using the StemMACS MSC Expansion Media Kit XF, human (# 130-104-182).

#### 2.2.2 Dissociation of fresh tissue

1. Wash umbilical cord in an appropriate buffer or cell culture medium, e.g., MACS Tissue Storage Solution or DMEM.  
▲ **Note:** For subsequent cell culture soak umbilical cord in 70 % ethanol for 30 seconds, then cut off the ends, and wash at least once in an appropriate buffer. Remove as much blood as possible.
2. Add Buffer L (volume according to the table in step 5) into the gentleMACS C Tube.
3. Cut off a piece of the umbilical cord (weight according to the table in step 5) and transfer it into the C Tube.
4. Cut the tissue in small pieces of 2–4 mm using scissors which reach the bottom of the C Tube.

5. Add enzymes into the C Tube according to the following table:

Tissue [g]	Buffer L [mL]	Enzyme D [μL]	Enzyme P [μL]	Enzyme B [μL]	Enzyme A [μL]
0.5	2.2	100	62.5	4	10
1	4.4	200	125	8	20
2-4	8.8	400	250	16	40

▲ **Note:** Enzyme A, B, and D can be premixed before addition into the C Tube. Do not premix any of the enzymes with Enzyme P.

6. Tightly close C Tube beyond the first resistance.  
If using the heating function of the gentleMACS Octo Dissociator with Heaters, attach the C Tube upside down onto the sleeves of the gentleMACS Octo Dissociator with Heaters. Run the gentleMACS Program **37C\_h\_UCDK\_1** and continue with step 10.
7. Incubate sample in a water bath at 37 °C for 3 hours.
8. Attach C Tube upside down onto the sleeve of the gentleMACS Dissociator.  
▲ **Note:** It has to be ensured that the sample material is still located in the area of the rotor/stator before starting the gentleMACS program.
9. Run the gentleMACS program **h\_cord\_01**.
10. After termination of the program, detach C Tube from the gentleMACS Dissociator.
11. Perform a short centrifugation step to collect the sample material at the tube bottom.
12. Apply sample to a MACS SmartStrainer (100 μm).  
▲ **Note:** Dissociated tissue can be removed from the closed C Tube by pipetting through the septum-sealed opening in the center of the cap of the C Tube. Use ART 1000 REACH 1000 μL pipette tips.
13. Wash the filter with DMEM. Use twice as much DMEM as volume in tube.  
▲ **Note:** For maximum cell recovery, rinse the used C Tube with wash buffer before transfer to the filter.
14. Discard the filter, invert tube after closing and centrifuge sample at 300×g for 15 minutes at room temperature.
15. Aspirate supernatant completely and resuspend cells with medium or an appropriate buffer to the required volume for further applications.
16. Process cells immediately for further applications.  
▲ **Note:** MSCs, endothelial cells, and leukocytes can be distinguished after dissociation by the surface markers CD90, CD73, CD34, and CD45 (MSCs are CD90<sup>+</sup>CD73<sup>+</sup>CD34<sup>+</sup>CD45<sup>-</sup>; endothelial cells are CD34<sup>+</sup>CD73<sup>+</sup>CD45<sup>-</sup>CD90<sup>-</sup>). If using the MSC Phenotyping Kit, be aware that CD105 is degraded during the dissociation step and that endothelial cells cannot be distinguished from leukocytes by this kit.  
▲ **Note:** For expansion of MSCs we recommend using the StemMACS MSC Expansion Media Kit XF, human (# 130-104-182).

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for local Miltenyi Biotec Technical Support contact information.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

gentleMACS, MACS, the Miltenyi Biotec logo, and StemMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.