

Introduction

Chimeric antigen receptors (CARs) have been used to direct autologous T cells to hematopoietic cancers with promising clinical results. The use of allogeneic T cells however comes with the risk of inducing graft-versus-host disease (GVHD). In this regard, CAR-expressing natural killer (CAR-NK) cells provide potential advantages: NK cells can be used in an allogeneic setting as they pose a low risk of GVHD. Moreover, besides the CAR-guided antitumor action, NK cells exert their activity

through their native receptors, potentially preventing cancer cell antigen escape. *Ex vivo* modification of NK cells to induce CAR expression requires isolation, activation, transduction, and expansion steps. To simplify procedures and widen the applicability for future clinical therapies, we are developing automated complex cell manufacturing processes in a closed system.

Methods

1 CD3⁺CD56⁺ NK cell isolation in the CliniMACS Prodigy®

We developed a fully automated procedure for the isolation of NK cells from leukapheresis products, the CliniMACS Prodigy CD3⁺ cell depletion/CD56⁺ cell enrichment process, enabling a two-step separation in a single tubing set

(CliniMACS Prodigy TS 310, fig. 1). The process is designed in distinct blocks, enabling the use of both separation steps individually. NKCE: natural killer cell engineering (see results section 3).




Figure 1

2 Enrichment of CD56⁺ cells and cultivation of NK cells in the CliniMACS Prodigy®

For cell cultivation and transduction, we developed the CliniMACS Prodigy CD56⁺ cell enrichment and cultivation process providing fully automated

CD56⁺ cell enrichment and subsequent expansion of NK cells. The entire process takes place in a single tubing set (CliniMACS Prodigy TS 520, fig. 2).




Figure 2

Results

1 CD3⁺CD56⁺ NK cell isolation in the CliniMACS Prodigy®

The CliniMACS Prodigy CD3⁺ cell depletion/CD56⁺ cell enrichment process resulted in a 4 log (average) depletion of T cells (fig. 3A) and a 98% purity of CD3⁺CD56⁺ cells. The resulting NK cell product

contained about 2x10⁸ NK cells (fig. 3B). NK cell recovery was in the range of 60% (fig. 3C) and the NK:T cell ratio amounted to 1,400 on average (fig. 3D).




Figure 3

2 Selective expansion of cytotoxic NK cells in NK MACS® GMP Medium

NK cells were expanded from PBMCs (n = 3) in the presence of 5% AB serum, 500 IU/mL IL-2, and 140 U/mL IL-15 for 14 days. Expansion rates were high for both NK MACS Medium (for research use) and NK MACS GMP Medium (fig. 4A). Both media supported selective NK cell growth, which resulted in a mean purity of approximately 80% NK cells

after 14 days of culture (fig. 4B). Starting NK cell frequency in PBMCs amounted to 15% on average. Expanded NK cells showed high cytotoxicity against K-562 cells at different effector-to-target cell ratios after 14 days of culture in both NK MACS Medium (research use) and NK MACS GMP Medium (fig. 4C).




Figure 4

3 CD3⁺CD56⁺ NK cell isolation and cultivation generates highly pure NK cells providing the option of genetic engineering

To provide highly pure genetically engineered NK cells, a process for natural killer cell engineering (NKCE) is currently being developed. In a first approach, CD3⁺ cell depletion will be performed in a separate process (fig. 1, part 1 of NKCE). The resulting cell product, which is depleted of CD3⁺ cells, will be the starting material for automated CD56⁺ cell enrichment, transduction, and expansion (fig. 2, part 2 of NKCE). Figure 5A exemplifies the results for

cell separation and cultivation of a sample from one donor. After separation, an average NK cell purity of 97% was achieved with a T cell depletion of 4.4–4.9 log (NK:T cell ratio of 1,949–5,591). Figure 5B shows the NK cell expansion over time, resulting in an average NK cell number of 2.4x10⁸ with a purity of >99.8% (fig. 5A). This process represents the basis for a fully automated platform enabling genetic engineering of NK cells.




Figure 5

Conclusion

- CD3⁺CD56⁺ NK cells can be enriched to high purity in the CliniMACS Prodigy.
- NK MACS GMP Medium enables high NK cell expansion rates for clinical settings.
- NK cells can be expanded at high purity (>99%) in the CliniMACS Prodigy.
- These fully automated processes provide the basis for natural killer cell engineering in a closed system. This will simplify CAR NK cell manufacturing procedures and therefore widen the applicability in future clinical applications.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS® GMP Products are for research use and *ex vivo* cell culture processing only, and are not intended for human *in vivo* applications. For regulatory status in the USA, please contact your local representative. MACS GMP Products are manufactured and tested under a quality system certified to ISO 13485 and are in compliance with relevant GMP guidelines. They are designed following the recommendations of USP <1043> on ancillary materials. The CliniMACS® System components, including Reagents, Tubing Sets, Instruments, and PBS/EDTA Buffer, are designed, manufactured and tested under a quality system certified to ISO 13485. In the EU, the CliniMACS System components are available as CE-marked medical devices for their respective intended use, unless otherwise stated. The CliniMACS Reagents and Biotin Conjugates are intended for *in vitro* use only and are not designated for therapeutic use or direct infusion into patients. The CliniMACS Reagents in combination with the CliniMACS System are intended to separate human cells. Miltenyi Biotec as the manufacturer of the CliniMACS System does not give any recommendations regarding the use of separated cells for therapeutic purposes and does not make any claims regarding a clinical benefit. For the manufacturing and use of target cells in humans the national legislation and regulations – e.g. for the EU the Directive 2004/23/EC (“human tissues and cells”), or the Directive 2002/98/EC (“human blood and blood components”) – must be followed. Thus, any clinical application of the target cells is exclusively within the responsibility of the user of a CliniMACS System. In the US, the CliniMACS CD34 Reagent System, including the CliniMACS Plus Instrument, CliniMACS CD34 Reagent, CliniMACS Tubing Sets TS and LS, and the CliniMACS PBS/EDTA Buffer, is FDA approved; all other products of the CliniMACS Product Line are available for use only under an approved Investigational New Drug (IND) application or Investigational Device Exemption (IDE). CliniMACS MicroBeads are for research use only and not for human therapeutic or diagnostic use. CliniMACS, CliniMACS Prodigy, MACS, and the MACS logo are registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. Copyright © 2018 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.