



# Sensitive analysis and isolation of ROR1<sup>+</sup> B cells for research on chronic lymphocytic leukemia

## Introduction

### The challenge of analyzing lymphoid cells in chronic lymphocytic leukemia research

Flow cytometric analysis of lymphoid cells from chronic lymphocytic leukemia (CLL) samples requires a complex panel of multiple markers. Most of these markers are expressed on both normal and malignant cells, and only characteristic marker expression profiles help distinguish between normal and malignant cells. Identification of new markers and the optimization of antibody panels for flow cytometric immunophenotyping would greatly enhance leukemia research.

Analysis of malignant B cells becomes especially difficult with samples containing only very low amounts of these cells, e.g., in minimal residual disease (MRD). Enrichment of malignant B cells from these samples would help overcome this obstacle. However, current techniques for the isolation of B cells do not allow for a distinction between malignant and normal cells. This drawback hampers the analysis of differential gene expression, for example.

### ROR1 as a B cell marker in CLL research


ROR1 is a member of the receptor tyrosine kinase-like orphan receptor (ROR) family and has characteristics of an oncofetal antigen. ROR1 signaling increases cell survival via the wnt pathway<sup>1</sup>. The protein is expressed on B cells from CLL blood samples, but not from blood of healthy donors.<sup>1,2,3</sup> ROR1 is also expressed on various solid tumors<sup>4</sup>, but not in major adult tissues apart from low levels in adipose tissue, on embryonic stem cells, and transiently at an early stage of B cell development<sup>5</sup>. Therefore, ROR1 is a potential marker to identify and isolate CLL B cells from blood samples for downstream analysis.

## Tools for the analysis and isolation of ROR1<sup>+</sup> B cells from CLL samples

### Analysis of ROR1<sup>+</sup> B cells in CLL samples

Miltenyi Biotec offers an Anti-ROR1 antibody, which increases sensitivity and specificity of flow cytometric B cell analysis in CLL research. It is available as a conjugate with PE, APC, or biotin. For functional studies, the antibody


is provided in an unconjugated form, without toxic preservatives or endotoxins. Figure 1 shows a typical immunophenotyping experiment with a CLL sample. PBMCs were labeled with a panel of eight MACS® Antibodies and analyzed by flow cytometry using the following setup: violet laser (405 nm): CD20-VioBlue®; blue laser (488 nm): CD43-FITC, CD5-PE, CD79b-PerCP-Vio700™, CD19-PE-Vio770™; red laser (633 nm): Anti-ROR1-APC, CD3-APC-Vio770.



**Figure 1: Identification of ROR1<sup>+</sup> cells and other immune cells in PBMCs from a CLL sample.** Cells were labeled with CD20-VioBlue, CD43-FITC, CD5-PE, CD79b-PerCP-Vio700, CD19-PE-Vio770, Anti-ROR1-APC, and CD3-APC-Vio770, and analyzed on the MACSQuant Analyzer 10. CD19<sup>+</sup>ROR1<sup>+</sup> B cells are shown in magenta, CD19<sup>+</sup>ROR1<sup>-</sup> B cells in light blue, T cells in yellow, and NK cells in black.

## Isolation of ROR1<sup>+</sup> cells from CLL samples

For the isolation of ROR1<sup>+</sup> B cells, Miltenyi Biotec has developed the Anti-ROR1 MicroBead Kit. Using this kit, ROR1<sup>+</sup> cells are indirectly magnetically labeled with an Anti-ROR1-PE antibody and Anti-PE MicroBeads. Subsequently, the labeled ROR1<sup>+</sup> cells are isolated manually using LS or MS Columns, or automatically with the autoMACS<sup>®</sup> Pro Separator. The isolated ROR1<sup>+</sup> B cells are then ready for any downstream analysis. ROR1<sup>+</sup> cells can be enriched to high purities and numbers suitable for further research processing – even from samples with very low ROR1<sup>+</sup> cell frequencies, e.g., in MRD (fig. 2).



**Figure 2: Enrichment of ROR1<sup>+</sup> cells from PBMCs of a donor with MRD.** ROR1<sup>+</sup> cells were isolated using the Anti-ROR1 MicroBead Kit, either manually (A) or automatically using the autoMACS Pro Separator (B). Cells were fluorescently labeled with Anti-ROR1-PE and CD19-APC and analyzed by flow cytometry. Purities of the manually and automatically isolated cells amounted to 94% and 98%, and yields amounted to 45% and 97%, respectively.

## Conclusion

- Anti-ROR1 antibodies enhance the sensitivity and specificity of flow cytometric immunophenotyping of CLL samples.
- Anti-ROR1 antibodies allow the identification and enumeration of B cells for CLL research.
- The Anti-ROR1 MicroBead Kit enables the enrichment of ROR1<sup>+</sup> CLL B cells – even from samples containing only very low amounts of these cells, for example, in MRD.
- Isolation of ROR1<sup>+</sup> cell populations to high purities facilitates research on the specific features and functions of CLL B cells.

## References

1. Fukuda, T. *et al.* (2008) Antisera induced by infusions of autologous Ad-CD154-leukemia B cells identify ROR1 as an oncofetal antigen and receptor for Wnt5a. *Proc. Acad. Natl. Sci USA* 105: 3047–3052.
2. Baskar, S. *et al.* (2012) Targeting malignant B cells with an immunotoxin against ROR1. *MAbs* 4: 349–361.
3. Baskar, S. *et al.* (2008) Unique cell surface expression of receptor tyrosinase kinase ROR1 in human B-cell chronic lymphocytic leukemia. *Clin. Cancer Res.* 14: 396–404.
4. Zhang, S. *et al.* (2012) The onco-embryonic antigen ROR1 is expressed by a variety of human cancers. *Am. J. Pathol.* 181: 1903–1910.
5. Hudecek, M. *et al.* (2010) The B-cell tumor-associated antigen ROR1 can be targeted with T cells modified to express a ROR1-specific chimeric antigen receptor. *Blood* 116: 4532–4541.

MACS Product*	Order no.
Anti-ROR1-PE, human	130-098-317
Anti-ROR1-APC, human	130-098-320
Anti-ROR1-Biotin, human	130-098-312
Anti-ROR1 pure, human	130-098-243
Anti-ROR1 MicroBead Kit, human	130-103-929
MACSQuant Analyzer 10	130-096-343
autoMACS Pro Separator – Starter Kit	130-092-545

\* Products are for research use only.

Miltenyi Biotec's portfolio comprises thousands of antibodies. Select the antigen and conjugate you need at [www.miltenyibiotec.com/antibodies](http://www.miltenyibiotec.com/antibodies)

Miltenyi Biotec's Custom Antibody Design Service enables you to tailor single and multicolor antibody conjugates, as well as multicolor antibody panels that perfectly meet your specific research needs. Learn more about it at [www.miltenyibiotec.com/customantibodies](http://www.miltenyibiotec.com/customantibodies)

For more information on the principle of MACS MicroBead Technology, see [www.miltenyibiotec.com/microbead](http://www.miltenyibiotec.com/microbead)



Miltenyi Biotec provides products and services worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, the MACS logo, Vio, VioBlue, Vio700, and Vio770 are registered trademarks or trademarks of Miltenyi Biotec GmbH. Copyright © 2014 Miltenyi Biotec GmbH. All rights reserved.