



Miltenyi Biotec

FeraTrack™ Direct

MRI contrast agent for *in vitro* labeling

5 tracking experiments

130-104-185

25 tracking experiments

130-104-186

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
 - 1.3 Physico-chemical properties
 - 1.4 Quality control
 - 1.5 Requirements
2. Protocol
 - 2.1 General advices
 - 2.2 Preparation of cells
 - 2.3 Intracellular magnetic labeling
3. References

1. Description

This product is for research use only.

Components	<p>1×200 µL FeraTrack Direct Contrast Particles MRI cell tracking agent (superparamagnetic iron oxide [SPIO] nanoparticles) or 5×200 µL FeraTrack Direct Contrast Particles MRI cell tracking agent (superparamagnetic iron oxide [SPIO] nanoparticles)</p>
Capacity	<p>For 5×10⁶ total cells, up to 5 tracking experiments or for 25×10⁶ total cells, up to 25 tracking experiments.</p>
Product format	<p>FeraTrack Direct Contrast Particles are supplied in a sterile PBS/EDTA buffer solution with an iron concentration of 45 mM.</p>
Appearance	<p>FeraTrack Direct Contrast Particles are a black to reddish-brown aqueous solution of dextran-coated SPIO nanoparticles.</p>
Storage	<p>Store protected from light at 2–8 °C. The expiration date is indicated on the vial label.</p>

For laboratory and animal research use only. **Warning: Not for human or animal therapeutic or diagnostic use. Make sure to comply with all laws and regulations governing research on animals.**

1.1 Background information

Magnetic resonance imaging (MRI) is the most frequently used technique for serial *in vivo* cell tracking applications due to high resolution of soft tissues, which makes it especially useful for imaging of the brain, muscles, or the heart.¹ In order to improve detectability of transplanted cells and produce a strong contrast against surrounding tissue, intracellular labeling of cells with iron oxide particles before transplantation has been described.² FeraTrack Direct Contrast Particles are SPIO nanoparticles specifically formulated for direct uptake from cell culture medium without the need of transfection reagent. Uptake of FeraTrack Direct Contrast Particles was proven by Prussian Blue staining, MRI, and electron microscopy. Biocompatibility was assured by cell expansion and differentiation assays of intracellularly labeled cells. Once within cells, SPIOs can induce decreased signal intensity on T1, T2, and T2*-weighted images.³

1.2 Applications

FeraTrack Direct Contrast Particles are optimized for *in vitro* labeling of cell lines (e.g. NIH-3T3, Jurkat), primary cells (e.g. T cells), and stem cells (e.g. mesenchymal stem cells) independent of an additional transfection reagent.

1.3 Physico-chemical properties

Compound: dextran-coated ferumoxide (Fe₃O₄).
Mean particle size: 150 nm (hydrodynamic diameter, intensity weighted).
Dosage: 100 µg Fe / 1×10⁶ cells.

1.4 Quality control

FeraTrack Direct Contrast Particles are tested for absence of microbial contamination by Thioglycate and Caso-Bouillon sterility testing.

1.5 Requirements

- Polystyrene cell culture dish, e.g., 6-well plate
- Standard CO₂ incubator
- Cell culture medium
- Phosphate-buffered saline (PBS)

2. Protocol

2.1 General advices

- ▲ Read the entire protocol before starting.
- ▲ Ensure sterile handling of FeraTrack Direct Contrast Particles and cells.
- ▲ Removal of serum and antibiotics from medium is optional. FeraTrack Direct Contrast Particles can be taken up from cell culture media containing up to 10% serum.
- ▲ Standard animal handling procedures and local regulations must be followed in case of *in vivo* application of labeled cells.

2.2 Preparation of cells

1. Adherent cells: plate $0.5\text{--}1 \times 10^6$ cells per 9 cm^2 cell culture dish (6-well plate) in 2 mL culture medium one day before labeling to ensure 90% confluency at the time of labeling.
2. Suspension cells: plate $0.5\text{--}1 \times 10^6$ suspension cells per 9 cm^2 cell culture dish (6-well plate) in 2 mL culture medium at least 10 minutes prior to labeling.

2.3 Intracellular magnetic labeling

1. Add 40 μL of the FeraTrack Direct Contrast Particles dropwise to the cells.
2. Incubate cells at 37 °C for 4–6 hours in a CO₂ incubator.
3. After 4–6 hours wash cells twice with PBS.
4. Add cultivation medium to labeled cells for subsequent cultivation steps or harvest cells for alternative downstream applications.

▲ **Note:** Cell type-dependent, efficiency of intracellular labeling may be increased by adaption of FeraTrack Direct Contrast Particles volume and incubation time.

3. References

1. Rogers, W. J. *et al.* (2006) Technology inside: *in vivo* cell tracking by use of MRI. *Nat. Clin. Pract. Cardiovasc. Med.* 10: 554–562.
2. Hoehn, M. *et al.* (2007) Cell tracking using magnetic resonance imaging. *J. Physiol.* 584: 25–30.
3. Bulte, J. W. and Kraitchman, D. L. (2004) Iron oxide MR contrast agents for molecular and cellular imaging. *NMR Biomed.* 17: 484–499.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

FeraTrack, MACS, and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.