



Miltenyi Biotec

Find the fittest
MACS[®] ART Annexin V System

Improve the quality of spermatozoa preparations

Be specific, be successful

The presence of apoptotic cells in sperm preparations for assisted reproduction is known to decrease specimen quality. Unfortunately, commonly used techniques for preparing sperm may not remove apoptotic spermatozoa with the required specificity.

That's why the MACS® ART Annexin V System was designed to selectively remove apoptotic spermatozoa from sperm preparations. Even early apoptotic spermatozoa with an intact cellular membrane are efficiently depleted, to yield spermatozoa fractions for use in subsequent operations in Assisted Reproductive Technology (ART).

Select only the highest quality for your next sperm preparation.

Researchers have reported that purified sperm populations show higher overall quality in terms of:


- Morphology¹
- Sperm Deformity Index¹
- Motility^{2,3}
- DNA fragmentation
- Mitochondria membrane potential integrity^{3,4,5}
- Apoptosis
- Cryosurvival rates^{2,4}

Scientist obtained a sperm population with:

- Improved sperm fertilization potential (sperm capacitation index, hamster oocyte penetration assay)^{6,12}
- Improved Induced Acrosome Reaction Test (IART)^{7,8}
- Increased sperm chromatin decondensation^{9,11}
- Improved cleavage rates after fertilization¹⁰

References

- 1 Aziz, N. *et al.* (2007) *Human Reproduction* 22: 1413–1419.
- 2 Said, T. *et al.* (2005) *Fertility and Sterility* 83: 1442–1446.
- 3 De Vantéry Arrighi, C. *et al.* (2009) *Reproductive Biology Endocrinology* 7: 1.
- 4 Grunewald, S. *et al.* (2006) *Cell and Tissue Banking* 7: 99–104.
- 5 Paasch, U. *et al.* (2005) *Asian Journal of Andrology* 7: 61–69.
- 6 Said, T. *et al.* (2006) *Biology of Reproduction* 74: 530–537.
- 7 Lee, T. *et al.* (2010) *Human Reproduction* 25: 839–846.
- 8 Grunewald, S. *et al.* (2006) *Ann. N. Y. Acad. Sci.* 1090: 138–146.
- 9 Grunewald, S. *et al.* (2009) *Fertility and Sterility* 92: 572–577.
- 10 Dirican, E. K. *et al.* (2008) *Journal of Assisted Reproduction and Genetics* 25: 375–381.
- 11 Delbes, G. *et al.* (2013) *Journal of Andrology* 1: 689–706.
- 12 Sheikh, A. *et al.* (2013) *Journal of Andrology* 1: 845–849.



Specific depletion of apoptotic spermatozoa

The MACS ART Annexin V System is intended for *in vitro* depletion of apoptotic human spermatozoa. The system selectively depletes unwanted cells using the ability of Annexin V to recognize the antigen EPS (externalized phosphatidylserine) in the plasma membrane of apoptotic cells.

Phosphatidylserine-positive cells are depleted from fresh, cryopreserved, or otherwise manipulated semen samples. The yielded spermatozoa fraction is suitable for use in Assisted Reproductive Technology.

How it works

Apoptotic cell depletion begins with magnetic labeling of apoptotic cells by the MACS ART Annexin V Reagent. The labeled cells are then passed through a separation column located in a magnetic field. Unwanted cells are selectively retained in the column. Living spermatozoa are not labeled by the reagent, so they pass through the column and are collected for later use.




Figure 1: Translocation of PS from the intracellular to the extracellular surface of the plasma membrane during apoptosis.





Figure 2: Isolation of PS-exposing cells using the MACS ART Annexin V System.



MACS® ART Annexin V System

Depletion of apoptotic spermatozoa with the MACS ART Annexin V System is simple and fast. It adapts easily into laboratories to enhance classical sperm parameters and lower sperm apoptotic markers. Additionally, it can be combined with other sperm preparation techniques for procedural flexibility.

Components of the MACS ART Annexin V System include the MACS ART Annexin V Reagent, MACS ART Binding Buffer, and MACS ART MS Column. All are free of human and animal-derived components. Product specifications, such as endotoxin content or sterility, are confirmed by batch-specific certificates of analysis to ensure consistency of quality.

The system has been tested in a Mouse Embryo Assay (MEA), which is widely used to assess the toxicity of materials used in human assisted reproduction and to detect embryo-toxic materials. The MACS ART Annexin V System shows no embryo toxicity.




Figure 3: The depletion of apoptotic spermatozoa with the MACS ART Annexin V System is simple and fast.

Product	Content	Order no.
MACS ART Annexin V Reagent	1×0.1 mL reagent	200-070-502
MACS ART Annexin V Reagent	6×0.1 mL reagent	200-070-503
MACS ART Binding Buffer	6×10 mL buffer	200-070-505
MACS ART MS Column	6 columns	200-070-500
MACS ART Separation Unit	1 unit	200-070-501
MACS MultiStand	1 piece	130-042-303
MACS ART Annexin V Complete Kit	6×0.1 mL reagent, 6×10 mL buffer, 6 columns	200-070-507



Easy, fast, and specific depletion of unwanted apoptotic spermatozoa

- Improve quality of sperm preparations
- Obtain desired population within 30 minutes
- Proven technology for optimal results



Miltenyi Biotec

Germany/Austria

Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macsde@miltenyi.com

USA/Canada

Miltenyi Biotec Inc.
2303 Lindbergh Street
Auburn, CA 95602, USA
Phone 800 FOR MACS
Phone +1 530 888 8871
Fax +1 877 591 1060
macsau@miltenyi.com

Australia

Miltenyi Biotec
Australia Pty. Ltd.
Unit 11, 2 Eden Park Drive
Macquarie Park, NSW 2113
Australia
Phone +61 2 8877 7400
Fax +61 2 9889 5044
macsau@miltenyi.com

Benelux

Miltenyi Biotec B.V.
Sandfortdreef 17
2333 ZZ Leiden
The Netherlands
macsnl@miltenyi.com
**Customer service
The Netherlands**
Phone 0800 4020120
Fax 0800 4020100

Customer service Belgium

Phone 0800 94016
Fax 0800 99626

Customer service Luxembourg

Phone 800 24971
Fax 800 24984

China

Miltenyi Biotec Technology &
Trading (Shanghai) Co., Ltd.
Room 401
No. 1077, Zhangheng Road
Pudong New Area
201203 Shanghai, P.R. China
Phone +86 21 6235 1005
Fax +86 21 6235 0953
macscn@miltenyi.com

France

Miltenyi Biotec SAS
10 rue Mercoeur
75011 Paris, France
Phone +33 1 56 98 16 16
Fax +33 1 56 98 16 17
macsfr@miltenyi.com

Italy

Miltenyi Biotec S.r.l.
Via Paolo Nanni Costa, 30
40133 Bologna
Italy
Phone +39 051 6 460 411
Fax +39 051 6 460 499
macsit@miltenyi.com

Japan

Miltenyi Biotec K.K.
NEX-Eitai Building 5F
16-10 Fuyuki, Koto-ku
Tokyo 135-0041, Japan
Phone +81 3 5646 8910
Fax +81 3 5646 8911
macsjp@miltenyi.com

Nordics and Baltics

Miltenyi Biotec Norden AB
Medicon Village
Scheeleorget 1
223 81 Lund
Sweden
macsse@miltenyi.com

Customer service Sweden

Phone 0200 111 800
Phone 80 20 30 10

Customer service Denmark

Phone 80 20 30 10

**Customer service
Norway, Finland, Iceland,
and Baltic countries**

Phone +46 46 280 72 80

Singapore

Miltenyi Biotec Asia Pacific Pte Ltd
438B Alexandra Road, Block B
Alexandra Technopark
#06-01
Singapore 119968
Phone +65 6238 8183
Fax +65 6238 0302
macssg@miltenyi.com

South Korea

Miltenyi Biotec Korea Co., Ltd.
Arigi Bldg. 8F
562 Nonhyeon-ro
Gangnam-gu
Seoul 06136, South Korea
Phone +82 2 555 1988
Fax +82 2 555 8890
macskr@miltenyi.com

Spain

Miltenyi Biotec S.L.
C/Luis Buñuel 2
Ciudad de la Imagen
28223 Pozuelo de Alarcón (Madrid)
Spain
Phone +34 91 512 12 90
Phone +34 91 512 12 91
macses@miltenyi.com

Switzerland

Miltenyi Biotec Swiss AG
Gibelinstrasse 27
4500 Solothurn
Switzerland
Phone +41 32 623 08 47
Fax +49 2204 85197
macsch@miltenyi.com

United Kingdom

Miltenyi Biotec Ltd.
Almac House, Church Lane
Bisley, Surrey GU24 9DR, UK
Phone +44 1483 799 800
Fax +44 1483 799 811
macsuk@miltenyi.com

www.miltenyibiotec.com

Cover image: Copyright © Steve Young/stock.adobe.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

This print material is not for use in the USA. The MACS ART Annexin V System is manufactured by Miltenyi Biotec B.V. & Co. KG according to an ISO 13485 certified quality system. These products are available in Europe as CE-marked medical devices. For availability in your country please contact your local representative. Not available for use in the USA. MACS and the Miltenyi Biotec logo are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.