

Isolation and cultivation of oligodendrocytes from adult mouse or rat brain

Contents

1. Description
 - 1.1 Background information
 - 1.2 Reagent and instrument requirements
2. Protocol
 - 2.1 Preparation of brain dissociation
 - 2.2 Magnetic labeling
 - 2.3 Magnetic separation
 - 2.4 Flow cytometric analysis
 - 2.5 Cell culture
 - 2.6 Immunocytochemical staining of cultured oligodendrocytes

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

1. Description

This product is for research use only.

1.1 Background information

Oligodendrocytes are the myelinating cells of the central nervous system. The main function of oligodendrocytes is to provide support and insulation to axons by creating the myelin sheath around axons. That reduces ion leakage and increases impulse speed, as saltatory propagation of action potentials occurs at the nodes of Ranvier.

This protocol has been developed to generate highly purified and viable oligodendrocytes from adult mouse or rat brain tissue. Brain tissue from mice or rats older than P7 is dissociated into single-cell suspensions using the Adult Brain Dissociation Kit. The extracellular matrix is enzymatically digested using the kit components, while the gentleMACS™ Octo Dissociator with Heaters is used for the mechanical dissociation steps during the on-instrument enzyme incubation. After the dissociation, the myelin and cell debris are removed using the Debris Removal Solution and is followed by an subsequent removal of erythrocytes using the Red Blood Cell Removal Solution. The Anti-O4 MicroBeads are used to isolate oligodendrocytes from the single-cell suspension.

1.2 Reagent and instrument requirements

- Dulbecco's phosphate-buffered saline (D-PBS) with calcium, magnesium, glucose, and pyruvate. Keep buffer cold (2–8 °C).

- D-PBS/BSA buffer: Prepare a solution containing D-PBS and 0.5% bovine serum albumin (BSA) by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with D-PBS. Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column.

Always use freshly prepared buffer. Do **not use** autoMACS Running Buffer or MACSQuant® Running Buffer as they contain a small amount of sodium azide that could affect the results.

▲ **Note:** BSA can be replaced by other proteins such as mouse or rat serum albumin, mouse or rat serum, or fetal bovine serum (FBS).

- MACS Columns and MACS Separators: O4⁺ cells can be enriched by using MS Columns. Positive selection can also be performed by using the autoMACS Pro.

Column	Max. number of labeled cells	Max. number of total cells	Separator
MS	10 ⁷	2×10 ⁷	MiniMACS, OctoMACS, VarioMACS, SuperMACS II
autoMACS	5×10 ⁷	10 ⁸	autoMACS Pro

▲ **Note:** Column adapters are required to insert certain columns into the VarioMACS™ or SuperMACS™ II Separators. For details refer to the respective MACS Separator data sheet.

For preparation of brain dissociation

- Adult Brain Dissociation Kit, mouse and rat (# 130-107-677)
- gentleMACS™ Octo Dissociator with Heaters (# 130-096-427)
- gentleMACS C Tubes (# 130-093-237, # 130-096-334)
- 35 mm diameter sterile petri dish
- Sterile scalpel
- Sterile forceps
- (Optional) ART® 1000 REACH™ pipet tips (Molecular BioProducts, Inc.) for removal of dissociated material from the closed C Tubes.
- MACS SmartStrainers (70 μm) (# 130-098-462)
- 15 mL and 50 mL tubes
- Centrifuge with swinging bucket rotor

For cell isolation and flow cytometric analysis

- Anti-O4 MicroBeads, human, mouse, rat (# 130-094-543)
- (Optional) Pre-Separation Filters (70 μm) (# 130-095-823)
- Anti-O4 antibodies, human, mouse, rat (clone O4 or clone REA576) conjugated to PE or APC. For more information about antibodies refer to www.miltenyibiotec.com/antibodies.
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD for flow cytometric exclusion of dead cells.
- (Optional) MACSQuant® Analyzer 10 (# 130-096-343)

For cell culture and immunocytochemical staining

- Double-distilled water (ddH₂O)
- Imaging Plate CG 1.5 (24 well) (# 130-098-263)
- MACS Neuro Medium (# 130-093-570) and MACS NeuroBrew[®]-21 (# 130-093-566)
- 200 mM L-Glutamine
- Poly-L-Lysine (0.01%)
- Penicillin/streptomycin
- Human PDGF-AA (# 130-093-977)
- Human FGF-2 (# 130-093-837)
- 2% paraformaldehyde (PFA) for the fixation
- FcR Blocking Reagent
- Anti-O4 pure, human, mouse, rat (# 130-115-810) and anti-mouse-IgM secondary antibody
- Staining buffer: Prepare a solution containing autoMACS[®] Running Buffer (# 130-091-221) with FcR Blocking Reagent, mouse (# 130-092-575) in a ratio 1:10, e.g., add 1 mL FcR Blocking Reagent to 9 mL autoMACS Running Buffer.
- Phosphate-buffered saline (PBS)
- autoMACS Running Buffer (# 130-091-221)

2. Protocol

2.1 Preparation of brain dissociation

▲ For subsequent cell separation and cultivation it is recommended to dissociate at least 800 mg of adult mouse or rat brain tissue.

▲ Volumes given below are for one adult mouse brain (max. 500 mg) in 1980 µL enzyme mix. When working with less than 500 mg, use the same volumes as indicated. When working with higher tissue quantities or brain tissue from adult rats, determine the weight and scale up all reagent volumes and total volumes accordingly. A maximum of 500 mg brain tissue per C Tube can be processed.

▲ A swinging bucket rotor is recommended for centrifugation, e. g., Heraeus[®] Multifuge 4KR by Thermo Fisher[®] Scientific.

- Enzyme P is ready to use. Prepare aliquots of appropriate volume to avoid repeated freeze-thaw-cycles. Store aliquots at -20 °C. This solution is stable for 6 months. Resuspend the lyophilized powder in the vial labeled Enzyme A with 1 mL Buffer A. Do not vortex. This solution should then be aliquoted and stored at -20 °C for later use. Avoid repeated freeze-thaw-cycles.
- Prepare enzyme mix 1 and enzyme mix 2 according to the table below.

Enzyme mix 1		Enzyme mix 2	
Enzyme P	Buffer Z	Buffer Y	Enzyme A
50 µL	1900 µL	20 µL	10 µL

Preparation of 1× Red Blood Cell Removal Solution

- Dilute the Red Blood Cell Removal Solution (10×) 1:10 with double-distilled water (ddH₂O), for example, dilute 0.1 mL of cold Red Blood Cell Removal Solution (10×) with 0.9 mL cold ddH₂O.
▲ **Note:** Do not use deionized water for dilution!
- Store the prepared 1× Red Blood Cell Removal Solution at 2–8 °C. Discard unused solution at the end of the day.

Preparation of cell culture dish

- Prepare the following medium: MACS Neuro Medium containing 2% MACS Neuro Brew-21, 1% penicillin/streptomycin and 0.5 mM L-Glutamine, 10 ng/mL Human PDGF-AA, and 10 ng/mL Human FGF-2.
- Coat the culture dish (24-well plate) with 0.01% Poly-L-Lysine overnight at 37 °C and wash three times with ddH₂O afterwards. Let the culture dish dry under sterile conditions.

2.1.1 Dissociation protocol

▲ For details on the use of the gentleMACS[™] Octo Dissociator with Heaters, refer to the user manual.

▲ A maximum of 500 mg mouse or rat brain tissue in 2 mL enzyme mix can be processed in one C Tube.

▲ For dissociation of small amount of neural tissue (<100 mg), refer to data sheet of Adult Brain Dissociation Kit, mouse and rat.

▲ For cell culture experiments subsequent to tissue dissociation, all steps should be performed under sterile conditions.

1. Remove the mouse or rat brain. Wash the brain in cold D-PBS.
2. Prepare the appropriate volume of enzyme mix 1 (refer to table in chapter 2.1) and transfer it into a gentleMACS C Tube.
3. Place the brain on a petri dish and cut it into 8 sagittal slices using a scalpel.
4. Transfer the tissue pieces into the C Tube containing 1950 µL of enzyme mix 1.
5. Transfer 30 µL of enzyme mix 2 into the C Tube.
6. Tightly close C Tube and attach it upside down onto the sleeve of the gentleMACS Octo Dissociator with Heaters.
7. Run the gentleMACS Program **37C_ABDK_01**.
8. After termination of the program, detach C Tube from the gentleMACS Octo Dissociator with Heaters.
9. (Optional) Centrifuge briefly to collect the sample at the bottom of the tube.
10. Resuspend sample and apply it to a MACS SmartStrainer (70 µm) placed on a 50 mL tube.

▲ **Note:** Moisten MACS SmartStrainer with buffer before use.

▲ **Note:** When upscaling the reagent volume and total volumes, increase also the number of MACS SmartStrainers (70 µm). One MACS SmartStrainer (70 µm) can be used for one adult mouse brain.

▲ **Note:** Dissociated tissue can be removed from the closed C Tube by pipetting through the septum-sealed opening in the center of the cap of the C Tube. Use ART 1000 REACH 1000 µL pipette tips.

▲ **Note:** Cells with a diameter >70 µm may be lost. To obtain these cells within the flow through, use a cell strainer with an appropriate mesh size.

- Apply 10 mL of cold D-PBS onto the MACS SmartStrainer (70 μm).
- Discard MACS SmartStrainer (70 μm) and centrifuge cell suspension at 300 \times g for 10 minutes at 4 $^{\circ}\text{C}$. Aspirate supernatant completely.
- Proceed to 2.1.2 for debris and red blood cell removal.

2.1.2 Debris and red blood cell removal

▲ Volumes given below are for the cell suspension from up to two adult mouse brains (max. 1 g) or up to 1 g of rat brain as starting material. When working with higher tissue quantities, scale up all reagent volumes accordingly.

▲ A maximum of cell suspension from two adult mouse brains (max. 1 g) or up to 1 g of rat brain can be processed in one 15 mL reagent tube.

▲ Always use pre-cooled buffers and solutions (4 $^{\circ}\text{C}$).

	Debris Removal Solution	D-PBS	Overlay (D-PBS)
1 brain (400–500 mg)	900 μL	3100 μL	4 mL
2 brains (800–1000 mg)	1800 μL	6200 μL	4 mL

- Resuspend cell pellet carefully with the appropriate volume of cold D-PBS according to the table above and transfer cell suspension to a 15 mL tube. Do not vortex.
- Add appropriate volume of cold Debris Removal Solution.
- Mix well.
- Overlay very gently with 4 mL of cold D-PBS.

▲ **Note:** Pipette very slowly to ensure that the D-PBS phase overlays the cell suspension and phases are not mixed.
- Centrifuge at 4 $^{\circ}\text{C}$ and 3000 \times g for 10 minutes with full acceleration and full brake.

▲ **Note:** If centrifuges give suboptimal centrifugation, the acceleration and brake can be reduced.
- Three phases are formed. Aspirate the two top phases completely and discard them.
- Fill up with cold D-PBS to a final volume of 15 mL.
- Gently invert the tube three times. Do not vortex!
- Centrifuge at 4 $^{\circ}\text{C}$ and 1000 \times g for 10 minutes with full acceleration and full brake. Aspirate supernatant completely.
- Resuspend cell pellet from up to two adult mouse brains carefully in 1 mL of cold 1 \times Red Blood Cell Removal Solution. Do not vortex.
- Incubate for 10 minutes in the refrigerator (2–8 $^{\circ}\text{C}$).
- Add 10 mL of cold D-PBS/BSA buffer.
- Centrifuge at 4 $^{\circ}\text{C}$ and 300 \times g for 10 minutes. Aspirate supernatant completely.
- Proceed to 2.2 for magnetic labeling.



2.2 Magnetic labeling

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will

prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ Volumes for magnetic labeling given below are for up to 10^7 total cells. When working with fewer than 10^7 cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10^7 total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ For optimal performance it is important to obtain a single-cell suspension before magnetic labeling. Pass cells through 70 μm nylon mesh (MACS SmartStrainer (70 μm), # 130-098-462) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

▲ The recommended incubation temperature is 2–8 $^{\circ}\text{C}$. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

- Resuspend cell pellet in 97.5 μL of PB buffer per up to 10^7 total cells.
- Add 2.5 μL of Anti-O4 MicroBeads.
- Mix well and incubate for 15 minutes in the dark in the refrigerator (2–8 $^{\circ}\text{C}$).
- Wash cells by adding 1 mL of PB buffer and centrifuge at 300 \times g for 5 minutes. Aspirate supernatant completely.
- Resuspend up to 10^7 cells in 500 μL of PB buffer.

▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
- (Optional) Take 20 μL for later flow cytometric analysis (original fraction).
- Proceed to magnetic separation (2.3).



2.3 Magnetic separation

▲ Choose an MS Column and an appropriate MACS Separator. For details refer to the table in section 1.4.

▲ Always wait until the column reservoir is empty before proceeding to the next step.

▲ Degas buffer before use, as air bubbles could block the column.

Magnetic separation with MS Columns

- Place column in the magnetic field of a suitable MACS Separator. For details refer to the respective MACS Column data sheet.
- (Optional) Place Pre-Separation Filter (70 μm) on top of the column to remove clumps which may clog the column.

▲ **Note:** Moisten Pre-Separation Filter with buffer before use.
- Prepare column by rinsing with 500 μL of PB buffer.
- Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
- Wash column with 3 \times 500 μL of PB buffer. Collect unlabeled cells that pass through and combine with the flow-through from step 3. This is the O4⁻ cell fraction.

▲ **Note:** Perform washing steps by adding buffer aliquots as soon as the column reservoir is empty.

- Remove column from the separator and place it on a suitable collection tube.
- Pipette 1 mL of PB buffer onto the column. Immediately flush out the magnetically labeled cells by firmly pushing the plunger into the column. This is the O4⁺ cell fraction.
- To increase the purity of O4⁺ cells, it is recommended to enrich the positive fraction over a second MS Column. Repeat the magnetic separation procedure as described in steps 1 to 6 by using a new column.
- Proceed to flow cytometric analysis (2.4).

Magnetic separation with the autoMACS[®] Pro Separator

▲ Refer to the respective user manual for instructions on how to use the autoMACS[®] Pro Separator.

▲ Use D-PBS/BSA buffer. Buffers used for operating the autoMACS Pro Separator should have a temperature of ≥10 °C.

- Prepare and prime the instrument.
- Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
- For a standard separation choose the following program:
Positive selection: Posseld2
Collect positive fraction in row C of the tube rack. This fraction represents the O4⁺ cells.
- (Optional) Collect negative fraction in row B of the tube rack. This fraction represents the O4⁻ cells.
- Proceed to flow cytometric analysis (2.4).

2.4 Flow cytometric analysis

▲ The recommended antibody dilution for labeling of cells is 1:11 for up to 10⁷ cells/100 µL of D-PBS/BSA buffer.

▲ Volumes given below are for up to 10⁷ nucleated cells. When working with fewer than 10⁷ cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10⁷ nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).

- (Optional) For analysis take 100 µL of positive and negative fraction. Include the 20 µL of the original fraction.
- Resuspend up to 10⁷ nucleated cells per 100 µL of buffer.
- Add 10 µL of Anti-O4-APC.
- Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).
▲ Note: Higher temperatures and/or longer incubation times may lead to nonspecific cell labeling. Working on ice requires increased incubation times.
- Wash cells by adding 1 mL of buffer and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
- Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry, e.g. using the MACSQuant[®] Analyzer 10.

2.5 Cell culture

- Plate 10⁵ cells in 50 µL of medium as a drop in the middle of

each well of a 24-well plate which has been coated overnight (refer to 2.1 “Preparation of cell culture dish”).

- Let the cells settle down for 30 minutes at 37 °C in the incubator.
- Carefully add 450 µL of prepared medium to each well.
- Maintain the culture by replacement of 50% of prepared medium every other day.

▲ Note: Cells loosely attach after the first days of cultivation. Therefore, a cultivation period of approximately 1 week is recommended.

2.6 Immunocytochemical staining of cultured oligodendrocytes

- Wash cells 3× with PBS.
- Fix cells with 2% PFA for 10 minutes at room temperature.
- Wash cells 3× with PBS.
▲ Note: Fixed cells can be stored in azide-containing buffer at 2–8 °C for up to 1 week.
- Add staining buffer and incubate for 10 minutes at room temperature.
- Discard staining buffer.
- Add Anti-O4 pure antibody in staining buffer to the cells with a final concentration of 5–10 µg/mL and incubate at room temperature in the dark for 10 minutes.
- Wash cells 3× with autoMACS Running Buffer.
- Add a corresponding secondary antibody (anti-mouse IgM), in staining buffer to the cells and incubate at room temperature in the dark for 10 minutes.
- Wash cells 3× with autoMACS Running Buffer.
▲ Note: For co-staining with additional antibodies repeat steps 6–9.
- Store cells in autoMACS Running Buffer.
- Cells are now ready for immunofluorescence microscopy.
▲ Note: Samples can be stored at 2–8 °C in the dark for up to one week.
▲ Note: When working with cells cultured on coverslips, the coverslips need to be mounted onto slides before imaging.

Refer to www.miltenyibiotec.com for all data sheets and protocols.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

autoMACS, gentleMACS, MACS, MACSQuant, MidiMACS, MiniMACS, MultiMACS, NeuroBrew, OctoMACS, QuadroMACS, SuperMACS, and VarioMACS are either registered trademarks or trademarks of Miltenyi Biotec GmbH. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. Thermo Fisher is a registered trademark of Thermo Fisher Scientific Inc. or its subsidiaries.

ART is a registered trademark and REACH is a trademark of Molecular BioProducts, Inc.

Copyright © 2017 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.