

## Contents

1. Description
  - 1.1 Background information
  - 1.2 Reagent and instrument requirements
2. Protocol
  - 2.1 Preparation of brain dissociation
  - 2.2 Magnetic labeling
  - 2.3 Magnetic separation
  - 2.4 Flow cytometric analysis
  - 2.5 Cell culture

## Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

## 1. Description

**This product is for research use only.**

### 1.1 Background information

Neurons are the basic information processing structures in the central nervous system (CNS). The information is processed and transmitted via electrical and chemical signals. Neurons display a remarkable heterogeneity in their morphology, polarity, location, function, electrophysical characteristics, and in their neurotransmitter production.

This protocol has been developed to generate highly purified and viable neurons from adult mouse brain tissue. Brain tissue from mice older than P7 is dissociated into single-cell suspensions using the Adult Brain Dissociation Kit. The extracellular matrix is enzymatically digested using the kit components, while the gentleMACS™ Dissociator with Heaters is used for the mechanical dissociation steps during the on-instrument enzyme incubation. After the dissociation, the myelin and cell debris are removed using the Debris Removal Solution and is followed by an subsequent removal of erythrocytes using the Red Blood Cell Removal Solution. The Neuron Isolation Kit, mouse is used to isolate neurons from the single-cell suspension.

### 1.2 Reagent and instrument requirements

- Dulbecco's phosphate-buffered saline (D-PBS) with calcium, magnesium, glucose, and pyruvate. Keep buffer cold (2–8 °C).
- D-PBS/BSA buffer: Prepare a solution containing D-PBS and 0.5% bovine serum albumin (BSA) by diluting MACS® BSA Stock Solution (# 130-091-376) 1:20 with D-PBS. Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column. Always use freshly prepared buffer. Do **not use** autoMACS® Running Buffer or MACSQuant® Running Buffer as they contain a small amount of sodium azide that could affect the results.
  - ▲ **Note:** BSA can be replaced by other proteins such as mouse serum albumin, mouse serum, or fetal bovine serum (FBS).
- MACS Columns and MACS Separators: neurons can be enriched by depletion using LS Columns. Depletion can also be performed by using the autoMACS Pro Separator.

Column	Max. number of labeled cells	Max. number of total cells	Separator
LS	$2 \times 10^7$	$4 \times 10^7$	MidiMACS, QuadroMACS, VarioMACS, SuperMACS II
autoMACS	$5 \times 10^7$	$10^8$	autoMACS Pro

▲ **Note:** Column adapters are required to insert certain columns into the VarioMACS™ or SuperMACS™ II Separators. For details refer to the respective MACS Separator data sheet.

### For preparation of brain dissociation

- Adult Brain Dissociation Kit, mouse and rat (# 130-107-677)
- gentleMACS™ Octo Dissociator with Heaters (# 130-096-427)
- gentleMACS C Tubes (# 130-093-237, # 130-096-334)
- 35 mm diameter sterile petri dish
- Sterile scalpel
- Sterile forceps
- (Optional) ART® 1000 REACH™ pipet tips (Molecular BioProducts, Inc.) for removal of dissociated material from the closed C Tubes.
- MACS SmartStrainers (70 µm) (# 130-098-462)
- 15 mL and 50 mL tubes
- Centrifuge with swinging bucket rotor

### For cell isolation and flow cytometric analysis

- Neuron Isolation Kit, mouse (# 130-115-389, # 130-115-390)
- (Optional) Pre-Separation Filters (70 µm) (# 130-095-823)
- Fluorochrome-conjugated antibodies for flow cytometric analysis, e.g., an astrocyte-specific antibody as Anti-ACSA-2-PE, an oligodendrocyte-specific antibody as Anti-O4-PE, or an microglia-specific antibody as CD11b-FITC. For more

information about antibodies refer to [www.miltenyibiotec.com/antibodies](http://www.miltenyibiotec.com/antibodies).

- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD for flow cytometric exclusion of dead cells.
- (Optional) MACSQuant® Analyzer 10 (# 130-096-343)

#### For cell culture

- Double-distilled water (ddH<sub>2</sub>O)
- Imaging Plate CG 1.5 (24 well) (# 130-098-263)
- MACS Neuro Medium (# 130-093-570) and MACS NeuroBrew®-21 (# 130-093-566)
- L-Glutamine
- Poly-L-Lysine (0.01%)
- Penicillin/streptomycin
- Human BDNF

## 2. Protocol

### 2.1 Preparation of brain dissociation

▲ For subsequent cell separation and cultivation it is recommended to dissociate at least 800 mg of adult mouse brain tissue.

▲ Volumes given below are for one adult mouse brain (max. 500 mg) in 1980 µL enzyme mix. When working with less than 500 mg, use the same volumes as indicated. When working with higher tissue quantities scale up all reagent volumes and total volumes accordingly.

▲ A swinging bucket rotor is recommended for centrifugation, e. g., Heraeus® Multifuge 4KR by Thermo Fisher® Scientific.

- Enzyme P is ready to use. Prepare aliquots of appropriate volume to avoid repeated freeze-thaw-cycles. Store aliquots at -20 °C. This solution is stable for 6 months. Resuspend the lyophilized powder in the vial labeled Enzyme A with 1 mL Buffer A. Do not vortex. This solution should then be aliquoted and stored at -20 °C for later use. Avoid repeated freeze-thaw-cycles.
- Prepare enzyme mix 1 and enzyme mix 2 according to the table below.

Enzyme mix 1		Enzyme mix 2	
Enzyme P 50 µL	Buffer Z 1900 µL	Buffer Y 20 µL	Enzyme A 10 µL

#### Preparation of 1× Red Blood Cell Removal Solution

- Dilute the Red Blood Cell Removal Solution (10×) 1:10 with double-distilled water (ddH<sub>2</sub>O), for example, dilute 0.1 mL of cold Red Blood Cell Removal Solution (10×) with 0.9 mL cold ddH<sub>2</sub>O.

▲ Note: Do not use deionized water for dilution!

- Store the prepared 1× Red Blood Cell Removal Solution at 2–8 °C. Discard unused solution at the end of the day.

#### Preparation of cell culture dish

- Prepare the following medium: MACS Neuro Medium containing 2% MACS Neuro Brew-21, 1% Penicillin/Streptomycin and 0.5 mM L-Glutamine.

- Coat the culture dish (24-well plate) with 0.01% Poly-L-Lysine overnight at 37 °C and wash three times with ddH<sub>2</sub>O afterwards. Let the culture dish dry under sterile conditions.

### 2.1.1 Dissociation protocol

▲ For details on the use of the gentleMACS™ Octo Dissociator with Heaters, refer to the user manual.

▲ A maximum of one mouse brain (max. 500 mg) in 2 mL enzyme mix can be processed in one C Tube.

▲ For cell culture experiments subsequent to tissue dissociation, all steps should be performed under sterile conditions.

1. Remove the mouse brain. Wash the brain in cold D-PBS.
2. Prepare the appropriate volume of enzyme mix 1 (refer to table in chapter 2.1) and transfer it into a gentleMACS C Tube.
3. Place the brain on a petri dish and cut it into 8 sagittal slices using a scalpel.
4. Transfer the tissue pieces into the C Tube containing 1950 µL of enzyme mix 1.
5. Transfer 30 µL of enzyme mix 2 into the C Tube.
6. Tightly close C Tube and attach it upside down onto the sleeve of the gentleMACS Octo Dissociator with Heaters.
7. Run the gentleMACS Program **37C\_ABDK\_01**.
8. After termination of the program, detach C Tube from the gentleMACS Octo Dissociator with Heaters.
9. (Optional) Centrifuge briefly to collect the sample at the bottom of the tube.
10. Resuspend sample and apply it to a MACS SmartStrainer (70 µm) placed on a 50 mL tube.
  - ▲ Note: Moisten MACS SmartStrainer with buffer before use.
  - ▲ Note: When upscaling the reagent volume and total volumes, increase also the number of MACS SmartStrainers (70 µm). One MACS SmartStrainer (70 µm) can be used for one adult mouse brain.
  - ▲ Note: Dissociated tissue can be removed from the closed C Tube by pipetting through the septum-sealed opening in the center of the cap of the C Tube. Use ART 1000 REACH 1000 µL pipette tips.
  - ▲ Note: Cells with a diameter >70 µm may be lost. To obtain these cells within the flow through, use a cell strainer with an appropriate mesh size.
11. Apply 10 mL of cold (4 °C) D-PBS onto the MACS SmartStrainer (70 µm).
12. Discard MACS SmartStrainer (70 µm) and centrifuge cell suspension at 300×g for 10 minutes at 4 °C. Aspirate supernatant completely.
13. Proceed to 2.1.2 for debris and red blood cell removal.

### 2.1.2 Debris and red blood cell removal

▲ Volumes given below are for the cell suspension from up to two adult mouse brains. When working with higher tissue quantities, scale up all reagent volumes accordingly.

▲ A maximum of cell suspension from two adult mouse brains can be processed in one 15 mL reagent tube.

▲ Always use pre-cooled buffers and solutions (4 °C).

	Debris Removal Solution	D-PBS	Overlay (D-PBS)
1 brain (400–500 mg)	900 µL	3100 µL	4 mL
2 brains (800–1000 mg)	1800 µL	6200 µL	4 mL

1. Resuspend cell pellet carefully with the appropriate volume of cold D-PBS according to the table above and transfer cell suspension to a 15 mL tube. Do not vortex.
2. Add appropriate volume of cold Debris Removal Solution.
3. Mix well.
4. Overlay very gently with 4 mL of cold D-PBS.  
▲ **Note:** Pipette very slowly to ensure that the D-PBS phase overlays the cell suspension and phases are not mixed.
5. Centrifuge at 4 °C and 3000×g for 10 minutes with full acceleration and full brake.  
▲ **Note:** If centrifuges give suboptimal centrifugation, the acceleration and brake can be reduced.
6. Three phases are formed. Aspirate the two top phases completely and discard them.
7. Fill up with cold D-PBS to a final volume of 15 mL.
8. Gently invert the tube three times. Do not vortex!
9. Centrifuge at 4 °C and 1000×g for 10 minutes with full acceleration and full brake. Aspirate supernatant completely.
10. Resuspend cell pellet from up to two adult mouse brains carefully in 1 mL of cold 1× Red Blood Cell Removal Solution. Do not vortex.
11. Incubate for 10 minutes in the refrigerator (2–8 °C).
12. Add 10 mL of cold D-PBS/BSA buffer.
13. Centrifuge at 4 °C and 300×g for 10 minutes. Aspirate supernatant completely.
14. Proceed to 2.2 for magnetic labeling.



### 2.2 Magnetic labeling

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ Volumes for magnetic labeling given below are for up to 10<sup>7</sup> total cells. When working with fewer than 10<sup>7</sup> cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10<sup>7</sup> total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ For optimal performance it is important to obtain a single-cell suspension before magnetic labeling. Pass cells through 70 µm

nylon mesh (MACS SmartStrainer (70 µm), # 130-098-462) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

▲ The recommended incubation temperature is 2–8 °C. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

1. Resuspend cell pellet in 80 µL of D-PBS/BSA buffer per 10<sup>7</sup> total cells.
2. Add 20 µL of Non-Neuronal Cells Biotin-Antibody Cocktail.
3. Mix well and incubate for 5 minutes in the dark in the refrigerator (2–8 °C).
4. Wash cells by adding 1 mL of D-PBS/BSA buffer per 10<sup>7</sup> cells and centrifuge at 300×g for 5 minutes. Aspirate supernatant completely.
5. Resuspend cell pellet in 80 µL of D-PBS/BSA buffer per 10<sup>7</sup> cells.
6. Add 20 µL of Anti-Biotin MicroBeads.
7. Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).
8. Adjust volume to 500 µL per 10<sup>7</sup> cells with D-PBS/BSA buffer.  
▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
9. (Optional) Take 20 µL for later flow cytometric analysis (original fraction).
10. Proceed to magnetic separation (2.3).



### 2.3 Magnetic separation

▲ Choose an LS Column and an appropriate MACS Separator. For details refer to the table in section 1.4.

▲ Always wait until the column reservoir is empty before proceeding to the next step.

#### Magnetic separation with LS Columns

1. Place column in the magnetic field of a suitable MACS Separator. For details refer to the respective MACS Column data sheet.
2. (Optional) Place Pre-Separation Filter (70 µm) on top of the column to remove clumps which may clog the column.  
▲ **Note:** Moisten Pre-Separation Filter with buffer before use.
3. Prepare column by rinsing with 3 mL of D-PBS/BSA buffer.
4. Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
5. Wash column with 2×1 mL of D-PBS/BSA buffer. Collect unlabeled cells that pass through and combine with the flow-through from step 4. This is the target cell fraction (negative fraction).  
▲ **Note:** Perform washing steps by adding buffer aliquots only when the column reservoir is empty.
6. (Optional, if non-neuronal cells are needed) Remove column from the separator and place it on a suitable collection tube.

- Pipette 3 mL of D-PBS/BSA buffer onto the column. Immediately flush out the magnetically labeled cells by firmly pushing the plunger into the column. This fraction represents the non-target cells (positive fraction).
- Proceed to flow cytometric analysis (2.4).

#### Magnetic separation with the autoMACS® Pro Separator

▲ Refer to the respective user manual for instructions on how to use the autoMACS® Pro Separator.

▲ Use D-PBS/BSA buffer. Buffers used for operating the autoMACS Pro Separator should have a temperature of  $\geq 10$  °C.

- Prepare and prime the instrument.
- Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
- For a standard separation choose the following program:  
**Depletion: Depl05**  
Collect negative fraction in row B of the tube rack. This fraction represents the target cells.
- (Optional if non-neuronal cells are needed) Collect positive fraction in row C of the tube rack. This fraction represents the non-target cells.
- Proceed to flow cytometric analysis (2.4).

#### 2.4 Flow cytometric analysis

▲ The recommended antibody dilution for labeling of cells is 1:10 for up to  $10^6$  cells/50  $\mu$ L of D-PBS/BSA buffer.

▲ Volumes given below are for up to  $10^6$  nucleated cells. When working with fewer than  $10^6$  cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for  $2 \times 10^6$  nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).

- (Optional) For analysis take 100  $\mu$ L of positive and negative fraction. Include the 20  $\mu$ L of the original fraction.
- Centrifuge at  $300 \times g$  for 5 minutes. Aspirate supernatant completely.
- Resuspend up to  $10^6$  nucleated cells per 45  $\mu$ L of D-PBS/BSA buffer.  
▲ **Note:** If more staining antibodies than Anti-Biotin-PE shall be added to the sample adjust buffer volume, accordingly.
- Add 5  $\mu$ L of Anti-Biotin-PE.
- (Optional) Add antibodies specific for non-neuronal cells, for example, 5  $\mu$ L of Anti-ACSA-2 and/or 5  $\mu$ L of Anti-O4 and/or 5  $\mu$ L of CD11b antibody.
- Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).  
▲ **Note:** Higher temperatures and/or longer incubation times may lead to nonspecific cell labeling. Working on ice requires increased incubation times.
- Wash cells by adding 1 mL of D-PBS/BSA buffer and centrifuge at  $300 \times g$  for 10 minutes. Aspirate supernatant completely.
- Resuspend cell pellet in a suitable amount of D-PBS/BSA buffer for analysis by flow cytometry, e.g. using the MACSQuant® Analyzer 10.

#### 2.5 Cell culture

- Plate  $10^5$  cells in 50  $\mu$ L of prepared medium as a drop in the middle of each well of a 24-well plate which has been coated overnight (refer to 2.1 “Preparation of cell culture dish”).
- Let the cells settle down for 30 minutes at 37 °C in the incubator.
- Carefully add 450  $\mu$ L of prepared medium to each well.
- Take off the medium to remove non-attached and dead cells.
- (Optional) Treat neurons with 50 ng/mL Human BDNF in 500  $\mu$ L of prepared medium for 3 to 6 hours in the incubator. Remove BDNF-containing medium afterwards.
- Add 500  $\mu$ L of prepared medium.
- Maintain the culture by replacement of 50% of medium every other day.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols.

#### Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

autoMACS, gentleMACS, MACS, MACSQuant, MidiMACS, MiniMACS, MultiMACS, NeuroBrew, OctoMACS, QuadroMACS, SuperMACS, and VarioMACS are either registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. Thermo Fisher is a registered trademark of Thermo Fisher Scientific Inc. or its subsidiaries.

ART is a registered trademark and REACH is a trademark of Molecular BioProducts, Inc.

Copyright © 2018 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.