

## Contents

1. Description
  - 1.1 Background information
  - 1.2 Reagent and instrument requirements
2. Protocol
  - 2.1 Preparation of brain dissociation
  - 2.2 Magnetic labeling
  - 2.3 Magnetic separation
  - 2.4 Flow cytometric analysis
  - 2.5 Cell culture
  - 2.6 Immunocytochemical staining of cultured microglia

## Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

## 1. Description

This product is for research use only.

### 1.1 Background information

Microglia, originating from the embryonic yolk sac and populating the brain during early embryonic development, are the resident immune-effector cells in the central nervous system (CNS). The most characteristic feature of microglial cells is their rapid activation in response to injury, inflammation, neurodegeneration, infection, and brain tumors. In the mature CNS, microglia represent about 10–15 percent of all the cells in the brain. Activated microglia serve as the major antigen-presenting cells in the CNS. They are morphologically, immunophenotypically, and functionally related to cells of the monocyte/macrophage lineage.

This protocol has been developed to generate highly purified and viable microglia from adult mouse or rat brain tissue. Brain tissue from mice or rats older than P7 is dissociated into single-cell suspensions using the Adult Brain Dissociation Kit. The extracellular matrix is enzymatically digested using the kit components, while the gentleMACS™ Dissociator with Heaters is used for the mechanical dissociation steps during the on-instrument enzyme incubation. After the dissociation, the myelin and cell debris are removed using the Debris Removal Solution and is followed by subsequent removal of erythrocytes using the Red Blood Cell Removal Solution. The CD11b (Microglia) MicroBeads, mouse and human or the CD11b/c (Microglia) MicroBeads, rat are used to isolate microglia from the single-cell suspension.

### 1.2 Reagent and instrument requirements

- Dulbecco's phosphate-buffered saline (D-PBS) with calcium, magnesium, glucose, and pyruvate. Keep buffer cold (2–8 °C).
- Phosphate-buffered saline (PBS)
- PB buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2 and 0.5% bovine serum albumin (BSA) by diluting MACS® BSA Stock Solution (# 130-091-376) 1:20 with PBS. Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column. Always use freshly prepared buffer. Do **not use** autoMACS® Running Buffer or MACSQuant® Running Buffer as they contain a small amount of sodium azide that could affect the results.
- ▲ **Note:** BSA can be replaced by other proteins such as mouse or rat serum albumin, mouse or rat serum, or fetal bovine serum (FBS).
- MACS Columns and MACS Separators: microglia can be enriched by using MS or LS Columns. Positive selection can also be performed by using the autoMACS Pro Separator.

Column	Max. number of labeled cells	Max. number of total cells	Separator
MS	10 <sup>7</sup>	2×10 <sup>7</sup>	MiniMACS, OctoMACS, VarioMACS, SuperMACS II
LS	2×10 <sup>7</sup>	4×10 <sup>7</sup>	MidiMACS, QuadroMACS, VarioMACS, SuperMACS II
autoMACS	5×10 <sup>7</sup>	10 <sup>8</sup>	autoMACS Pro

▲ **Note:** Column adapters are required to insert certain columns into the VarioMACS™ or SuperMACS™ II Separators. For details refer to the respective MACS Separator data sheet.

### For preparation of brain dissociation

- Adult Brain Dissociation Kit, mouse and rat (# 130-107-677)
- gentleMACS Octo Dissociator with Heaters (# 130-096-427)
- gentleMACS C Tubes (# 130-093-237, # 130-096-334)
- 35 mm diameter sterile petri dish
- Sterile scalpel
- Sterile forceps
- (Optional) ART® 1000 REACH™ pipet tips (Molecular BioProducts, Inc.) for removal of dissociated material from the closed C Tubes.
- MACS SmartStrainers (70 μm) (# 130-098-462)
- 15 mL and 50 mL tubes
- Centrifuge with swinging bucket rotor

## For cell isolation and flow cytometric analysis

- CD11b (Microglia) MicroBeads, mouse (# 130-093-634, # 130-093-636)
- CD11b/c (Microglia) MicroBeads, rat (# 130-105-634, # 130-105-543)
- (Optional) Pre-Separation Filters (70 µm) (# 130-095-823)
- CD11b antibodies, human and mouse (clone M1/70.15.11.5), or CD45 antibodies, mouse (clone 30F11), or CD11b/c antibodies rat (clone REA325), or CD45 antibodies, rat (clone REA504) conjugated to PE or APC. For more information about antibodies refer to [www.miltenyibiotec.com/antibodies](http://www.miltenyibiotec.com/antibodies).
- (Optional) Propidium Iodide Solution (# 130-093-233) or 7-AAD Staining Solution (# 130-111-568) for flow cytometric exclusion of dead cells.
- (Optional) MACSQuant<sup>®</sup> Analyzer 10 (# 130-096-343)

## For cell culture and immunocytochemical staining

- Double-distilled water (ddH<sub>2</sub>O)
- Imaging Plate CG 1.5 (24 well) (# 130-098-263)
- DMEM with stable glutamine
- 200 mM L-glutamine
- Poly-L-lysine (0.01%)
- Penicillin/streptomycin
- FBS (fetal bovine serum)
- 2% paraformaldehyde (PFA) for the fixation
- CD11b pure, human and mouse (# 130-115-811) and anti-rat IgG2bk secondary antibody or CD68 pure, mouse (# 130-115-808) and anti-rat IgG2a secondary antibody
- FcR Blocking Reagent, mouse (# 130-092-575)
- autoMACS<sup>®</sup> Running Buffer (# 130-091-221)
- Staining buffer: Prepare a solution containing autoMACS Running Buffer (# 130-091-221) with FcR Blocking Reagent, mouse (# 130-092-575) in a ratio 1:10, e.g., add 1 mL FcR Blocking Reagent to 9 mL autoMACS Running Buffer.
- Phosphate-buffered saline (PBS)
- (Optional) 0.2% TRITON<sup>™</sup> X-100 in PBS

## 2. Protocol

### 2.1 Preparation of brain dissociation

▲ For subsequent cell separation and cultivation it is recommended to dissociate at least 800 mg of adult mouse or rat brain tissue.

▲ Volumes given below are for one adult mouse brain (max. 500 mg) in 1980 µL enzyme mix. When working with less than 500 mg, use the same volumes as indicated. When working with higher tissue quantities, scale up all reagent volumes and total volumes accordingly. A maximum of 500 mg brain tissue per C Tube can be processed.

▲ A swinging bucket rotor is recommended for centrifugation, e. g., Heraeus<sup>®</sup> Multifuge 4KR by Thermo Fisher<sup>®</sup> Scientific.

- Enzyme P is ready to use. Prepare aliquots of appropriate volume to avoid repeated freeze-thaw-cycles. Store

aliquots at -20 °C. This solution is stable for 6 months. Resuspend the lyophilized powder in the vial labeled Enzyme A with 1 mL Buffer A. Do not vortex. This solution should then be aliquoted and stored at -20 °C for later use. Avoid repeated freeze-thaw-cycles.

- Prepare enzyme mix 1 and enzyme mix 2 according to the table below.

Enzyme mix 1		Enzyme mix 2	
Enzyme P 50 µL	Buffer Z 1900 µL	Buffer Y 20 µL	Enzyme A 10 µL

### Preparation of 1× Red Blood Cell Removal Solution

- Dilute the Red Blood Cell Removal Solution (10×) 1:10 with double-distilled water (ddH<sub>2</sub>O), for example, dilute 0.1 mL of cold Red Blood Cell Removal Solution (10×) with 0.9 mL cold ddH<sub>2</sub>O.

▲ **Note:** Do not use deionized water for dilution!

- Store the prepared 1× Red Blood Cell Removal Solution at 2–8 °C. Discard unused solution at the end of the day.

### Preparation of cell culture dish

- Prepare the following medium: DMEM containing 10% FBS, 1% penicillin/streptomycin, and 2 mM L-glutamine.
- Coat the culture dish (24-well plate) with 0.01% Poly-L-lysine overnight at 37 °C and wash three times with ddH<sub>2</sub>O afterwards.

#### 2.1.1 Dissociation protocol

▲ For details on the use of the gentleMACS<sup>™</sup> Octo Dissociator with Heaters, refer to the user manual.

▲ A maximum of 500 mg mouse or rat brain tissue in 2 mL enzyme mix can be processed in one C Tube.

▲ For cell culture experiments subsequent to tissue dissociation, all steps should be performed under sterile conditions.

1. Remove the mouse or rat brain. Wash the brain in cold D-PBS.
2. Prepare the appropriate volume of enzyme mix 1 (refer to table in chapter 2.1) and transfer it into a gentleMACS C Tube.
3. Place the brain on a petri dish and cut it into 8 sagittal slices using a scalpel.
4. Transfer the tissue pieces into the C Tube containing 1950 µL of enzyme mix 1.
5. Transfer 30 µL of enzyme mix 2 into the C Tube.
6. Tightly close C Tube and attach it upside down onto the sleeve of the gentleMACS Octo Dissociator with Heaters.
7. Run the gentleMACS Program **37C\_ABDK\_01**.
8. After termination of the program, detach C Tube from the gentleMACS Octo Dissociator with Heaters.
9. (Optional) Centrifuge briefly to collect the sample at the bottom of the tube.
10. Resuspend sample and apply it to a MACS<sup>®</sup> SmartStrainer (70 µm) placed on a 50 mL tube.

▲ **Note:** Moisten MACS SmartStrainer with buffer before use.

▲ **Note:** When upscaling the reagent volume and total volumes, increase also the number of MACS SmartStrainers (70 µm). One MACS SmartStrainer (70 µm) can be used for one adult mouse brain.

▲ **Note:** Dissociated tissue can be removed from the closed C Tube by pipetting through the septum-sealed opening in the center of the cap of the C Tube. Use ART 1000 REACH 1000 µL pipette tips.

▲ **Note:** Cells with a diameter >70 µm may be lost. To obtain these cells within the flow through, use a cell strainer with an appropriate mesh size.

11. Apply 10 mL of cold (4 °C) D-PBS onto the MACS® SmartStrainer (70 µm).
12. Discard MACS SmartStrainer (70 µm) and centrifuge cell suspension at 300×g for 10 minutes at 4 °C. Aspirate supernatant completely.
13. Proceed to 2.1.2 for debris and red blood cell removal.

### 2.1.2 Debris and red blood cell removal

▲ Volumes given below are for the cell suspension from up to two adult mouse brains or up to 1 g of rat brain as starting material. When working with higher tissue quantities, scale up all reagent volumes accordingly.

▲ A maximum of cell suspension from two adult mouse brains or up to 1 g of rat brain can be processed in one 15 mL reagent tube.

▲ Always use pre-cooled buffers and solutions (4 °C).

	Debris Removal Solution	D-PBS	Overlay (D-PBS)
1 brain (400–500 mg)	900 µL	3100 µL	4 mL
2 brains (800–1000 mg)	1800 µL	6200 µL	4 mL

▲ **Note:** In case of very small amount of tissue (< 100 mg) cell debris removal can be performed in a 5 mL reagent tube using 450 µL of Debris Removal Solution, 1550 µL of D-PBS for resuspension of the cell pellet, and 2000 µL D-PBS for overlay.

1. Resuspend cell pellet carefully with the appropriate volume of cold D-PBS according to the table above and transfer cell suspension to a 15 mL tube. Do not vortex.
2. Add appropriate volume of cold Debris Removal Solution.
3. Mix well.
4. Overlay very gently with 4 mL of cold D-PBS.  
▲ **Note:** Pipette very slowly to ensure that the D-PBS phase overlays the cell suspension and phases are not mixed.
5. Centrifuge at 4 °C and 3000×g for 10 minutes with full acceleration and full brake.  
▲ **Note:** If centrifuges give suboptimal centrifugation, the acceleration and brake can be reduced.
6. Three phases are formed. Aspirate the two top phases completely and discard them.
7. Fill up with cold D-PBS to a final volume of 15 mL.
8. Gently invert the tube three times. Do not vortex!
9. Centrifuge at 4 °C and 1000×g for 10 minutes with full acceleration and full brake. Aspirate supernatant completely.
10. Resuspend cell pellet from up to two adult mouse brains carefully in 1 mL of cold 1× Red Blood Cell Removal Solution. Do not vortex.
11. Incubate for 10 minutes in the refrigerator (2–8 °C).
12. Add 10 mL of cold PB buffer.
13. Centrifuge at 4 °C and 300×g for 10 minutes. Aspirate supernatant completely.

14. Proceed to 2.2 for magnetic labeling.



## 2.2 Magnetic labeling

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ Volumes for magnetic labeling given below are for up to 10<sup>7</sup> total cells. When working with fewer than 10<sup>7</sup> cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10<sup>7</sup> total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ For optimal performance it is important to obtain a single-cell suspension before magnetic labeling. Pass cells through 70 µm nylon mesh (MACS SmartStrainer (70 µm), # 130-098-462) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

▲ The recommended incubation temperature is 2–8 °C. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

### Mouse microglia cells

1. Resuspend cell pellet carefully in 90 µL of cold PB buffer per 10<sup>7</sup> total cells by pipetting slowly up and down. Do not vortex.
2. Add 10 µL of CD11b (Microglia) MicroBeads, human and mouse.
3. Mix well and incubate for 15 minutes in the dark in the refrigerator (2–8 °C).
4. Wash cells by adding 1 mL of cold PB buffer per 10<sup>7</sup> cells and centrifuge at 300×g for 5 minutes. Aspirate supernatant completely.
5. Resuspend up to 10<sup>7</sup> cells in 500 µL of PB buffer.  
▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
6. (Optional) Take 20 µL for later flow cytometric analysis (original fraction).
7. Proceed to magnetic separation (2.3).

### Rat microglia cells

1. Resuspend cell pellet carefully in 80 µL of cold PB buffer per 10<sup>7</sup> total cells by pipetting slowly up and down. Do not vortex.
2. Add 20 µL of CD11b/c (Microglia) MicroBeads, rat.
3. Mix well. Do not vortex. Incubate for 15 minutes in the dark in the refrigerator (2–8 °C).
4. Wash cells by adding 1 mL of cold PB buffer per 10<sup>7</sup> cells and centrifuge at 300×g for 5 minutes. Aspirate supernatant completely.
5. Resuspend up to 10<sup>7</sup> cells in 500 µL of PB buffer.  
▲ **Note:** For higher cell numbers, scale up buffer volume accordingly.
6. (Optional) Take 20 µL for later flow cytometric analysis (original fraction).
7. Proceed to magnetic separation (2.3).



## 2.3 Magnetic separation

- ▲ Choose an MS or LS Column and an appropriate MACS® Separator. For details refer to the table in section 1.2.
- ▲ Always wait until the column reservoir is empty before proceeding to the next step.

### Magnetic separation with MS or LS Columns

1. Place column in the magnetic field of a suitable MACS Separator. For details refer to the respective MACS Column data sheet.
2. (Optional) Place Pre-Separation Filter (70 µm) on top of the column to remove clumps which may clog the column.

▲ **Note:** Moisten Pre-Separation Filter with buffer before use.

3. Prepare column by rinsing with the appropriate amount of PB buffer:

MS: 500 µL      LS: 3 mL

4. Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
5. Wash column with the appropriate amount of buffer. Collect unlabeled cells that pass through and combine with the flow-through from step 3. This is the non-target cell fraction.

MS: 3×500 µL      LS: 3×3 mL

▲ **Note:** Perform washing steps by adding buffer aliquots only when the column reservoir is empty.

6. Remove column from the separator and place it on a suitable collection tube.
7. Pipette the appropriate amount of buffer onto the column. Immediately flush out the magnetically labeled cells by firmly pushing the plunger into the column. This is the target cell fraction (positive fraction).

MS: 1 mL      LS: 5 mL

8. To increase the purity of microglia cells, it is recommended to enrich the positive fraction over a second MS or LS Column. Repeat the magnetic separation procedure as described in steps 1 to 6 by using a new column.
9. Proceed to flow cytometric analysis (2.4).

### Magnetic separation with the autoMACS® Pro Separator

- ▲ Refer to the respective user manual for instructions on how to use the autoMACS® Pro Separator.
- ▲ Use PBS buffer. Buffers used for operating the autoMACS Pro Separator should have a temperature of ≥10 °C.

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.
3. For a standard separation choose the following program:

#### Positive selection of mouse microglia cells: Possel

Collect positive fraction in row C of the tube rack. This fraction represents the target cells.

#### Positive selection of rat microglia cells: Posseld2

Collect positive fraction in row C of the tube rack. This fraction represents the target cells.

4. (Optional) Collect negative fraction in row B of the tube rack. This fraction represents the non-target cells.
5. Proceed to flow cytometric analysis (2.4).

## 2.4 Flow cytometric analysis

### Mouse microglia cells

▲ The recommended antibody dilution for labeling of cells is 1:10 for up to 10<sup>7</sup> cells/100 µL of buffer.

▲ Volumes given below are for up to 10<sup>7</sup> nucleated cells. When working with fewer than 10<sup>7</sup> cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10<sup>7</sup> nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).

1. (Optional) For analysis take 100 µL of positive and negative fraction. Include the 20 µL of the original fraction.
2. Resuspend up to 10<sup>7</sup> nucleated cells per 100 µL of buffer.
3. Add 10 µL of CD11b-PE or CD11b-APC.
4. Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).

▲ **Note:** Higher temperatures and/or longer incubation times may lead to nonspecific cell labeling. Working on ice requires increased incubation times.

5. Wash cells by adding 1 mL of buffer and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
6. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry, e.g. using the MACSQuant® Analyzer 10.

### Rat microglia cells

▲ The recommended antibody dilution for labeling of cells and subsequent analysis by flow cytometry is 1:10 for up to 10<sup>6</sup> cells/50 µL of buffer.

▲ Volumes given below are for up to 10<sup>6</sup> nucleated cells. When working with fewer than 10<sup>6</sup> cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10<sup>6</sup> nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).

1. (Optional) For analysis take 100 µL of positive and negative fraction. Include the 20 µL of the original fraction.
2. Resuspend up to 10<sup>6</sup> nucleated cells per 45 µL of buffer.
3. Add 5 µL of CD11b/c-PE or CD11b/c-APC.
4. Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).

▲ **Note:** Higher temperatures and/or longer incubation times may lead to nonspecific cell labeling. Working on ice requires increased incubation times.

5. Wash cells by adding 1 mL of buffer and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
6. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry, e.g. using the MACSQuant Analyzer 10.

## 2.5 Cell culture

1. Plate  $10^5$  cells in 50  $\mu$ L of prepared medium as a drop in the middle of each well of a 24-well plate which has been coated overnight (refer to 2.1 “Preparation of cell culture dish”).
4. Let the cells settle down for 30 minutes at 37 °C in the incubator.
5. Carefully add 450  $\mu$ L of prepared medium to each well.
6. Maintain the culture by replacement of 50% of prepared medium every other day.  
**▲ Note:** Cells loosely attach after the first days of cultivation. Therefore, a cultivation period of approximately 1 week is recommended.

## 2.6 Immunocytochemical staining of cultured microglia

1. Wash cells 3× with PBS.
2. Fix cells with 2% PFA for 10 minutes at room temperature.
3. Wash cells 3× with PBS.  
**▲ Note:** When working with CD68 antibodies, add 0.2% TRITON™ X-100 in PBS, incubate for 10 minutes at room temperature, and wash cells 3× with autoMACS® Running Buffer.  
**▲ Note:** Fixed cells can be stored in azide-containing buffer at 2–8 °C for up to 1 week.
4. Add staining buffer and incubate for 10 minutes at room temperature.
5. Discard staining buffer.
6. Add primary antibody in staining buffer to the cells with a final concentration of 5–10  $\mu$ g/mL.  
When working with CD11b antibodies: Incubate at room temperature in the dark for 10 minutes.  
When working with CD68 antibodies: Incubate at 2–8 °C overnight in the dark.
7. Wash cells 3× with autoMACS Running Buffer.
8. Add a corresponding secondary antibody in staining buffer to the cells.  
When working with CD11b antibodies: Incubate at room temperature in the dark for 10 minutes.  
When working with CD68 antibodies: Incubate at room temperature in the dark for 1 hour.
9. Wash cells 3× with autoMACS Running Buffer.  
**▲ Note:** For co-staining with additional antibodies repeat steps 6–9.
10. Store cells in autoMACS Running Buffer.
11. Cells are now ready for immunofluorescence microscopy.  
**▲ Note:** Samples can be stored at 2–8 °C in the dark for up to one week.  
**▲ Note:** When working with cells cultured on coverslips, the coverslips need to be mounted onto slides before imaging.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols.

## Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

autoMACS, gentleMACS, MACS, the MACS logo, MACSQuant, MidiMACS, MiniMACS, OctoMACS, QuadroMACS, SuperMACS, and VarioMACS are either registered trademarks or trademarks of Miltenyi Biotec GmbH and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

HERAEUS is a registered trademark of Heraeus Holding GmbH licensed to Thermo Fisher Scientific. Thermo Fisher is a registered trademark of Thermo Fisher Scientific Inc. or its subsidiaries.

ART is a registered trademark and REACH is a trademark of Molecular BioProducts, Inc.

Copyright © 2017 Miltenyi Biotec GmbH and/or its affiliates. All rights reserved.