



Multi Tissue Dissociation Kit 2

Preservation of cell surface epitopes

Background

The Multi Tissue Dissociation Kit 2 enables the thorough dissociation of various tough tissues to enhance the yield of cells found at low frequencies after dissociation, such as endothelial cells. Preservation of cell surface epitopes is essential for several downstream applications such as cell separation, flow analysis, and sorting. The dissociation process involves enzymatic digestion steps. Here we tested how the enzymes of the Multi Tissue Dissociation Kit 2 affect the integrity of key epitopes serving for the identification of major leukocyte populations, endothelial cells, and epithelial cells in mouse and human. Our results show that many epitopes are well preserved during enzymatic treatment.

Method

For mouse tissue, preservation of cell surface epitopes was tested on mouse spleens, lungs, blood, or perfused livers, either treated with the Multi Tissue Dissociation Kit 2 enzymes or incubated in buffer without enzymes (30 min at 37 °C). For human tissue, preservation of cell surface epitopes was tested on PBMCs, monocyte-derived dendritic cells, HUVECs, or SKBR3 cells either treated with the Multi Tissue Dissociation Kit 2 enzymes or incubated in buffer without enzymes (30 min at 37 °C).

The enzymatically treated and untreated samples were stained with each of the fluorochrome-conjugated MACS® Antibodies detecting cell surface epitopes and subsequently analyzed by flow cytometry. A decrease in fluorescence signal intensity of the fluorochrome conjugates in the treated cells indicated effects on epitope integrity.

Results, mouse

Antibody detecting cell surface epitope	Antibody clone	Antibody detecting cell surface epitope	Antibody clone	Antibody detecting cell surface epitope	Antibody clone
Anti-F4/80	REA126	CD16/32	REA377	CD106	429 (MCCAM.A)
Anti-Gr-1	RB6-8C5	CD19	6D5	CD107a	1D4B
Anti-LGR5	DA04-10E8.9	CD25	REA568	CD127	A7R34
Anti-Ly-6C	1G7.G10	CD31	390	CD140a	APA5
Anti-Ly-6G	REA526	CD40	FGK45.5	CD140b	APB5
Anti-MHC II	M5/114.15.2	CD43	REA364	CD205	NLDC-145
Anti-NK1.1	PK136	CD45	30F11.1	CD309	Avas12
Anti-NKp46	29A1.4.9	CD45R	RA3-6B2	CD326 (EpcAM)	caa7-9G8
Anti-PDCA-1	JF05-1C2.4.1	CD49b	REA541		
Anti-Siglec-H	551.3D3	CD51	mEF-SK4		
Anti-Ter119	Ter-119	CD61	2C9.G2(HMβ3-1)		
CD3e	REA606	CD62L	MEL-14-H2.100		
CD4	REA604	CD80	16-10A1		
CD8a	REA601	CD81	EAT2		
CD11b	REA592	CD90.2	30-H12		
CD11c	N418	CD105	MJ7/18		

Results, human

Antibody detecting cell surface epitope	Antibody clone	Antibody detecting cell surface epitope	Antibody clone	Antibody detecting cell surface epitope	Antibody clone
anti-CLA	HECA-452	CD28	15E8	CD133	AC133
anti-HLA-DR	AC122	CD31	AC128	CD138	B-B4
anti-KIR2D	NKVFS1	CD33	AC104.3E3	CD141 (BDCA-3)	AD5-14H12
anti-LGR5	DA03-22H2.8	CD34	AC136	CD144	REA199
anti-NKp80	4A4.D10	CD36	AC106	CD146	ME-9F1
anti-PTK7 (CCK-4)	188B	CD38	IB6	CD158a/h (KIR2DL1/DS1)	11PB6
anti-Slan (M-DC8)	DD1	CD39	MZ18-23C8	CD158b (KIR2DL2/DL3)	DX27
anti-TCR α/β	BW242/412	CD40	HB14	CD158e (KIR3DL1)	DX9
anti-TCR γ/δ	11F2	CD44	DB105	CD158e/k (KIR3DL1/DL2)	5.331
CD1c (BDCA-1)	AD5-8E7	CD45	5B1	CD161 (NKR-P1A)	191B8
CD2	LT2	CD45RA	T6D11	CD163	GHI/61.1
CD3	BW264/56	CD45RO	UCHL1	CD197 (CCR7)	FR11-11E8
CD4	VIT4	CD49d	MZ18-24A9	CD203c	FR3-16A11
CD4	M-T466	CD49e	NKI-SAM1	CD209	DCN47.5
CD6	M-T411	CD49f	GoH3	CD226 (DNAM-1)	DX11
CD8	BW135/80	CD56	AF12-7H3	CD266	ITEM-4
CD10	97C5	CD57	TB03	CD271	ME20.4 1.H4
CD11b	M1/70.15.11.5	CD62L	145/15	CD294 (CRTH2)	BM16
CD11c	MJ4-27G12.4.6	CD63	H5C6	CD303 (BDCA-2)	AC144
CD14	TÜK4	CD66abce	TET2	CD314 (NKG2D)	BAT221
CD15	VIMC6	CD73	AD2	CD324	67A4
CD16	VEP13	CD84	MZ18-21F6	CD326 (EpCAM)	HEA-125
CD19	LT19	CD90	DG3	CD335 (NKp46)	9E2
CD20	LT20	CD95	DX2	CD337 (NKp30)	AF29-4D12
CD24	32D12	CD101	BB27		
CD25	4E3	CD105	43A4E1		
CD26	FR10-11G9	CD123	AC145		
CD27	M-T271	CD127	MB15-18C9		

Legend
 ■ Stable
 ■ Moderate sensitivity*
 ■ Strong sensitivity*

* Reducing the concentration of one of the kit components (Enzyme P) can aid in preserving sensitive epitopes. For more information please contact us.



Miltenyi Biotec GmbH | Phone +49 2204 8306-0 | Fax +49 2204 85197 | macs@miltenyibiotec.de | www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS is a registered trademark of Miltenyi Biotec GmbH. Copyright © 2018 Miltenyi Biotec GmbH. All rights reserved.