

### Contents

1. Description
  - 1.1 Principle of a suppression assay using the MSC Suppression Inspector
  - 1.2 Background information
  - 1.3 Applications
  - 1.4 Reagent and instrument requirements
2. Protocol
  - 2.1 Sample preparation
    - 2.1.1 Preparation of cells
    - 2.1.2 Preparation of MSC Suppression Inspector
  - 2.2 Stimulation and suppression assay
3. Examples of a CFSE- and a tritium-based suppression assay using the MSC Suppression Inspector
4. References

### 1. Description

**Components** 2.5 mL MSC Suppression Inspector, human:  $5 \times 10^7$  Anti-Biotin MACSiBead™ Particles pre-loaded with biotinylated CD2, CD3, and CD28 antibodies.


**Product format** MSC Suppression Inspector is supplied in an azide-free suspension.

**Storage** Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

#### 1.1 Principle of a suppression assay using the MSC Suppression Inspector

Mesenchymal stem cells (MSCs) are often functionally analyzed *in vitro* by a so-called suppression assay. For this purpose, MSCs are co-cultured with  $CD4^+CD25^-$  or  $CD4^+$  responder T cells (Tresp) at different ratios in the presence of a polyclonal stimulus, in this case the MSC Suppression Inspector. Tresp cells alone show a proliferative response. Co-culture of MSCs with Tresp cells results in reduced proliferation of Tresp cells. Cell proliferation is determined by  $^3H$ -thymidine incorporation but can also be detected by carboxyfluorescein succinimidyl ester (CFSE) staining. The suppression assay is performed with a dilution series ranging from a ratio of 1:1 to 8:1 of Tresp cells:MSCs as outlined in tables 1 and 3. As additional control, Tresp and MSCs are cultured alone with and without the MSC Suppression Inspector. The dilution series is carried out in triplicate to achieve significant results. All volumes given in the protocol are calculated for one assay.

### Principle of the assay



## 1.2 Background information

MSCs are fibroblast-like plastic-adherent cells that can be isolated from a variety of tissues, such as bone marrow or adipose tissue. During the last few years the attention of scientists was redirected away from the multipotentiality of MSCs towards their possibility for immunomodulation. It was observed that bone marrow derived MSCs suppress T-cell proliferation.<sup>1,2</sup> This function of MSCs can be analyzed using the MSC Suppression Inspector which contains an optimal T cell stimulation reagent for a MSC suppression assay. The MSC Suppression Inspector consists of Anti-Biotin MACSiBead Particles that are pre-loaded with biotinylated CD2, CD3, and CD28 antibodies.

## 1.3 Applications

- Functional characterization of human MSCs by *in vitro* suppression assays.

## 1.4 Reagent and instrument requirements

- (Optional) CD4<sup>+</sup>CD25<sup>+</sup> Regulatory T Cell Isolation Kit, human (# 130-091-301), CD4<sup>+</sup>CD25<sup>+</sup>CD127<sup>dim/-</sup> Regulatory T Cell Isolation Kit II, human (# 130-094-775).
- Cell culture medium, for example, RPMI 1640 (# 130-091-440) supplemented with 10% AB serum, X-VIVO 15™ (Cambrex), or X-VIVO 15™ supplemented with 5% AB serum.
  - ▲ **Note:** 2-Mercaptoethanol (0.01 mM) can be added to preserve cell viability in case of rapid cell growth.
- 96-well culture plates (flat bottom).
- <sup>3</sup>H-thymidine.
- CFSE: 5(6)-carboxyfluorescein diacetate *N*-succinimidyl ester.
- Humidified incubator.
- (Optional) NH Expansion Medium, human (# 130-091-680).
- (Optional) CytoMix – MSC, human (# 130-093-552).
- (Optional) CD271 MicroBead Kit (PE) (# 130-092-819) or CD271 MicroBead Kit (APC) (# 130-092-283).
- (Optional) Anti-MSCA-1 (W8B2) MicroBead Kit, human (# 130-093-583)

## 2. Protocol

▲ This protocol has been developed for a tritium-based suppression assay. For a special protocol for a CFSE-based suppression assay refer to [www.macs-stemcells.com/downloads](http://www.macs-stemcells.com/downloads).

### 2.1 Sample preparation

▲ All steps in the protocol have to be performed under aseptic conditions. In this protocol one MACSiBead Particle per cell (bead-to-cell ratio 1 : 1) is used for stimulation.

Ratio Tresp cells : MSCs	Tresp cells	MSCs	MSC Suppression Inspector (amount of MACSiBead Particles)
1:0	5 × 10 <sup>4</sup>	–	5 × 10 <sup>4</sup>
0:1	–	5 × 10 <sup>4</sup>	5 × 10 <sup>4</sup>
1:1	5 × 10 <sup>4</sup>	5 × 10 <sup>4</sup>	10 × 10 <sup>4</sup>
2:1	5 × 10 <sup>4</sup>	2.5 × 10 <sup>4</sup>	7.5 × 10 <sup>4</sup>
4:1	5 × 10 <sup>4</sup>	1.3 × 10 <sup>4</sup>	6.3 × 10 <sup>4</sup>
8:1	5 × 10 <sup>4</sup>	0.6 × 10 <sup>4</sup>	5.6 × 10 <sup>4</sup>
Control 1:0	5 × 10 <sup>4</sup>	–	–
Control 0:1	–	5 × 10 <sup>4</sup>	–
Total cells/ MACSiBeads	3 × 10 <sup>5</sup>	2 × 10 <sup>5</sup>	4 × 10 <sup>5</sup>
Total cells/ MACSiBeads for 1 assay (triplicates)	9 × 10 <sup>5</sup>	6 × 10 <sup>5</sup>	12 × 10 <sup>5</sup>

**Table 1:** Number of responder T cells (Tresp), mesenchymal stem cells (MSCs), and MSC Suppression Inspector (MACSiBead Particles) per well.

### 2.1.1 Preparation of cells

#### Determine the concentration and the total number of responder T cells (Tresp)

▲ Start with CD4<sup>+</sup>CD25<sup>-</sup> or CD4<sup>+</sup> responder T cells isolated under aseptic conditions, e.g., with the CD4<sup>+</sup>CD25<sup>+</sup> Regulatory T Cell Isolation Kit, human (# 130-091-301). For details concerning Treg isolation refer to the respective data sheet.

Human PBMCs	
<b>Depletion of non-CD4<sup>+</sup> cells</b>	<ol style="list-style-type: none"> <li>1. Indirect magnetic labeling of non-CD4<sup>+</sup> cells with CD4<sup>+</sup> T Cell Biotin-Antibody Cocktail and Anti-Biotin MicroBeads.</li> <li>2. Magnetic separation using an LD Column or an autoMACS Column (program "Depl05").</li> </ol>
Pre-enriched CD4 <sup>+</sup> cells (flow-through fraction)	
<b>Depletion of CD4<sup>+</sup>CD25<sup>+</sup> regulatory T cells (Tregs)</b>	<ol style="list-style-type: none"> <li>1. Direct magnetic labeling of CD25<sup>+</sup> T cells with CD25 MicroBeads.</li> <li>2. Magnetic separation using two MS Columns or an autoMACS Column (program "Posseld2").</li> </ol>
CD4 <sup>+</sup> CD25 <sup>-</sup> Tresp (flow-through fraction, first column) CD4 <sup>+</sup> CD25 <sup>+</sup> Tregs (eluted cells, second column)	

**Table 2:** Isolation of Tresp with the CD4<sup>+</sup>CD25<sup>+</sup> Regulatory T Cell Isolation Kit, human.

### Preparation of mesenchymal stem cells (MSCs)

Use MSCs after isolation from human tissue, e.g., bone marrow or adipose tissue.

When using frozen MSCs it is suggested to use the cells after 2–3 days of cultivation.

1. Determine the concentration and the total number of MSCs and Tresp cells. For one assay, as outlined in table 1,  $9 \times 10^5$  Tresp cells and  $6 \times 10^5$  MSCs are needed.
2. Transfer required volumes of cell suspension to suitable tubes.
3. Add 5–10 volumes culture medium to the cells and centrifuge at  $300 \times g$  for 10 minutes. Aspirate supernatant completely.
4. Resuspend the Tresp cells ( $9 \times 10^5$ ) in 1800  $\mu\text{L}$  of medium and the MSCs ( $6 \times 10^5$ ) in 1200  $\mu\text{L}$ . The concentration of the cell suspensions is now  $5 \times 10^5$  cells/mL.
5. Pipette the appropriate volumes of MSCs and Tresp cell suspension in a 96-well culture plate. Refer to table 3 for the respective volumes.

#### 2.1.2 Preparation of MSC Suppression Inspector

1. Resuspend MSC Suppression Inspector thoroughly and transfer 60  $\mu\text{L}$  to a suitable tube.

▲ Note: Concentration of MSC Suppression Inspector is  $2 \times 10^7$  MACSiBead Particles per mL.

2. Add 0.3–0.6 mL of culture medium and centrifuge at  $300 \times g$  for 5 minutes. Aspirate supernatant completely.
3. Resuspend MSC Suppression Inspector in 120  $\mu\text{L}$  of culture medium. The reagent is now ready to use.

▲ Note: Concentration of prepared MSC Suppression Inspector is  $1 \times 10^7$  MACSiBead Particles per mL.

Ratio Tresp cells: MSCs	Tresp cells ( $5 \times 10^5$ cells/mL)	MSCs ( $5 \times 10^5$ cells/mL)	MSC Suppression Inspector ( $1 \times 10^7$ MACSiBead particles/mL)	Culture medium
1:0	100 $\mu\text{L}$	–	5 $\mu\text{L}$	105 $\mu\text{L}$
0:1	–	100 $\mu\text{L}$	5 $\mu\text{L}$	105 $\mu\text{L}$
1:1	100 $\mu\text{L}$	100 $\mu\text{L}$	10 $\mu\text{L}$	–
2:1	100 $\mu\text{L}$	50 $\mu\text{L}$	7.5 $\mu\text{L}$	53 $\mu\text{L}$
4:1	100 $\mu\text{L}$	25 $\mu\text{L}$	6.5 $\mu\text{L}$	79 $\mu\text{L}$
8:1	100 $\mu\text{L}$	12.5 $\mu\text{L}$	6.0 $\mu\text{L}$	92 $\mu\text{L}$
Control 1:0	100 $\mu\text{L}$	–	–	110 $\mu\text{L}$
Control 0:1	–	100 $\mu\text{L}$	–	110 $\mu\text{L}$
Total volume	600 $\mu\text{L}$	387.5 $\mu\text{L}$	40 $\mu\text{L}$	654 $\mu\text{L}$
Total volume for 1 assay (triplicates)	1800 $\mu\text{L}$	1200 $\mu\text{L}$	120 $\mu\text{L}$	approx. 2 mL

**Table 3:** Pipetting scheme for one assay with a total volume of 210  $\mu\text{L}$  per well using cell suspensions that contain  $5 \times 10^5$  cells/mL.


### 2.2 Stimulation and suppression assay

1. Resuspend the prepared MSC Suppression Inspector thoroughly and add the required amount to the wells (bead-to-cell ratio 1:1). For a detailed pipetting scheme refer to table 3.  
▲ Note: The bead-to-cell ratio refers to the total cell number per well.
2. Fill up wells to a total volume of 210  $\mu\text{L}$  with culture medium (refer to table 3).
3. Incubate at 37 °C and 5–7%  $\text{CO}_2$  for 4–5 days.
4. Add 1  $\mu\text{Ci}$   $^3\text{H}$ -thymidine to each well and incubate at 37 °C and 5–7%  $\text{CO}_2$  for 16 hours.
5. Measure  $^3\text{H}$ -thymidine incorporation, e.g., by using a liquid scintillation counter.

### 3. Examples of a CFSE- and a tritium-based suppression assay using the MSC Suppression Inspector

#### Tritium-based

MSCs were isolated from human bone marrow and culture expanded with NH Expansion Medium (# 130-091-680). After two passages MSCs were co-cultured with  $\text{CD4}^+\text{CD25}^-$  responder T cells at different ratios. For T cell stimulation, the MSC Suppression Inspector was added to the culture. As controls, MSCs and  $\text{CD4}^+\text{CD25}^-$  responder T cells alone were cultured without any stimulus. Proliferation of T cells was determined by  $^3\text{H}$ -thymidine incorporation.  $^3\text{H}$ -thymidine was added for 16 hours after 5 days of culture.




Tresp show high proliferation after stimulation with the MSC Suppression Inspector (A). When adding MSCs Tresp proliferation is suppressed dramatically (B). Unstimulated Tresp show no proliferation (D). MSCs alone show little proliferation with or without stimulation (C and E).

### CFSE-based

MSCs were isolated from human bone marrow either by plastic adherence (PA-MSCs) or by CD271 isolation (CD271-MSCs) using CD271 MicroBead Kit (APC) (# 130-092-283). Both PA-MSCs and CD271-MSC were culture expanded with NH Expansion Medium (# 130-091-680).

After two passages MSCs were co-cultured with CFSE-labeled CD4<sup>+</sup>CD25<sup>-</sup> responder T cells. For T cell stimulation, the MSC Suppression Inspector was added to the cultures. As control, MSCs and CD4<sup>+</sup>CD25<sup>-</sup> Tresp were cultured without the MSC Suppression Inspector. Cells were harvested after 5 days and the percentage of proliferating Tresp was measured as CFSE dye dilution analyzed by flow cytometry using the MACSQuant<sup>®</sup> Analyzer.



Almost 90% of all Tresp proliferate after stimulation with the MSC Suppression Inspector (A). When adding PA-MSCs Tresp proliferation is suppressed to a level of about 28% (B). The addition of CD271<sup>+</sup> MSCs suppresses the Tresp proliferation to a level of 20% (D). MSCs alone show no proliferation with or without stimulation (C, E, G, and H). Unstimulated Tresp show a proliferation rate of 30% (F).

### 4. References

1. Di Nicola, M. *et al.* (2002) Human bone marrow stromal cells suppress T-lymphocyte proliferation induced by cellular or nonspecific mitogenic stimuli. *Blood* 99: 3838–3843.
2. Bartholomew, A. *et al.* (2002) Mesenchymal stem cells suppress lymphocyte proliferation *in vitro* and prolong skin graft survival *in vivo*. *Exp. Hematol.* 30: 42–48.

All protocols and data sheets are available at [www.miltenyibiotec.com](http://www.miltenyibiotec.com).

### Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS, the MACS logo, and MACSQuant are registered trademarks and MACSiBead is a trademark of Miltenyi Biotec GmbH.

X-VIVO 15 is a registered trademark of Cambrex or its subsidiaries.

Copyright © 2011 Miltenyi Biotec GmbH. All rights reserved.