

Excerpt from MACS&more Vol 15 – 2/2013


Enrichment, stimulation, and viral transduction of naive and central memory CD8⁺ T cells under GMP conditions for translational research towards the development of adoptive cell therapy of cancer patients

Anna Casati¹, Azam Varghaei-Nahvi², Steven A. Feldman², Mario Assenmacher¹, Steven A. Rosenberg², Mark E. Dudley², and Alexander Scheffold³

¹ Miltenyi Biotec GmbH, Bergisch Gladbach, Germany

² Surgery Branch, National Cancer Institute, National Institute of Health, Bethesda, MD, USA

³ Department of Rheumatology and Clinical Immunology, Charité, Universitätsmedizin Berlin, Berlin, Germany CCS Istituto Neurologico Carlo Besta, Milan, Italy



Miltenyi Biotec

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. CliniMACS, gentleMACS, MACS, MACSiBead, MACSiMAG, the MACS logo, and TransAct are registered trademarks or trademarks of Miltenyi Biotec GmbH. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only. Copyright © 2013 Miltenyi Biotec GmbH. All rights reserved.

Enrichment, stimulation, and viral transduction of naive and central memory CD8⁺ T cells under GMP conditions for translational research towards the development of adoptive cell therapy of cancer patients

Anna Casati¹, Azam Varghaei-Nahvi², Steven A. Feldman², Mario Assenmacher¹, Steven A. Rosenberg², Mark E. Dudley², and Alexander Scheffold³

¹ Miltenyi Biotec GmbH, Bergisch Gladbach, Germany

² Surgery Branch, National Cancer Institute, National Institute of Health, Bethesda, MD, USA

³ Department of Rheumatology and Clinical Immunology, Charité, Universitätsmedizin Berlin, Berlin, Germany

Introduction

Adoptive immunotherapy involving the transfer of genetically engineered T cells can mediate tumor regression in patients with metastatic melanoma or other forms of cancer.¹⁻⁴ However, the efficacy of T cell immunotherapy depends on a number of critical steps. Autologous T cells collected from the patients are activated, genetically modified, expanded, and finally infused back into the patient to unfold their anti-tumor activity.

Peripheral blood T cells are heterogeneous with regard to their function and differentiation states, i.e., their self-renewal capacity and persistence. In particular terminally differentiated effector memory T (TEM) cells show poor persistence and only weak capacity to mediate tumor regression¹⁻⁵. It has been shown that engineering of T cells can improve their anti-tumor reactivity and *in vivo* persistence.^{6,7} Moreover, naive T cells (TN) show a higher efficacy in adoptive cell transfer in a cancer setting than memory T cells.⁸⁻¹⁰ Among memory T cells, central memory T cells (TCM) were shown to exert a higher anti-tumor activity than TEM cells in mice^{11,12}, non-human primate¹³, and humanized mouse models¹⁴. However, the translation of this principle into therapeutic strategies is hampered by the lack of appropriate protocols for the clinical-scale isolation of CD8⁺ TN and TCM cells, as well as GMP-compliant reagents and reliable processes

for cell activation and transduction.

One of the current standard procedures for cell activation involves soluble CD3 antibodies in combination with, for example, autologous antigen-presenting cells (APCs) for cross-linking. This protocol is highly dependent on the quality of the cross-linking cells and requires large amounts of viral supernatants, which is a major cost factor in this process. The other standard T cell activation process involves CD3 and CD28 antibodies immobilized on solid surfaces, such as cell culture vessels or large cell-sized particles.^{8,13} These methods have major drawbacks, as they are time-consuming and difficult to standardize for GMP compliance. Moreover, large particles need to be removed prior to infusion of T cells into the patient.

We developed a reliable method for the enrichment of CD8⁺ TN and TCM cells under GMP conditions. We also introduced a new flexible nanomatrix conjugated with CD3 and CD28 antibodies (TransAct™ Reagent), which allows easier handling for the activation of T cells in a closed system under GMP conditions. The activated TN and TCM cells isolated from blood of melanoma patients were expanded and transduced efficiently. TN and TCM cells activated with TransAct Reagent or CD3/CD28-conjugated cell-sized MACSiBead™ Particles showed only minimal effector differentiation.¹⁵

Materials and methods

Enrichment of human CD8⁺ TN and TCM cells

For all experiments fresh leukapheresis samples from patients with metastatic melanoma were used. Patients were enrolled in clinical trials approved by the National Cancer Institute Institutional Review Board. PBMCs were prepared by density gradient centrifugation (Ficoll™ Hypaque™) and washed twice in PBS. For each enrichment, 2–4×10⁹ PBMCs were used.

TN cells were enriched with a three-step procedure. Firstly, PBMCs were incubated with CliniMACS® CD4 MicroBeads on a rotator for 30 min at room temperature (RT). Cells were washed by centrifuging and resuspending in 100 mL of buffer A, i.e., CliniMACS PBS/EDTA Buffer supplemented with 0.5% human albumin (Plasbumin®-25 25%, USP; Talecris Biotherapeutics), and transferred to a collection bag. Magnetically labeled CD4⁺ cells were depleted using the CliniMACS Plus Instrument (program Depletion 2.1) and the CliniMACS Tubing Set LS. Secondly, the unlabeled target cells were incubated with CliniMACS CD62L MicroBeads for 15 min in cold buffer A and washed. Target cells were enriched using the CliniMACS Plus Instrument (program Enrichment 3.2) and the CliniMACS Tubing Set. Finally, the target fraction was labeled with a CD45RA-Biotin, human antibody and Anti-

Biotin MACSiBead Particles. Cells were further enriched using the MACSiMAG™ Separator.¹⁵ TCM cells were enriched with a two-step procedure. PBMCs were incubated with CliniMACS CD4 MicroBeads and CliniMACS CD45RA MicroBeads (in cold buffer A) on a rotator for 30 min at RT. Subsequently, the cells were washed, resuspended in 100 mL of buffer A, and transferred to a collection bag. Magnetically labeled cells were depleted using the CliniMACS Plus Instrument (program Depletion 2.1) and the CliniMACS Tubing Set LS. The unlabeled target cells were incubated with CliniMACS CD62L MicroBeads for 15 min in cold buffer A and washed. Target cells were enriched using the CliniMACS Plus Instrument (program Enrichment 3.2) and the CliniMACS Tubing Set.¹⁵

Flexible nanomatrix for T cell stimulation (TransAct Reagent)

The TransAct Reagent was provided by Miltenyi Biotec. TransAct Reagent is a colloidal reagent consisting of nanoscale iron oxide crystals embedded into a biocompatible polysaccharide matrix with an overall diameter of about 100 nm. Antibodies against CD3 (clone OKT3) and CD28 (clone 15E8) are covalently attached to the matrix. The matrix can be produced under cGMP conditions, sterilized by filtration and unbound reagent can easily be removed from the cell suspension by centrifugation of the cells. In the experiments described here, the reagent was used at 0.5 µg/mL of both antibodies conjugated in a 1:1 ratio. Anti CD3 and anti CD28 TransAct Reagents can also be used when conjugated on separate matrices and mixed in variable ratios and concentrations to achieve optimal T cell activation.

T cell stimulation and transduction

Three different methods were used for stimulation of the enriched T cell subsets: 1) CD3 antibody coated on the cell culture dish (OKT3; 10 µg/mL) in combination with soluble CD28 antibody (15E8; 1 µg/mL), referred to as CD3/CD28, 2) MACSiBead Particles loaded with CD3, CD28, and CD2 antibodies, and 3) TransAct Reagent. Interleukin 2 (IL-2; 300 IU/mL) was present in the cell culture for all three approaches.¹⁵ T cells were transduced using the mF5 retrovirus encoding an anti-MART-1 T cell receptor (TCR) comprising constant regions of the mouse TCR α and

β chains in addition to human TCR CDR antigen-combining regions¹⁶. Transductions were performed as described.¹⁷

Flow cytometry

To analyze the phenotypes of the T cell subsets, cells were labeled with fluorochrome-conjugated antibodies against CD8, CD4, CD45RA, CD62L, CD127, IFN- γ , IL-2, TNF- α (all from Miltenyi Biotec), CCR7, CD27, CD45RO (BD® Biosciences), and MART-1 tetramer (Beckman-Coulter; Allele HLA-A 0201, peptide MART-1, ELAGIGILTV). For the analysis of intracellular cytokine production, cells were cocultured with the 624 (MART-1⁺HLA-A2⁺) melanoma cell line in the presence of brefeldin A. After 6 hours cells were labeled with the CD8 antibody, fixed, and permeabilized with the Inside Stain Kit (Miltenyi Biotec). Subsequently, cells were analyzed by flow cytometry.

Statistical analysis

Groups were compared using the two-tailed unpaired t-test. The error bars in the graphs represent the standard error of means.

Results

Clinical-scale enrichment of human CD8⁺ T_N and T_{CM} cells

Human CD8⁺ T cells represent a heterogeneous population of T_N, T_{CM}, and T_{EM} cells, which can be distinguished by their expression of CD45RA and CD62L¹⁸ (table 1). We used these phenotypic characteristics to isolate T_N and T_{CM} cells. T_N cells were enriched in a three-step procedure, involving CD4⁺ cell depletion, positive selection of CD62L⁺ cells, and a second positive selection of CD45RA⁺ cells. Purity and recovery are shown in table 1. The remaining CD8⁺ cells contained in the target fraction were mainly NK and B cells.

T_{CM} cells were enriched in a two-step procedure, i.e., simultaneous depletion of CD4⁺ and CD45RA⁺ cells and subsequent

positive selection of CD62L⁺ cells. The unlabeled non-target fraction from the CD62L⁺ cell enrichment step was enriched in CD8⁺CD45RA⁻CD62L⁻ cells and represented the T_{EM} cells in our *in vitro* experiments.

Phenotypic characterization of CD8⁺ T_N and T_{EM} cells

The CD8⁺ cell subsets are characterized by differential expression of a number of additional markers. As CD62L and CCR7 are generally co-expressed in T_N and T_{CM} cells, we analyzed the enriched cell populations for CCR7 surface expression (fig. 1) and found that 70±5.5% of the T_N cells and 54±7% of the T_{CM} cells expressed CCR7. In contrast, only 13.2±2.4% of the T_{EM} cell fraction expressed CCR7 (see fig. 1 for representative dot plots). A large proportion of the T_N cell fraction (62.5±5.7%) expressed CD127 (IL-7R), whereas expression in T_{CM} and T_{EM} cells amounted to only 30±3.7% and 30.7±7.3%, respectively. Expression of CD27 is known to positively correlate with replicative capacity and tumor response^{19,20}. We found that CD27 expression was significantly higher in T_N and T_{CM} cells than in T_{EM} cells (88±1.8% and 80.5±5.3 vs. 62±5.4%). In contrast, granzyme A expression was higher in T_{EM} cells (85±2.2%) compared to T_N and T_{CM} cells (16.7±3.6% and 42.5±6.6%). This is an agreement with the notion that lymph node homing capacity and the acquisition of effector function are coordinately regulated^{18,21}. Taken together, these data indicate that the GMP-compliant procedures used for enrichment of T_N and T_{CM} cells yield populations with distinct, typical phenotypes (fig. 1).

Transduction of T cell subsets with antigen-specific TCR

The standard protocols for *in vitro* activation of T cells involve autologous antigen-presenting cells (APCs), which are required for cross-linking of soluble CD3 antibodies

T cell subset	Phenotype	Purity	Recovery
Naive T (T _N) cells	CD8 ⁺ CD45RA ⁺ CD62L ⁺	75.2±15.5%	44.6±16.7%
Central memory T (T _{CM}) cells	CD8 ⁺ CD45RA ⁻ CD62L ⁺	84±9.6%	60.9±23.8%
Effector memory T (T _{EM}) cells	CD8 ⁺ CD45RA ⁻ CD62L ⁻	n.d.	n.d.

Table 1 Phenotypes of human CD8⁺ T cell subsets. Purities (among CD8⁺ cells) and recoveries after enrichment with the CliniMACS System are shown. Numbers indicate means ± SEM (n=6). Data were adapted from reference 15.




Figure 1 Enriched TN, TCM, and TEM cells show distinct phenotypes. Enriched T cell subsets and PBMCs were analyzed for expression of CCR7, CD45RO, CD62L, CD127, CD27, and granzyme A (GranzA). Representative dot plots are shown. Numbers indicate the percentages of cells in the respective gates. *0.01 < P < 0.05; **0.005 < P < 0.01; ***P < 0.005. Data were adapted from reference 15.

and costimulatory signals. Since APCs are removed during enrichment of the T cell subsets, an alternative approach was required to effectively stimulate the enriched T cells. We tested three methods: 1) CD3 antibodies coated onto the cell culture vessel plus soluble CD28 antibody, 2) CD3 and CD28 antibodies linked to MACSiBead Particles, and 3) TransAct Reagent, a novel GMP-compatible reagent, which consists of CD3 and CD28 antibodies immobilized on a novel flexible nanomatrix. In contrast to MACSiBead Particles, the TransAct Reagent can be effectively removed from the culture by centrifugation. Thus, the reagent can be handled similarly to soluble antibodies.

CD8⁺ TN, TCM, and TEM cells were stimulated for 60 hours and then transduced with a retroviral vector encoding the anti-MART-1 TCR. Ten days after transduction, the transduction efficiency was similar for all three T cell subsets, regardless of the stimulation protocol, as determined by MART-1 tetramer staining (data not shown). We analyzed the transduced cell subsets for expression of CCR7 and CD62L after 13–15 days of expansion. CD62L was expressed at equally high levels on TN and TCM cells, whereas TEM cells showed a lower cell surface expression of this marker (fig. 2). CCR7 was highly expressed on TN cells, in contrast to TCM cells, which showed reduced surface expression, and TEM cells, which were CCR7 negative. TN cells stimulated with TransAct Reagent or MACSiBead Particles showed a higher frequency of CCR7⁺MART-1 tetramer⁺CD8⁺ cells (57.22±6.8% and 58.26±10.1%) compared to cells stimulated with CD3/CD28 (30.5±8.4%). This indicates that the phenotype, in particular expression of




Figure 2 Expression of CD62L and CCR7 in MART-1 tetramer⁺ TN, TCM, and TEM cells. The three T cell subsets were enriched from PBMCs and stimulated with CD3/CD28 antibodies, MACSiBead Particles, or TransAct Reagent in the presence of IL-2. PBMCs were stimulated with soluble CD3 and CD28 antibodies. Cells were transduced 60 hours after stimulation to express MART-1-specific TCR. Cells were analyzed by flow cytometry 13–15 days after stimulation. Data indicate frequencies of (A) MART-1 tetramer⁺CD62L⁺ cells and (B) MART-1 tetramer⁺CCR7⁺ cells among CD8⁺ cells. Data were adapted from reference 15.

the homing receptor CCR7, can be modulated by the method of cell stimulation.

Taken together, our results show that the enriched T cell subsets can be efficiently stimulated, transduced, and expanded over two weeks. However, expression patterns in expanded cells can differ from the initial cell population, and marker expression varies dependent on the stimulation method, which may have an effect on the cells' behavior *in vivo*.

Functional analysis of T cell subsets

To test whether the MART-1 TCR in transduced cells was functional, we re-stimulated the T cell subsets with a MART-1⁺HLA-A2⁺ cell line and analyzed the production of IL-2, TNF- α , and IFN- γ . The percentage of CD8⁺IL-2⁺ T_N cells was higher when they had been stimulated with MACSiBead Particles (13 \pm 4.2%) or TransAct Reagent (7.3 \pm 1.3%) compared to CD3/CD28 (3.8 \pm 0.7%) (fig. 3B). Similarly, T_N cells stimulated with MACSiBead Particles showed a higher percentage of TNF- α -producing cells (42.2 \pm 4.2%) than CD3/CD28-stimulated T_N cells (22 \pm 5.6%) (fig. 3D). The percentage of IFN- γ -producing T cells was similar for all three methods of stimulation (fig. 3C). There was a slight trend towards increased percentages of IFN- γ ⁺ cells in T_{CM} and T_{EM} cells in comparison to T_N cells.

Compared to CD3/CD28-stimulated T_N cells, the cells stimulated with MACSiBead Particles or TransAct Reagent showed an increase in IL-2, equivalent or increased TNF- α levels, and equivalent or reduced IFN- γ levels, which indicates less progress towards terminal differentiation.

Conclusion

- We developed a reliable and easy-to-perform GMP-compliant method for the enrichment of CD8⁺ T_N and T_{CM} cells from leukapheresis products.
- We developed a novel reagent for effective T cell stimulation, TransAct Reagent, which can be sterile filtered and added to cell cultures at defined concentrations. Excess reagent can be easily removed from the cells by centrifugation.
- The enriched T cell subsets can be stimulated and efficiently transduced to produce MART-1-specific cells.
- Following stimulation of T cells with TransAct Reagent or MACSiBead Particles,




Figure 3 Cytokine production by T_N, T_{CM}, and T_{EM} cells following re-stimulation with MART-1⁺ melanoma cells. T cells were stimulated using the different protocols as indicated, and transduced to express a MART-1-specific TCR. Transduced T cell subsets and PBMCs were cultured for a total of 13–15 days and then washed to remove IL-2. Cells were cultured for two more days and re-stimulated with a MART-1⁺HLA-A2⁺ melanoma cell line for 6 hours. Cytokine production was determined by intracellular staining. (A) IL-2 expression in CD8⁺ T_N cells. Numbers indicate percentages of IL-2⁺ cells among CD8⁺ cells. (B) Frequencies of CD8⁺IL-2⁺, CD8⁺IFN- γ ⁺, and CD8⁺TNF- α ⁺ cells among CD8⁺MART-1 tetramer⁺ cells. *0.01 < P < 0.05; **0.005 < P < 0.01; ***P < 0.005. Data were adapted from reference 15.

expression of CD62L and CCR7 is preserved, in particular on T_N cells, which may confer improved *in vivo* persistence to the cells.

- The enrichment and expansion procedure yields a high number of cells providing a solid basis for future translation into clinical studies.



Dissociate and discover

gentleMACS™ Dissociators

Get a good start with automated tissue dissociation

Walk-away processing

Automated heaters enable enzymatic dissociation at 37 °C

Flexible

High yields of viable single cells from virtually any tissue

Time-saving

Run up to eight programs in parallel or independently

Available for research use only in the United States.

Safe

Closed system, sterile sample handling

Reproducible

Standardized tissue dissociation and homogenization using optimized procedures

► gentlemacs.com

References

- Restifo, N.P. *et al.* (2012) *Nat. Rev. Immunol.* 12: 269–281.
- Brenner, M.K. and Heslop, H.E. (2010) *Curr. Opin. Immunol.* 22: 251–257.
- Johnson, L.A. *et al.* (2009) *Blood* 114: 535–546.
- Morgan, R.A. *et al.* (2006) *Science* 314: 126–129.
- Morgan, R.A. *et al.* (2010) *Mol. Ther.* 18: 843–851.
- Lamers, C.H. *et al.* (2011) *Blood* 117: 72–82.
- Johnson, L.A. *et al.* (2006) *J. Immunol.* 177: 6548–6559.
- Hinrichs, C.S. *et al.* (2009) *Proc. Natl. Acad. Sci. USA* 106: 17469–17474.
- Gattinoni, L. *et al.* (2005) *J. Clin. Invest.* 115: 1616–1626.
- Hinrichs, C.S. *et al.* (2011) *Blood* 117: 808–814.
- Klebanoff, C.A. *et al.* (2005) *Proc. Natl. Acad. Sci. USA* 102: 9571–9576.
- Klebanoff, C.A. *et al.* (2004) *Proc. Natl. Acad. Sci. USA* 101: 1969–1974.
- Berger, C. *et al.* (2008) *J. Clin. Invest.* 118: 294–305.
- Wang, X. *et al.* (2011) *Blood* 117: 1888–1898.
- Casati, A. *et al.* (2013) *Cancer Immunol. Immunother.* 62: 1563–1573.
- Bialer, G. *et al.* (2010) *J. Immunol.* 184: 6232–6241.
- Bunnell, B.A. *et al.* (1995) *Proc. Natl. Acad. Sci. USA* 92: 7739–7743.
- Sallusto, F. *et al.* (1999) *Nature* 401: 708–712.
- Huang, J. *et al.* (2006) *J. Immunol.* 176: 7726–7735.
- Huang, J. *et al.* (2005) *J. Immunother.* 28: 258–267.
- Geginat, J. *et al.* (2003) *Blood* 101: 4260–4266.

MACS Product	Order no.
CliniMACS Plus Instrument	151-01
CliniMACS CD4 MicroBeads	130-030-401
CliniMACS CD45RA MicroBeads	130-020-003
CliniMACS CD62L MicroBeads	130-020-006
CliniMACS PBS/EDTA Buffer	700-25
CliniMACS Tubing Set	161-01
CliniMACS Tubing Set LS	162-01
MACSiMAG Separator	130-092-168
Anti-Biotin MACSiBead Particles – cell culture grade	130-092-357
MACS GMP TransAct CD3/CD28 Kit	Coming soon

The CliniMACS System components: Reagents, Tubing Sets, Instruments and PBS/EDTA Buffer, are manufactured and controlled under an ISO 13485 certified quality system. In Europe, the CliniMACS System components are available as CE-marked medical devices. In the USA, the CliniMACS System components including the CliniMACS Reagents are available for use only under an approved Investigational New Drug (IND) application or Investigational Device Exemption (IDE). CliniMACS MicroBeads are for research use only and not for human therapeutic or diagnostic use. The MACSiMAG Separator and Anti-Biotin MACSiBead Particles are for research use only. The MACS GMP TransAct CD3/CD28 Kit is for research use and *ex vivo* cell culture processing only, and is not intended for human *in vivo* applications. For regulatory status in the USA, please contact your local representative. Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use.