

Enrichment of pluripotent stem cell derived neural crest stem cells and further differentiation to peripheral neurons

Background

Neural crest stem cells (NCSC) can be induced from pluripotent stem cells by synergistic action of two inhibitors of SMAD signaling, Noggin and SB431542^{1,2} or alternatively, by a single small molecule dorsomorphin blocking the signaling of several (TGF- β) superfamily receptors.³

CD271⁺ neural crest stem cells can be selected ten days post-induction. They show expression of typical neural crest markers (HNK1, AP2) and they can be differentiated to homogenous peripheral neurons.⁴

CD271, also known as LNGFR (low-affinity nerve growth factor receptor), NGFR (nerve growth factor receptor), or p75NTR (neurotrophin receptor), belongs to the tumor necrosis factor receptor superfamily. CD271 was initially described to be expressed on cells of the nervous system and was suggested to be involved in the development, survival and differentiation of neural cells.⁵ CD271 can be found in the central and peripheral nervous system on autonomic and sensory neurons⁶ as well as on glial cells, including oligodendrocytes⁷, astrocytes⁸, Schwann cells^{9,10}, and NCSCs¹¹.

Materials and methods

Note: For reagents, solutions, and materials required for the removal of feeder cells or the isolation of TRA-1-60⁺ pluripotent stem cells and NCSCs, please refer to the data sheets of the Feeder Removal MicroBeads, mouse or Anti-Tra-1-60 MicroBead Kit, human and Neural Crest Stem Cell MicroBeads, human, respectively.

Reagents and solutions

PEB buffer: prepare a medium containing Dulbecco's Phosphate-Buffered Saline (DPBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA.

NCSC differentiation medium: prepare a solution containing 80% DMEM/F12, 1x non-essential amino acids (100 \times), 1 mM L-glutamine, 0.1 mM β -Mercaptoethanol, 20% KnockOut[™] Serum Replacement, and 1% Penicillin/Streptomycin. Before use, add fresh 1 μ M Stemolecule Dorsomorphin. Only for initial plating of the iPS cells, add fresh 2 μ M Stemolecule Thiazovivin before use.

N-2 medium: prepare a medium containing 100% DMEM/F12, 1 mM L-glutamine, 1 \times N-2 Supplement (100 \times), 0.1% MACS[®] NeuroBrew-21 w/o VitaminA, 20 μ g/mL insulin, and 1% Penicillin/Streptomycin. Before use, add fresh 10 ng/mL Human FGF-2, premium grade and 10 ng/mL Human EGF, premium grade.

Peripheral neurons differentiation medium: prepare a medium containing 50% N-2 medium (without FGF-2 and EGF), 50% MACS Neuro Medium, and 1% MACS NeuroBrew-21 w/o VitaminA. Before use, add fresh 100 ng/mL cAMP, 10 ng/mL BDNF, 20 ng/mL NGF- β , and 100 μ M ascorbic acid.

- Anti-Tra-1-60 MicroBead Kit, human (# 130-095-816) or Feeder Removal MicroBeads, mouse (# 130-095-531)
- Neural Crest Stem Cell MicroBeads, human (# 130-097-127)
- BD Matrigel[™], hESC-qualified Matrix (BD, Cat.#354277)
- Poly-L-ornithin
- Laminin
- Fibronectin
- N-2 Supplement (100 \times) (PAA)
- KnockOut[™] Serum Replacement (Life Technologies)
- Stemolecule Dorsomorphin (# 130-095-574)
- Stemolecule Thiazovivin (# 130-095-568)
- MACS NeuroBrew-21 w/o VitaminA (# 130-097-263)
- MACS Neuro Medium (# 130-093-570)
- Human FGF-2, premium grade (# 130-093-564)
- Human EGF, premium grade (# 130-097-749)
- Human BDNF (# 130-096-285)
- Human NGF- β (# 130-093-971)
- Ascorbic acid
- cAMP
- Insulin

Coating of cell culture plate with BD Matrigel[™] for NCSC differentiation

Coat cell culture plate with BD Matrigel[™] according to the manufacturer's instructions. Briefly, thaw BD Matrigel[™] over night at 4 $^{\circ}$ C. Use cold pipettes. Resuspend 150 μ L BD Matrigel[™] in 12 culture medium, e.g., DMEM or DMEM/F12 and pipette it onto the plate. Store plate over night at 4 $^{\circ}$ C. Remove the medium from the plate just before usage.

Coating of cell culture plate for cultivation and differentiation of NCSCs to peripheral neurons

1. Cover the plate with 100 µg/mL poly-L-ornithin in PBS. Incubate over night at 37 °C.
2. Remove the poly-L-ornithin/PBS from the plate
3. Wash two times with DPBS.
4. Cover the plate with 10 µg/mL laminin and 5 µg/mL fibronectin in PBS. Incubate over night at 4 °C.
5. Remove PBS shortly before usage.

Differentiation of induced pluripotent stem (iPS) cells to NCSCs

Note: Put BD Matrigel™ coated plate to room temperature one hour prior usage.

Harvest iPS cells cultured on feeder cells. Enrich human iPS cells using the Anti-Tra-1-60 MicroBead Kit, human or remove murine feeder cells using Feeder Removal MicroBeads, mouse. For detailed information on the protocols, please refer to the respective data sheets.

Note: To maintain the cells performance it is mandatory to perform all washing steps and the elution from the column with NCSC differentiation medium supplemented with 2 µM Stemolecule Thiazovivin. Optionally, collect cells by centrifugation after elution from the column. Resuspend cell pellet in NCSC differentiation medium supplemented with 2 µM Stemolecule Thiazovivin.

1. Determine the cell number and adjust the volume accordingly.
2. Aspirate the BD Matrigel™ solution from the plate.
3. Plate cells at density of 2.5×10^4 – 5×10^4 cells/cm² in NCSC differentiation medium containing 2 µM Stemolecule Thiazovivin.
4. Culture cells in NCSC differentiation medium for ten days. Renew medium every day.

Enrichment of CD271⁺ NCSCs

Harvest of NCSCs differentiated from induced pluripotent cells (iPS). Enrich CD271⁺ NCSCs using Neural Crest Stem Cell MicroBeads, human. For detailed information on the protocol, please refer to the respective data sheet.

Note: To increase the purity of CD271⁺ cells, it is highly recommended to enrich the eluted fraction, containing CD271⁺ neural crest stem cells, over a second MS Column.

Cultivation of CD271⁺ NCSCs

1. Collect the isolated CD271⁺ NCSCs by centrifugation. Resuspend cell pellet in N-2 medium.
2. Determine the cell number and adjust the volume accordingly.
3. Remove the PBS containing the laminin and fibronectin from the plate.

4. Plate cells at density of 5×10^4 – 1×10^5 cells/cm².
5. Cultivate NCSCs in N-2 medium. Renew medium every other day. Passage cells on day 4–7. For passaging, cells can be trypsinized according to standard procedures.

Differentiation towards peripheral neurons

1. When NCSCs are confluent, change medium to peripheral neuron differentiation medium.
2. Renew medium every 2–3 days.
3. The differentiation of NCSCs to peripheral neuron takes approximately 21 days.
4. After differentiation peripheral neurons can be further processed or analyzed.




Figure 1: Peripheral neurons differentiated from enriched CD271⁺ NCSCs. Cells were stained for peripherin. Nuclei were visualized using 4',6-diamidino-2-phenylindole (DAPI).

Publications

1. Chambers, S.M. *et al.* (2009) Highly efficient neural conversion of human ES and iPS cells by dual inhibition of SMAD signaling. *Nat. Biotechnol.* 27 (3): 275–280.
2. Lee, G. *et al.* (2010) Derivation of neural crest cells from human pluripotent stem cells. *Nat. Protoc.* 5 (4): 688–701.
3. Zhou, J. *et al.* (2010) High-efficiency induction of neural conversion in human ESCs and human induced pluripotent stem cells with a single chemical inhibitor of transforming growth factor beta superfamily receptors. *Stem Cells* 28 (10): 1741–1750.
4. Schreiner, C. *et al.* (2011). Sequential magnetic enrichment of TRA-1-60⁺ PSCs and CD271⁺ (p75) NCSCs enhances peripheral neuron differentiation. Poster available at www.macs-stemcells.com/downloads.
5. Thomson, T. M. *et al.* (1988) Expression of human nerve growth factor receptor on cells derived from all three germ layers. *Exp. Cell Res.* 174 (2): 533–539.
6. Kashiba, H. *et al.* (1995) Coexpression of trk family members and low-affinity neurotrophin receptors in rat dorsal root ganglion neurons. *Brain Res. Mol. Brain Res.* 30 (1): 158–164.
7. Casaccia-Bonnel, P. *et al.* (1996) Death of oligodendrocytes mediated by the interaction of nerve growth factor with its receptor p75. *Nature* 383 (6602): 716–719.
8. Rudge, J. S. *et al.* (1994) Neurotrophic factor receptors and their signal transduction capabilities in rat astrocytes. *Eur J. Neurosci.* 6 (5): 693–705.
9. DiStefano, P. S. and Johnson, E. M. Jr. (1988) Identification of a truncated form of the nerve growth factor receptor. *Proc. Natl. Acad. Sci. USA* 85 (1): 270–274.
10. Vroemen, M. and Weidner, N. (2003) Purification of Schwann cells by selection of p75 low affinity nerve growth factor receptor expressing cells from adult peripheral nerve. *J. Neuroscience Methods* 12(2): 135–143.
11. Chalazonitis, A. *et al.* (1998) Age-dependent differences in the effects of GDNF and NT-3 on the development of neurons and glia from neural crest-derived precursors immunoselected from the fetal rat gut: Expression of GFRα-1 *in vitro* and *in vivo*. *Dev. Biol.* 204 (2): 385–406.



Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS is a registered trademark of Miltenyi Biotec GmbH. BD and BD Matrigel are trademarks of Becton, Dickinson and Company. KnockOut is a trademark of Life Technologies Corporation. Copyright © 2013 Miltenyi Biotec GmbH. All rights reserved.