

**Find the fittest**  
**MACS® ART Annexin V System**



## Easy, fast, and specific depletion of unwanted apoptotic spermatozoa

- Improve quality of sperm preparations
- Obtain desired population within 30 minutes
- Proven technology for optimal results

## Improve the quality of spermatozoa preparations

### Be specific, be successful

The presence of apoptotic cells in sperm preparations for assisted reproduction is known to decrease specimen quality. Unfortunately, commonly used techniques for preparing sperm may not remove apoptotic spermatozoa with the required specificity.

That's why the MACS® ART Annexin V System was designed to selectively remove apoptotic spermatozoa from sperm preparations. Even early apoptotic spermatozoa with an intact cellular membrane are efficiently depleted, to yield spermatozoa fractions for use in subsequent operations in Assisted Reproductive Technology (ART).

Select only the highest quality for your next sperm preparation.

### Researchers have reported that purified sperm populations show higher overall quality in terms of:

- Morphology<sup>1</sup>
- Sperm Deformity Index<sup>1</sup>
- Motility<sup>2,3</sup>
- DNA fragmentation
- Mitochondria membrane potential integrity<sup>3,4,5</sup>
- Apoptosis
- Cryosurvival rates<sup>2,4</sup>


### Scientist obtained a sperm population with:

- Improved sperm fertilization potential (sperm capacitation index, hamster oocyte penetration assay)<sup>6,12</sup>
- Improved Induced Acrosome Reaction Test (IART)<sup>7,8</sup>
- Increased sperm chromatin decondensation<sup>9,11</sup>
- Improved cleavage rates after fertilization<sup>10</sup>

## References

- 1 Aziz, N. *et al.* (2007) *Human Reproduction* 22: 1413–1419.
- 2 Said, T. *et al.* (2005) *Fertility and Sterility* 83: 1442–1446.
- 3 De Vantéry Arrighi, C. *et al.* (2009) *Reproductive Biology Endocrinology* 7: 1.
- 4 Grunewald, S. *et al.* (2006) *Cell and Tissue Banking* 7: 99–104.
- 5 Paasch, U. *et al.* (2005) *Asian Journal of Andrology* 7: 61–69.
- 6 Said, T. *et al.* (2006) *Biology of Reproduction* 74: 530–537.
- 7 Lee, T. *et al.* (2010) *Human Reproduction* 25: 839–846.
- 8 Grunewald, S. *et al.* (2006) *Ann. N. Y. Acad. Sci.* 1090: 138–146.
- 9 Grunewald, S. *et al.* (2009) *Fertility and Sterility* 92: 572–577.
- 10 Dirican, E. K. *et al.* (2008) *Journal of Assisted Reproduction and Genetics* 25: 375–381.
- 11 Delbes, G. *et al.* (2013) *Journal of Andrology* 1: 689–706.
- 12 Sheikh, A. *et al.* (2013) *Journal of Andrology* 1: 845–849.






## Specific depletion of apoptotic spermatozoa

The MACS ART Annexin V System is intended for *in vitro* depletion of apoptotic human spermatozoa. The system selectively depletes unwanted cells using the ability of Annexin V to recognize the antigen EPS (externalized phosphatidylserine) in the plasma membrane of apoptotic cells.


Phosphatidylserine-positive cells are depleted from fresh, cryopreserved, or otherwise manipulated semen samples. The yielded spermatozoa fraction is suitable for use in Assisted Reproductive Technology.

### How it works


Apoptotic cell depletion begins with magnetic labeling of apoptotic cells by the MACS ART Annexin V Reagent. The labeled cells are then passed through a separation column located in a magnetic field. Unwanted cells are selectively retained in the column. Living spermatozoa are not labeled by the reagent, so they pass through the column and are collected for later use.



**Figure 1:** Translocation of PS from the intracellular to the extracellular surface of the plasma membrane during apoptosis.



**Figure 2:** Isolation of PS-exposing cells using the MACS ART Annexin V System.




## MACS® ART Annexin V System

Depletion of apoptotic spermatozoa with the MACS ART Annexin V System is simple and fast. It adapts easily into laboratories to enhance classical sperm parameters and lower sperm apoptotic markers. Additionally, it can be combined with other sperm preparation techniques for procedural flexibility.

Components of the MACS ART Annexin V System include the MACS ART Annexin V Reagent, MACS ART Binding Buffer, and MACS ART MS Column. All are free of human and animal-derived components. Product specifications, such as endotoxin content or sterility, are confirmed by batch-specific certificates of analysis to ensure consistency of quality.

The system has been tested in a Mouse Embryo Assay (MEA), which is widely used to assess the toxicity of materials used in human assisted reproduction and to detect embryo-toxic materials. The MACS ART Annexin V System shows no embryo toxicity.



**Figure 3:** The depletion of apoptotic spermatozoa with the MACS ART Annexin V System is simple and fast.

Product	Content	Order no.
MACS ART Annexin V Reagent	1×0.1 mL reagent	200-070-502
MACS ART Annexin V Reagent	6×0.1 mL reagent	200-070-503
MACS ART Binding Buffer	6×10 mL buffer	200-070-505
MACS ART MS Column	6 columns	200-070-500
MACS ART Separation Unit	1 unit	200-070-501
MACS MultiStand	1 piece	130-042-303
MACS ART Annexin V Complete Kit	6×0.1 mL reagent, 6×10 mL buffer, 6 columns	200-070-507



## Miltenyi Biotec

### Germany/Austria/ Switzerland

Miltenyi Biotec GmbH  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone +49 2204 8306-0  
Fax +49 2204 85197  
macs@miltenyibiotec.de

### USA/Canada

Miltenyi Biotec Inc.  
2303 Lindbergh Street  
Auburn, CA 95602, USA  
Phone 800 FOR MACS  
Phone +1 530 888 8871  
Fax +1 530 888 8925  
macs@miltenyibiotec.com

### Australia

Miltenyi Biotec  
Australia Pty. Ltd.  
Unit 16A, 2 Eden Park Drive  
Macquarie Park NSW 2113  
Australia  
Phone +61 2 8877 7400  
Fax +61 2 9889 5044  
macs@miltenyibiotec.com.au

### Benelux

Miltenyi Biotec B.V.  
Schipholweg 68 H  
2316 XE Leiden  
The Netherlands  
macs@miltenyibiotec.nl

### Customer service

**The Netherlands**  
Phone 0800 4020120  
Fax 0800 4020100

### Customer service Belgium

Phone 0800 94016  
Fax 0800 99626

### Customer service Luxembourg

Phone 800 24971  
Fax 800 24984

### China

Miltenyi Biotec Technology &  
Trading (Shanghai) Co., Ltd.  
Room 2309  
No. 319, Xianxia Road  
Changning District  
200051 Shanghai, P.R. China  
Phone +86 21 62351005  
Fax +86 21 62350953  
macs@miltenyibiotec.com.cn

### France

Miltenyi Biotec SAS  
10 rue Mercoeur  
75011 Paris, France  
Phone +33 1 56 98 16 16  
Fax +33 1 56 98 16 17  
macs@miltenyibiotec.fr

### Italy

Miltenyi Biotec S.r.l.  
Via Persicetana, 2/D  
40012 Calderara di Reno (BO)  
Italy

Phone +39 051 6 460 411  
Fax +39 051 6 460 499  
macs@miltenyibiotec.it

### Japan

Miltenyi Biotec K.K.  
Nittsu-Eitai Building 5F  
16-10 Fuyuki, Koto-ku,  
Tokyo 135-0041, Japan  
Phone +81 3 5646 8910  
Fax +81 3 5646 8911  
macs@miltenyibiotec.jp

### Nordics and Baltics

Miltenyi Biotec Norden AB  
Scheelevägen 17  
223 70 Lund  
Sweden

macs@miltenyibiotec.se

### Customer service Sweden

Phone 0200-111 800  
Fax 046-280 72 99

### Customer service Denmark

Phone 80 20 30 10  
Fax +46 46 280 72 99

### Customer service

**Norway, Finland, Iceland,  
and Baltic countries**

Phone +46 46 280 72 80  
Fax +46 46 280 72 99

### Singapore

Miltenyi Biotec Asia Pacific Pte Ltd.  
100 Beach Road  
#28-06 to 28-08 Shaw Tower  
Singapore 189702  
Phone +65 6238 8183  
Fax +65 6238 0302  
macs@miltenyibiotec.com.sg

### South Korea

Miltenyi Biotec Korea  
Petra building  
659 Yeoksam-dong  
Gangnam-gu  
Seoul, South Korea  
Phone +82 2 555 1988  
Fax +82 2 555 8890  
macs@miltenyibiotec.co.kr

### Spain

Miltenyi Biotec S.L.  
C/Luis Buñuel 2  
Ciudad de la Imagen  
28223 Pozuelo de Alarcón (Madrid)  
Spain  
Phone +34 91 512 12 90  
Fax +34 91 512 12 91  
macs@miltenyibiotec.es

### United Kingdom

Miltenyi Biotec Ltd.  
Almac House, Church Lane  
Bisley, Surrey GU24 9DR, UK  
Phone +44 1483 799 800  
Fax +44 1483 799 811  
macs@miltenyibiotec.co.uk

[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

Miltenyi Biotec provides products and services worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

This print material is not for use in the USA. The MACS ART Annexin V System is manufactured by Miltenyi Biotec GmbH and tested under an ISO 13485 certified quality system. These products are available in Europe as CE-marked medical devices. For availability in your country, please contact your local sales representative. MACS is a registered trademark of Miltenyi Biotec GmbH. Copyright © 2014 Miltenyi Biotec GmbH. All rights reserved