

## Contents

5 MACS D Columns for SuperMACS and SuperMACS II plus accessories. Sterile packed. Capacity: up to  $10^9$  magnetically labeled cells.

## Applications

- Depletion of large numbers of cells labeled with MACS MicroBeads.
- Magnetic separation of biological material labeled with MACS MicroBeads.

## Storage of MACS Separation Columns

Store columns dry, protected from light. Do not use after expiry date.

## Instrument and Reagent Requirement

Magnetic cell separator SuperMACS or SuperMACS II.  
Column Adapter for CS and D Columns in combination with SuperMACS II.

MACS MicroBeads for magnetic labeling of the cells.

70 % ethanol in dd  $H_2O$ .

Buffer: phosphate buffered saline supplemented with 2 mM EDTA and 0.5 % bovine serum albumin.

50 ml syringe.


## How to Use MACS Depletion Column Type D

D Columns have been developed for depletion of large amounts of magnetically labeled human and animal cells out of a heterogeneous suspension with SuperMACS or SuperMACS II. They can be used to separate different biological material including plant cells, bacteria, viruses, protozoa or cell organelles. Various running buffers may be used with the depletion column. The suitability of a specific buffer has to be tested experimentally. By passing the material through the column, the suspension is depleted for magnetically labeled material.


## Protocol for Cell Separation Using D Columns

### Assembly and Preparation of the Depletion Column

3-way stopcock positions



### Assembly of the D Column



- Remove yellow cap from separation column and attach 3-way stopcock to column at port A.
- Fill a syringe with **70 % ethanol** and attach to port B of the 3-way stopcock.
- Turn 3-way stopcock to position "fill".
- Move the gap of the SuperMACS to approximately 3 cm width using the side wheel. Place assembled column in MACS separator with 3-way stopcock in adjustment and secure with the lever.
- When using in combination with SuperMACS II, insert column into the mounted Column Adapter for CS and D Columns and move Column Adapter into the magnetic field by turning the handle (for details, see SuperMACS II data sheet).

- Fill the column upright from the bottom with 70 % ethanol from the syringe until the solution reaches the reservoir.
- Turn the 3-way stopcock to position “run” and rinse column by filling from the top with buffer. Allow 70 % ethanol to run into the column until it reaches the white filter device. Then add fresh buffer. Rinse with 500 ml of buffer. Do not let the column “run dry” at any time during procedure.
- Choose a flow resistor (22G; 21G or 20G, see table below) and cut off the tip of the plastic sheath with the pliers supplied with the SuperMACS Starting Kit. Leave the plastic sheath in place for safety and attach the flow resistor to port C of the 3-way stopcock.

Caution: Never remove plastic sheath from flow resistor.

Flow resistor	Flow rates in ml/min for buffer
22G	3.5
21G	4.0
20G	8.0
Note: The lower the flow rate, the higher is the depletion efficiency of magnetically labeled cells.	

- Fill the syringe with **buffer** and attach to port B of the 3-way stopcock. Leave the syringe attached during separation, except when refilling.
- Turn the stopcock to position “rinse” and flush the air from the flow resistor. Turn the stopcock in position “closed”. The column is now ready for separation.

### Depletion of Cells

- Apply magnetically labeled cell suspension in appropriate volume of buffer (up to  $10^8$  cells per 500  $\mu$ l) onto the column with the flow resistor attached. Turn the stopcock to position “run” and allow cell suspension to penetrate the matrix.
- Wash the column with 200 ml of buffer. Collect effluent as negative fraction that is depleted of the magnetically labeled cells.

### Important Notes

- ▲ Never use a flow resistor without plastic cover.
- ▲ Recommended buffer is PBS supplemented with 2 mM EDTA and 0.5 % BSA. Other buffers have to be tested for their flow conditions. Do not use buffers of too high viscosity.
- ▲ Fill column properly avoiding air bubbles. Small air bubbles underneath the filter can be neglected.
- ▲ Column must not “run dry” at any time during the procedure.
- ▲ Do not store column after filling.
- ▲ Columns can be cooled just prior to use by passing 2–3 column volumes of ice-cold buffer through it.
- ▲ The number of cells loaded onto the column can be up to  $10^8$  per 500  $\mu$ l.

- ▲ When working with anti-coagulated blood or buffy coat, dilute before separation with buffer.
- ▲ Do not use separation columns in combination with magnetic beads other than MACS MicroBeads. Magnetic forces in the column are very high and may damage biological material, if other beads are used.
- ▲ If the buffer does not flow well, there may be an air bubble in the capillary of the flow resistor. Switch the stopcock to “rinse”, flush the needle with buffer from the side syringe, switch the stopcock back to “run” and continue the run. The needle may also become blocked by cell clumps, which means the needle has to be replaced.
- ▲ If necessary, it is possible to pass your cells through the same column a second time.

### Technical Specifications

- Matrix volume: 43 ml; reservoir volume: 20 ml.
- Typical capacity: e.g.  $10^9$  retained lymphocytes and up to  $10^{11}$  total applied cells.
- Recommended sample size for leukocytes:  $10^7$ - $10^9$  magnetically labeled cells in  $10^9$ - $10^{11}$  total cells. Sample concentration: up to  $10^8$  leukocytes/500  $\mu$ l cell suspension.

Refer to [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

## Legal notices

### Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com), as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

### Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for the most up-to-date information on Miltenyi Biotec products.

### Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at [www.miltenyibiotec.com](http://www.miltenyibiotec.com) for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

### Trademarks

MACS, the Miltenyi Biotec logo, and SuperMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.