




Miltényi Biotec

IL-12 Secretion Assay – Cell Enrichment and Detection Kit (PE) human

For 50 tests with 10^7 cells

Order no. 130-092-122



Unless otherwise specifically indicated, Miltényi Biotec products and services are for research use only and not for diagnostic or therapeutic use.



Miltényi Biotec

Miltényi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68, 51429 Bergisch Gladbach, Germany
Phone +49 2204 8306-0, Fax +49 2204 85197
macsde@miltényi.com, www.miltényibiotec.com

140-001-451.03

Contents

1. Description


Contents

1. Description
 - 1.1 Principle of the IL-12 Secretion Assay
 - 1.2 Background information
 - 1.3 Applications
 - 1.4 Reagent and instrument requirements
2. Protocol overview
3. Experimental setup
 - 3.1 Controls
 - 3.2 Counterstaining of cytokine-secreting cells
 - 3.3 Two-color cytokine analysis
 - 3.4 Detection without prior enrichment
4. Protocol for the IL-12 Secretion Assay
 - 4.1 Cell preparation
 - 4.2 *In vitro* stimulation
 - 4.3 Cytokine Secretion Assay
 - 4.4 Magnetic labeling
 - 4.5 Magnetic separation
5. Detection and analysis of IL-12-secreting cells
6. References

7. Appendix: Flask and dish sizes for *in vitro* stimulation

1. Description

Components	<p>1 mL IL-12 Catch Reagent, human: anti-IL-12 monoclonal antibody (mouse IgG1) conjugated to cell surface-specific monoclonal antibody (mouse IgG2a).</p> <p>1 mL IL-12 Detection Antibody, human: anti-IL-12 monoclonal antibody (mouse IgG1) conjugated to PE (R-phycoerythrin).</p> <p>1 mL Anti-PE MicroBeads: MicroBeads conjugated to monoclonal anti-PE antibodies (isotype: mouse IgG1).</p>
Size	For 50 tests with 10^7 cells.
Product format	All components are supplied in buffer containing stabilizer and 0.05% sodium azide.
Storage	Store protected from light at 2–8 °C. Do not freeze. The expiration dates are indicated on the vial labels.



1.1 Principle of the IL-12 Secretion Assay

IL-12-secreting cells are analyzed and isolated using the IL-12 Secretion Assay starting from whole blood, peripheral blood mononuclear cells (PBMCs), or other leukocyte-containing single-cell preparations.

The cells are restimulated for a certain period of time with the desired stimulus. Subsequently, an IL-12-specific Catch Reagent is attached to the cell surface of all leukocytes. The cells are then incubated for a short time at 37 °C to allow cytokine secretion. The secreted IL-12 binds to the IL-12 Catch Reagent on the positive, secreting cells. These cells are subsequently labeled with a second IL-12-specific antibody, the IL-12 Detection Antibody conjugated R-phycoerythrin (PE) for sensitive detection by flow cytometry.

The IL-12-secreting cells can now be magnetically labeled with Anti-PE MicroBeads and enriched over a MACS® Column which is placed in the magnetic field of a MACS Separator. The magnetically labeled

cells are retained within the MACS Column while the unlabeled cells run through. After removing the column from the magnetic field, the magnetically retained cells can be eluted as the positively selected cell fraction, enriched for cytokine-secreting cells. The cells can now be used for cell culture or analysis. Since viable cells are analyzed, non-specific background can be minimized by dead cell exclusion. This provides highest sensitivity of analysis.

1.2 Background information

Interleukin 12 (IL-12) consists of a p35 subunit and a p40 subunit that form the biologically active heterodimer IL-12p70. It is produced mainly by monocytes and macrophages, but also by dendritic cells, primarily in response to stimulation with bacterial or viral antigens.

IL-12p70 is produced, e.g., in response to LPS with costimulation of T helper cells through CD40L or IFN- γ .¹ Moreover, monocytes can be induced to synthesize IL-12p70 without costimulation of T cells, through combined stimulation with toll-like receptor (TLR) TLR-4 and TLR-8 ligands.² IL-12 displays various effects on both T and NK cells, playing an essential role in the interaction between innate and adaptive immunity.

The IL-12 Cytokine Secretion Assay is designed to detect and enrich cells, secreting the active heterodimer IL-12p70.

Quantitative analysis of IL-12-secreting leukocytes can provide important information on the natural course of immune responses. Enrichment of the IL-12-secreting leukocytes using the MACS

Technology increases the sensitivity of analysis and enables further functional characterization, culturing of the cells, as well as downstream applications.

1.3 Applications

- Detection and enrichment of viable IL-12-secreting leukocytes.
- Detection and enrichment of IL-12-secreting cells for enumeration and phenotypic analysis as well as for culturing and functional characterization.
- Monitoring and analysis of the stimulatory or regulatory functions of cytokine-secreting monocytes, macrophages, and dendritic cells.

1.4 Reagent and instrument requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C). Degas buffer before use, as air bubbles could block the column.
- Culture medium, e.g., RPMI 1640 (# 130-091-440), containing 5% human serum, like autologous or AB serum (do not use BSA or FCS because of non-specific stimulation!).


- Propidium Iodide Solution (# 130-093-233) or 7-AAD for flow cytometric exclusion of dead cells.
- Dead Cell Removal Kit (# 130-090-101) for the depletion of dead cells.
- Refrigerated centrifuge (2–8 °C).
- Rotation device for tubes: MACSmix™ Tube Rotator (# 130-090-753).
- MACS Columns and MACS Separators:

Column	Max. number of labeled cells	Max. number of total cells	Separator
MS	10 ⁷	2 × 10 ⁸	MiniMACS, OctoMACS, VarioMACS, SuperMACS
LS	10 ⁸	2 × 10 ⁹	MidiMACS, QuadroMACS, VarioMACS, SuperMACS
autoMACS	2 × 10 ⁸	4 × 10 ⁹	autoMACS Pro, autoMACS

▲ Note: Column adapters are required to insert certain columns into VarioMACS™ Separator or SuperMACS™ Separator. For details, see MACS Separator data sheets.

- (Optional) Cell stimulation reagents, for example, LPS and R-848 for restimulation of human cells.
- (Optional) Staining reagents such as CD14-FITC (# 130-080-701), CD14-APC (# 130-091-243), CD11c-FITC (# 130-092-410), or CD11c-APC (# 130-092-412). For more information about fluorochrome-conjugated antibodies see www.miltenyibiotec.com.
- (Optional) Pre-Separation Filters (# 130-041-407) to remove cell clumps.

2. Protocol overview



3. Experimental setup

3.1 Controls

Negative control

For accurate detection of IL-12-secreting cells, a negative control sample should always be included. This will provide information about IL-12 secretion unrelated to the *in vitro* stimulation but, for example, due to ongoing *in vivo* immune responses. The control sample should be treated exactly the same as the antigen-stimulated sample except for the addition of stimulus, or by using a control stimulus.

3.2 Counterstaining of cytokine-secreting cells

The IL-12-secreting cells are stained with PE-conjugated IL-12 Detection Antibody. To identify cells of interest, counterstaining for cells with, for example, CD14-FITC (# 130-080-701), CD14-APC (# 130-091-243), CD11c-FITC (# 130-092-410), or CD11c-APC (# 130-092-412) is important.

▲ Do not use tandem conjugates of phycoerythrin as they may also be recognized by the Anti-PE MicroBeads.

▲ The samples should be stained with Propidium Iodide Solution (# 130-093-233) or 7-AAD prior to acquisition, to exclude dead cells from analysis. This will reduce non-specific background staining and increase sensitivity.

3.3 Two-color cytokine analysis

IL-12-secreting cells can be analyzed simultaneously for, e.g., IL-10 production by two-color cytokine analysis combining the IL-12 Secretion Assay – Detection Kit (PE) with the IL-10 Secretion Assay – Detection Kit (APC) (# 130-090-761). Detailed protocols are included in the data sheets of the Cytokine Secretion Assay – Detection Kits (APC) and are available at www.miltenyibiotec.com/protocols.

3.4 Detection without prior enrichment

(Optional) If the sample contains more than 0.01–0.1% of IL-12-secreting cells, you may be able to analyze IL-12-secreting cells without prior enrichment (see also: IL-12 Secretion Assay – Detection Kit (PE), # 130-092-124). The assay can also be performed directly starting from whole blood. For details on the procedure refer to www.miltenyibiotec.com/protocols.

4. Protocol for the IL-12 Secretion Assay

4.1 Cell preparation

To detect and isolate cytokine-secreting cells, best results are achieved by starting the assay with fresh PBMCs, or with other leukocyte-containing single-cell preparations from tissues or cell lines. Alternatively, frozen cell preparations can be used.

▲ **Note:** PBMCs may be stored overnight. The cells should be resuspended and incubated in culture medium as described in 4.2 step 2. However, the stimulus should not be added to the culture until the next day.

▲ **Note:** Remove platelets after density gradient separation. Resuspend cell pellet, fill tube with buffer, and mix. Centrifuge at 200×g for 10–15 minutes at 20 °C. Carefully remove supernatant.

Special protocols for whole blood: You can start the IL-12 Secretion Assay directly from whole blood. For details on the procedure refer to www.miltenyibiotec.com/protocols.

4.2 In vitro stimulation

▲ Always include a negative control in the experiment (see 3.1). A positive control may also be included.

▲ Do not use media containing any non-human proteins, such as BSA or FCS, as they lead to non-specific stimulation.



Protocol for *in vitro* stimulation

1. Wash cells by adding medium, centrifuge at 300×g for 10 minutes.
2. Resuspend cells in culture medium, containing 5% human serum, adjust to 10^6 cells/mL and 5×10^5 cells/cm² (see 7. Appendix A: Flask and dish sizes for *in vitro* stimulation).
3. Add stimulus or control reagent, for example:
LPS (TLR-4 agonist) and R-848 (TLR-7/-8 agonist): LPS (e.g. 1 µg/mL) and R-848 (e.g. 2.5 µg/mL) for 4 hours at 37 °C, 5–7% CO₂.
For comparison of different experiments, the stimulation time should always be the same.
4. Collect cells carefully by using a cell scraper. Rinse the dish with cold buffer. Check microscopically for any remaining cells, if necessary, rinse the dish again.

4.3 Cytokine Secretion Assay

General considerations

▲ The assay is optimized for cell samples containing <5% of total IL-12-secreting cells. If ≥5% of IL-12-secreting cells are expected, it is necessary to dilute the cells further during the cytokine secretion period, and therefore a larger test tube will be needed (see table below).

The dilution prevents non-specific staining of cells not secreting IL-12 during this period.

- ▲ For each test with 10^7 total cells, prepare:
100 mL of cold buffer (2–8 °C)
100 µL of cold medium (2–8 °C)
20 mL (or 100 mL; see table below) of warm medium (37 °C).

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling (exception: warm medium during secretion period).

▲ Volumes given below are for up to 10^7 total cells. When working with fewer than 10^7 cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10^7 total cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ Do not remove supernatant by decanting. This will lead to cell loss and incorrect incubation volumes. Aspirate supernatant.

▲ Dead cells may bind non-specifically to antibodies. Therefore, when working with cell preparations containing large amounts of dead cells, they should be removed before starting the IL-12 Secretion Assay, e.g., by density gradient centrifugation or by using the Dead Cell Removal Kit (# 130-090-101).



Labeling cells with IL-12 Catch Reagent

1. Use 10^7 total cells in a 15 mL closable tube per sample.
2. Wash cells by adding 10 mL of cold buffer, centrifuge at 300×g for 10 minutes at 2–8 °C, aspirate supernatant completely.
3. Resuspend cell pellet in 80 µL of cold medium per 10^7 total cells.
4. Add 20 µL of IL-12 Catch Reagent per 10^7 total cells, mix well, and incubate for 5 minutes on ice.



IL-12 secretion period

1. Add warm (37 °C) medium to dilute the cells according to the following table:

Expected number of IL-12-secreting cells	Dilution	Amount of medium to add per 10^7 total cells
<5%	5×10^5 cells/mL	20 mL
≥5%	$\leq 10^5$ cells/mL	100 mL

▲ **Note:** For frequencies of cytokine-secreting cells >20% the cells need to be further diluted, e.g., by a factor of 5.

▲ **Note:** Choose an appropriate tube to allow for addition of cold medium to stop the secretion period.

2. Incubate cells in closed tube for 45 minutes at 37 °C under slow continuous rotation using the MACSmix Tube Rotator (# 130-090-753), or turn tube every 5 minutes to resuspend settled cells.

▲ **Note:** During this step it is crucial to prevent contact of cells to avoid cross contamination with cytokines.



Labeling cells with IL-12 Detection Antibody

1. Put the tube on ice.
2. Wash the cells by filling up the tube with cold buffer and centrifuge at 300×g for 10 minutes at 2–8 °C. Aspirate supernatant completely.
▲ **Note:** If the volume of the cell suspension was higher than the volume of the added buffer, then repeat the wash step.
3. Resuspend cell pellet in 80 µL of cold buffer per 10^7 total cells.
4. Add 20 µL of IL-12 Detection Antibody per 10^7 total cells.
5. (Optional) Add staining antibodies, e.g., 10 µL of CD14-FITC (# 130-080-701), CD14-APC (# 130-091-243), CD11c-FITC (# 130-092-410), or CD11c-APC (# 130-092-412).
6. Mix well and incubate for 10 minutes on ice.
7. Add 10 mL of cold buffer and centrifuge at 300×g for 10 minutes at 2–8 °C. Aspirate supernatant completely.

4.4 Magnetic labeling**Magnetic labeling with Anti-PE MicroBeads**

1. Resuspend cell pellet in 80 μ L of cold buffer per 10^7 total cells.
2. Add 20 μ L of Anti-PE MicroBeads per 10^7 total cells, mix well, and incubate for 15 minutes in the refrigerator (2–8 °C).
3. Wash cells by adding 10 mL of cold buffer, centrifuge at 300 \times g for 10 minutes at 2–8 °C. Aspirate supernatant completely.
4. Resuspend cell pellet in 500 μ L of cold buffer. For higher cell numbers than 5×10^7 use a dilution of 10^8 cells/mL.
5. (Optional) Take an aliquot of the fraction before enrichment for flow cytometric analysis and cell count.
6. Proceed to magnetic separation (see 4.5).

4.5 Magnetic separation**Magnetic separation using MS or LS Columns**

- ▲ Choose an appropriate MACS Column and MACS Separator according to the number of total cells. For details see table in section 1.4.

- ▲ Always wait until the column reservoir is empty before proceeding to the next step.

- ▲ To achieve highest purities, perform two consecutive column runs.

1. Place MS or LS Column in the magnetic field of a suitable MACS Separator. For details see the respective MACS Column data sheet.
2. Prepare column by rinsing with the appropriate amount of buffer:

MS: 500 μ L LS: 3 mL

3. Apply cell suspension onto the column. Collect flow-through containing unlabeled cells.
4. Wash column with the appropriate amount of buffer. Collect unlabeled cells that pass through and combine with the effluent from step 3.

MS: 3 \times 500 μ L LS: 3 \times 3 mL

5. Remove column from the separator and place it on a suitable collection tube.

▲ **Note:** To perform a second column run, you may elute the cells directly from the first onto the second, equilibrated column instead of a collection tube.

6. Pipette 1 mL of buffer onto the column. Immediately flush out the magnetically labeled cells by firmly pushing the plunger into the column.

MS: 1 mL LS: 5 mL

7. To increase purity of IL-12-secreting cells, enrich the eluted fraction over a second column. Repeat the magnetic separation procedure as described in steps 1 to 6 by using a new column.

▲ **Note:** For subsequent cell culture, the cells can also be eluted with medium. If part of the cells are analyzed by flow cytometry, the medium should not contain phenol red.

8. Proceed to analysis (see section 5.), cell culture, or other subsequent experiments.

**Magnetic separation with the autoMACS Pro Separator or the autoMACS Separator**

▲ Refer to the respective user manual for instructions on how to use the autoMACS Pro Separator or the autoMACS Separator.

▲ Buffers used for operating the autoMACS Pro Separator or the autoMACS Separator should have a temperature of ≥ 10 °C.

▲ Program choice depends on the isolation strategy, the strength of magnetic labeling, and the frequency of magnetically labeled cells. For details refer to the section describing the cell separation programs in the respective user manual.

Magnetic separation with the autoMACS Pro Separator

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube in row A of the tube rack and the fraction collection tubes in rows B and C.

3. For a standard separation choose the following program:

Positive selection: "Posseld"

Collect positive fraction in row C of the tube rack.

Magnetic separation with the autoMACS Separator

1. Prepare and prime the instrument.
2. Apply tube containing the sample and provide tubes for collecting the labeled and unlabeled cell fractions. Place sample tube at the uptake port and the fraction collection tubes at port neg1 and port pos2.

3. For a standard separation choose the following program:

Positive selection: "Posseld"

Collect positive fraction from outlet port pos2.

5. Detection and analysis of IL-12-secreting cells

▲ Add Propidium Iodide Solution (# 130-093-233) or 7-AAD to a final concentration of 0.5 μ g/mL just prior to acquisition to exclude dead cells from flow cytometric analysis. Incubating with propidium iodide for longer periods will affect the viability of the cells. Do not fix the cells when using Propidium Iodide Solution or 7-AAD.

▲ For optimized sensitivity, an appropriate number of viable cells has to be acquired from the stimulated sample as well as from the control sample.

- Acquire 2×10^5 viable cells from each sample.

To illustrate the analysis, we describe the detection of IL-12-secreting cells using the IL-12 Secretion Assay. The detailed description, including how to set gates, should serve as a model for the analysis of your own sample.


1. 10^7 human PBMCs have been stimulated for 4 hours with LPS (1 $\mu\text{g}/\text{mL}$) and R-848 (2.5 $\mu\text{g}/\text{mL}$) or have been left untreated.
2. The IL-12 Secretion Assay was performed on the stimulated and the unstimulated sample.
3. Counterstaining of cells was performed using CD14-FITC.
4. Dead cells were stained with Propidium Iodide Solution, which was added just prior to flow cytometric analysis to a final concentration of 0.5 $\mu\text{g}/\text{mL}$.
5. 200,000 viable cells were acquired by flow cytometry, from the stimulated as well as from the unstimulated sample.
6. A leukocyte gate based on forward and side scatter properties (FSC/SSC) was activated prior to further gating to exclude debris (A).
7. Dead cells were excluded according to PI-staining in a fluorescence 2 (PE) versus fluorescence 3 (PI) plot (B).

The dead cell exclusion is crucial for the analysis of rare cells, as dead cells may bind non-specifically to antibodies or MicroBeads. This could lead to false positive events.


8. Analysis of secreted IL-12 (PE) versus CD14-FITC staining of viable leukocytes is displayed (C).

A. Leukocyte gate using FSC versus SSC

Before enrichment




After enrichment




B. Dead cell exclusion

Before enrichment




After enrichment




C. Counterstaining of IL-12-secreting cells using CD14-FITC

Stimulated Sample

Before enrichment




After enrichment




Unstimulated control sample

Before enrichment



After enrichment



6. References

1. Krug, A. *et al.* (2001) Toll-like receptor reveals CpG DNA as a unique microbial stimulus for plasmacytoid dendritic cells which synergizes with CD40 ligand to induce high amounts of IL-12. *Eur. J. Immunol.* 31: 3026–3037. [1215]
2. Bekeredjian-Ding, I. *et al.* (2006) T Cell-Independent, TLR-Induced IL-12p70 Production in Primary Human Monocytes. *J. Immunol.* 176: 7438–7446. [8765]

7. Appendix: Flask and dish sizes for *in vitro* stimulation

For *in vitro* stimulation (see 4.2) the cells should be resuspended in culture medium, containing 5% of human serum, at a dilution of 10^6 cells/mL. The cells should be plated at a density of 5×10^5 cells/cm². Both the dilution and the cell density are important to assure optimum stimulation.

The following table lists culture plate, dish and flask sizes suitable for different cell numbers. It also indicates the appropriate amount of medium to add.

Total cell number	Medium volume to add	Culture plate	Well diameter
0.15×10^6	0.15 mL	96 well	0.64 cm
0.50×10^6	0.50 mL	48 well	1.13 cm
1.00×10^6	1.00 mL	24 well	1.60 cm
2.00×10^6	2.00 mL	12 well	2.26 cm
5.00×10^6	5.00 mL	6 well	3.50 cm
Total cell number	Medium volume to add	Culture dish	Dish diameter
4.5×10^6	4.5 mL	small	3.5 cm
10.0×10^6	10.0 mL	medium	6 cm
25.0×10^6	25.0 mL	large	10 cm
50.0×10^6	50.0 mL	extra large	15 cm
Total cell number	Medium volume to add	Culture flask	Growth area
12×10^6	12 mL	50 mL	25 cm ²
40×10^6	40 mL	250 mL	75 cm ²
80×10^6	80 mL	720 mL	162 cm ²
120×10^6	120 mL	900 mL	225 cm ²

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, MACSmix, MidiMACS, the Miltenyi Biotec logo, MiniMACS, OctoMACS, QuadroMACS, SuperMACS, and VarioMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.