

Tandem Signal Enhancer human

0.3 mL
1 mL

Order no. 130-099-887
Order no. 130-099-888

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
 - 1.3 Recommended antibody dilution
 - 1.4 Reagent requirements
2. Protocol
 - 2.1 Sample preparation
 - 2.2 Use of Tandem Signal Enhancer and antibodies for the labeling of human PBMCs
 - 2.3 Use of Tandem Signal Enhancer and antibodies for the labeling of human whole blood
3. Examples of immunofluorescent staining on whole blood using the Tandem Signal Enhancer

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

1. Description

This product is for research use only.

Components	0.3 mL of Tandem Signal Enhancer or 1 mL of Tandem Signal Enhancer
Capacity	300 µL: 30 tests of PBMCs and whole blood or up to 3×10^8 total cells 1 mL: 100 tests of PBMCs and whole blood or up to 10^9 total cells
Product format	All products are supplied in buffer containing stabilizer and 0.05% sodium azide.
Storage	Store protected from light at 2–8 °C. Do not freeze. The expiration date is indicated on the vial label.

1.1 Background information

When used for labeling human peripheral blood mononuclear cells (PBMCs) or whole blood tandem-dye-conjugated antibodies often show non-specific binding, especially to monocytes. The Tandem Signal Enhancer strongly reduces non-specific binding of tandem conjugates. Accordingly, negative cells show decreased background staining, which results in higher stain indices for target cell populations. The additional use of FcR Blocking Reagent is not required.

1.2 Applications

- Blocking of tandem-dye-conjugated antibodies binding to monocytes and non-target cells. Suitable for tandem dyes, such as PerCP-Vio® 700, PerCP-Cy® 5.5, PE-Cy5.5, PE-Vio 770, PE-Cy7, APC-Vio 770, APC-Cy7, APC-H7, ECD, PE-Vio 615, BD Horizon™ PE-CF594, PE/Dazzle™ 594, or PE-eFluor® 610.

1.3 Recommended antibody dilution

The recommended Tandem Signal Enhancer dilution is **1:11 for up to 10^7 cells/100 µL** of buffer for labeling of cells and subsequent analysis by flow cytometry.

1.4 Reagent requirements

- Buffer: Prepare a solution containing phosphate-buffered saline (PBS), pH 7.2, 0.5% bovine serum albumin (BSA), and 2 mM EDTA by diluting MACS® BSA Stock Solution (# 130-091-376) 1:20 with autoMACS® Rinsing Solution (# 130-091-222). Keep buffer cold (2–8 °C).
▲ **Note:** EDTA can be replaced by other supplements such as anticoagulant citrate dextrose formula-A (ACD-A) or citrate phosphate dextrose (CPD). BSA can be replaced by other proteins such as human serum albumin, human serum, or fetal bovine serum (FBS). Buffers or media containing Ca^{2+} or Mg^{2+} are not recommended for use.
- Red Blood Cell Lysis Solution (10×) (# 130-094-183).
- (Optional) For labeling of human whole blood: MACSQuant® Storage Solution (# 130-092-748).
- (Optional) For labeling of human whole blood: fluorochrome-conjugated CD45 Antibody, anti-human. For antibodies for additional staining refer to www.miltenyibiotec.com/antibodies.
- (Optional) Propidium Iodide Solution (# 130-093-233) for flow cytometric exclusion of dead cells without fixation.
- (Optional) Fixation and Dead Cell Discrimination Kit (# 130-091-163) for cell fixation and flow cytometric exclusion of dead cells.

2. Protocol

2.1 Sample preparation

When working with anticoagulated peripheral blood or buffy coat, PBMCs should be isolated by density gradient centrifugation, for example, using Ficoll-Paque™.

▲ **Note:** To remove platelets after density gradient separation, resuspend cell pellet in buffer and centrifuge at 200×g for 10–15 minutes at 20 °C. Carefully aspirate supernatant. Repeat washing step.

When working with tissues or lysed blood, prepare a single-cell suspension using standard methods.

For details refer to the protocols section at www.miltenyibiotec.com/protocols.

▲ Work fast, keep cells cold, and use pre-cooled solutions. This will prevent capping of antibodies on the cell surface and non-specific cell labeling.

▲ Volumes given below are for up to 10^7 nucleated cells. When working with fewer than 10^7 cells, use the same volumes as indicated. When working with higher cell numbers, scale up all reagent volumes and total volumes accordingly (e.g. for 2×10^7 nucleated cells, use twice the volume of all indicated reagent volumes and total volumes).

▲ For optimal performance it is important to obtain a single-cell suspension before labeling. Pass cells through 30 μ m nylon mesh (Pre-Separation Filters (30 μ m), # 130-041-407) to remove cell clumps which may clog the column. Moisten filter with buffer before use.

▲ The recommended incubation temperature is 2–8 °C. Higher temperatures and/or longer incubation times may lead to non-specific cell labeling. Working on ice may require increased incubation times.

2.2 Use of Tandem Signal Enhancer and antibodies for the labeling of human PBMCs

1. Determine cell number.
2. Centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely.
3. Resuspend up to 10^7 nucleated cells per 90 μ L of buffer.
▲ **Note:** If staining with several antibodies (see step 5), adjust the volume of buffer accordingly. Recommended dilution is 1:11.
4. Add 10 μ L of Tandem Signal Enhancer.
5. Add staining antibodies as recommended by the manufacturer. Using MACS Antibodies add the antibody according to the respective data sheet.
▲ **Note:** If staining with several antibodies, adjust the volume of buffer accordingly. Recommended dilution is 1:11.
6. Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).
7. Wash cells by adding 1 mL of buffer and centrifuge at 300×g for 10 minutes. Aspirate supernatant completely.
8. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry or fluorescence microscopy.


2.3 Use of Tandem Signal Enhancer and antibodies for the labeling of human whole blood

1. Prepare a 4× lysis solution by diluting Red Blood Lysis Solution (10×) 1:2.5 with double-distilled water (ddH₂O).
2. Pipette all conjugates into a reaction vessel in a titer as recommended by the manufacturer. For a better enumeration always use CD45 antibodies for triggering.
3. Add 10 μ L of Tandem Signal Enhancer.
4. Adjust volume to 60 μ L with buffer.
▲ **Note:** If staining with several antibodies, adjust the volume of buffer accordingly. Recommended dilution is 1:11.
5. Add 50 μ L of whole blood.
▲ **Note:** Best suitable is whole blood supplemented with EDTA or Heparin.
6. Mix well and incubate for 10 minutes in the refrigerator (2–8 °C).

7. Add 250 μ L of 4× lysis solution and 750 μ L of MACSQuant Storage Solution.
8. Mix well and incubate for 10 minutes in the dark in the refrigerator (2–8 °C).
9. Use lysed blood directly (lyse – no wash) in flow cytometry. Alternatively, centrifuge cell suspension at 300×g for 10 minutes. Aspirate supernatant completely (lyse – wash).
10. Resuspend cell pellet in a suitable amount of buffer for analysis by flow cytometry or fluorescence microscopy.

3. Example of immunofluorescent staining on whole blood using the Tandem Signal Enhancer

Fresh lysed human whole blood was stained with CD45 antibodies for targeting on all CD45⁺ cells and analyzed using the MACSQuant Analyzer. Then, a gate was set on all lymphocytes and monocytes. The histogram shows staining with CD20-PE-Vio 770 without supplement (black line) and with the use of Tandem Signal Enhancer (purple line).



Conjugate	Stain index without Tandem Signal Enhancer	Stain index with Tandem Signal Enhancer
CD20-PE-Vio 770	37.4	87.2

For more examples please refer to the respective product page at www.miltenyibiotec.com/antibodies.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

autoMACS, MACS, MACSQuant, the Miltenyi Biotec logo, and Vio are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only.

Ficoll-Paque is a trademark of GE Healthcare companies.

Cy is a registered trademark of GE Healthcare UK Limited.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.