



Miltenyi Biotec

MACSPlex Cytokine Basic Kit

human and mouse

For up to 100 tests

Order no. 130-109-701



Miltenyi Biotec B.V. & Co. KG
Friedrich-Ebert-Straße 68, 51429 Bergisch Gladbach, Germany
Phone +49 2204 8306-0, Fax +49 2204 85197
macsde@miltenyi.com, www.miltenyibiotec.com

MACSPlex-500-014

Unless otherwise specifically indicated,
Miltenyi Biotec products and services are for research use only
and not for diagnostic or therapeutic use.

Content

1. Description

Content

1. Description
 - 1.1 Principle of MACSPlex Cytokine Assays
 - 1.2 Principle of the MACSPlex Cytokine Basic Kit
 - 1.3 Applications
 - 1.4 Reagent and instrument requirements
2. Protocols for assay performance
 - 2.1 Preparation of the MACSPlex Cytokine Standard
 - 2.2 Sample preparation
 - 2.3 MACSPlex Cytokine Assay
3. Flow cytometer set up
 - 3.1 Setup of the MACSQuant® Instrument
 - 3.2 Setup of other flow cytometers
4. Flow cytometric acquisition and data analysis using the MACSQuant® Express Mode
5. Troubleshooting

1. Description

This product is for research use only.

Components	2×100 mL MACSPlex Buffer
	0.35 mL MACSPlex Setup Beads
Size	1 MACSPlex Filter Plate
	2 adhesive foils
Product format	up to 100 tests
	MACSPlex Setup Beads are supplied in buffer containing stabilizer and 0.05% sodium azide. MACSPlex Buffer contains stabilizer and 0.09% sodium azide.
Storage	Store MACSPlex Setup Beads protected from light at 2–8 °C. Do not freeze.
	Store MACSPlex Buffer, MACSPlex Filter Plate, and adhesive foils at room temperature. The expiration dates are indicated on the labels.

1.1 Principle of MACSPlex Cytokine Assays

MACSPlex Cytokine Assays are designed for determining concentrations of soluble analytes in a single sample. The analysis is based on MACSPlex Capture Beads, which display defined fluorescence properties and can be identified using standard flow cytometry techniques.

MACSPlex Capture Beads within the single MACSPlex Cytokine Reagents Kit contain a fluorescently labeled bead population, coated with a specific antibody reacting with one of the analytes within the sample. MACSPlex Capture Beads from different reagents kits display different fluorescence intensities.

Samples containing unknown levels of analytes are incubated with the antibody-coated MACSPlex Capture Beads, and analytes bind to the specific antibody. Each detection reagent contains APC-conjugated antibodies specific for one analyte and is also added. Consequently, sandwich complexes are formed between the MACSPlex Capture Bead, the analyte, and the detection reagent. These complexes can be analyzed based on the fluorescence characteristics of both the MACSPlex Capture Bead and the detection reagent. Standards of known quantities of given analytes are provided with the MACSPlex Cytokine Standards and are used for the quantification of the analytes within the unknown samples.

The analysis of MACSPlex Cytokine Assays require a flow cytometer with a blue (488 nm) and a red (635 nm) laser, which are capable of detecting FITC, PE, and APC. In combination with the Express Modes of the MACSQuant® Flow Cytometers the MACSPlex Cytokine Assays are optimized for automated measurement. They simplify flow cytometric analysis via predefined experiment settings as well as acquisition and analysis templates. They apply an automated gating

strategy to populations of interest that will be automatically adjusted for each data file individually to achieve optimal results.




Figure 1.1: Principle of MACSPlex Cytokine Assays.

1.2 Principle of the MACSPlex Cytokine Basic Kit

The MACSPlex Cytokine Basic Kit contains the supplementary components necessary to perform a MACSPlex Cytokine assay.

It has to be combined with the MACSPlex Cytokine 12 Standard, human and the desired MACSPlex Cytokine Reagents Kits, human or with the MACSPlex Cytokine 10 Standard, mouse and the desired MACSPlex Cytokine Reagents Kits, mouse, respectively.

Do not mix human and mouse specificities in one experiment.

MACSPlex Cytokine Reagents Kits are designed for detection of one or maximal up to seven analytes of choice in a single sample.

1.3 Applications

The MACSPlex Cytokine Basic Kit in combination with the MACSPlex Cytokine Reagents Kits and MACSPlex Cytokine Standard has been designed for flow cytometric detection of one or up to seven analytes of choice in a single sample. The kit has been optimized for use with serum, plasma, and cell culture supernatants.

1.4 Reagent and instrument requirements

- MACSPlex Cytokine 12 Standard, human (# 130-106-197) or MACSPlex Cytokine 10 Standard, mouse (# 130-106-198).

▲ **Note:** Do not mix human and mouse specificities in one experiment.

- MACSPlex Cytokine Reagents Kits, human or MACSPlex Cytokine Reagents Kits, mouse of choice.

▲ **Note:** Do not mix human and mouse specificities in one experiment.

- Polypropylene or polystyrene reagent tubes for serial dilutions of the MACSPlex Cytokine Standards as well as for preparation, dilution, and storage of unknown samples.
- Vacuum manifold or centrifuge with adapters to accommodate microtiter plates.
- Orbital shaker for microtiter plates or tubes (frequency 450–1400 rpm).
- MACSQuant® Analyzer, MACSQuant Analyzer 10 (# 130-096-343), or other flow cytometers equipped with blue (488 nm) and red (635 nm) lasers able to discriminate FITC, PE, and APC fluorescence.
▲ **Note:** The MACSQuant VYB cannot be used.
- MACS® Chill 96 Rack (# 130-094-459), when using the MACSQuant Analyzer or MACSQuant Analyzer 10.
- MACSQuant Calibration Beads (# 130-093-607), when using the MACSQuant Analyzer or MACSQuant Analyzer 10.
- (Optional for human assays) Cell culture medium, e.g., TexMACS™ Medium (# 130-097-196), DendriMACS™ GMP Medium (# 170-076-302), or RPMI 1640 supplemented with human AB serum.
- (Optional for mouse assays) Cell culture medium, e.g., RPMI 1640 supplemented with mouse serum or fetal bovine serum (FBS).
- Disposable pipette tips or pipettes.
- Universal microplate lid to cover the filter plate during washing steps in the centrifuge.

2. Protocols for assay performance

- ▲ Use of disposable pipette tips and/or pipettes is recommended to avoid any potential contamination or cross-contamination of reagents or samples, which might invalidate test results.

- ▲ Avoid air bubbles.

Prepare standards and samples and transfer to MACSPlex Filter Plate	20–45 min
Prepare MACSPlex Cytokine Capture Beads and transfer to MACSPlex Filter Plate	5–7 min
Incubate plate	2 hours
Wash plate (2x)	5 min
Prepare MACSPlex Cytokine Detection Reagent and transfer to MACSPlex Filter Plate	5–7 min
Incubate plate	1 hour
Wash plate	2 min
Resuspend samples with MACSPlex Buffer	3 min
Acquire data using the Express Mode of the MACSQuantify Software	

Figure 2.1: Experimental overview for the assay using a MACSPlex Filter Plate.

2.1 Preparation of the MACSPlex Cytokine Standard

- ▲ Reconstitute and dilute MACSPlex Cytokine Standard with MACSPlex Buffer, or use the same media as is used for the dilution of the unknown sample.
- ▲ Only use freshly prepared MACSPlex Cytokine Standard solutions. Do not store or reuse reconstituted or diluted standards.
- ▲ Use polypropylene or polystyrene reagent tubes. Do not use glass vials.

For the human MACSPlex Cytokine assay prepare the MACSPlex Cytokine 12 Standard, human as described in the data sheet.

For the mouse MACSPlex Cytokine assay prepare the MACSPlex Cytokine 10 Standard, mouse as described in the data sheet.

2.2 Sample preparation

- ▲ Handle all blood components and biological material as potentially hazardous.
- ▲ If unknown samples are expected or known to contain levels >2000 pg/mL, it is recommended to dilute the samples to make sure the fluorescence values are within the dynamic range of the standard curve.

- ▲ Use polypropylene or polystyrene reagent tubes. Do not use glass vials for sample preparation, dilution, or storage.

Preparation of serum samples

1. Allow the blood to clot for at least 30 minutes.
2. Centrifuge at 10,000×g for 10 minutes at 4 °C.
3. Transfer serum into a new tube and dilute at least 1:8 with MACSPlex Buffer, i.e., add 25 µL of the undiluted sample to 175 µL of MACSPlex Buffer.
4. Proceed to section 2.3.

Preparation of plasma samples

- ▲ Use freshly drawn blood samples supplemented with EDTA as anticoagulant. Do not store blood samples longer than 30 minutes before performing the assay.

1. Centrifuge peripheral blood at 10,000×g for 10 minutes at 4 °C.
2. Transfer plasma into a new tube and dilute at least 1:8 with MACSPlex Buffer, i.e., add 25 µL of the undiluted sample to 175 µL of MACSPlex Buffer.
3. Proceed to section 2.3.

Preparation of cell culture supernatant samples

1. Centrifuge cell culture supernatant at 10,000×g for 10 minutes at 4 °C.

2. Transfer the supernatant into a new tube.
3. (Optional) Dilute with cell culture medium or MACSPlex Buffer.
4. Proceed to section 2.3.

Frozen samples of serum, plasma, or cell culture supernatant


1. Thaw the samples completely and mix well by vortexing.
2. Centrifuge at 10,000×g for 10 minutes at 4 °C to remove particulates.
3. Transfer serum or plasma into a new tube and dilute at least 1:8 with MACSPlex Buffer, i.e., add 25 µL of the undiluted sample to 175 µL of MACSPlex Buffer. Cell culture supernatants can be diluted optionally with cell culture medium or MACSPlex Buffer.
4. Proceed to section 2.3.

2.3 MACSPlex Cytokine Assay

- ▲ Run the assay at room temperature. Work fast and keep samples protected from light, for example, cover plate or tubes with aluminum foil, especially during incubation steps.
- ▲ Unknown samples should be run in replicates, for example, in duplicates or triplicates and in different dilutions (e.g., undiluted, 1:3, 1:10, 1:30, 1:100) to make sure the fluorescence values are within the dynamic range of the standard curve.


2.3.1 Protocol for the assay using the MACSPlex Filter Plate

Design your assay using two columns of the MACSPlex Filter Plate for the standards. Add each of the standard samples in duplicates next to each other. Standards should be run in order from the lowest concentration (blank control: 0 pg/mL) to the highest concentration (stock solution: 10,000 pg/mL for human or 10,000–50,000 pg/mL for mouse standards). Start with the unknown sample in the next column of the plate. For details, see figure 2.2 for human or figure 2.3 for mouse assay.



Well position	Sample	Dilution
A1/A2	Blank control	
B1/B2	MACSPlex Cytokine 12 Standard (3.2 pg/mL)	1:3125 (Dilution5; 1:5 ⁵)
C1/C2	MACSPlex Cytokine 12 Standard (16 pg/mL)	1:625 (Dilution4; 1:5 ⁴)
D1/D2	MACSPlex Cytokine 12 Standard (80 pg/mL)	1:125 (Dilution3; 1:5 ³)
E1/E2	MACSPlex Cytokine 12 Standard (400 pg/mL)	1:25 (Dilution2; 1:5 ²)
F1/F2	MACSPlex Cytokine 12 Standard (2 ng/mL)	1:5 (Dilution1; 1:5 ¹)
G1/G2	MACSPlex Cytokine 12 Standard (10 ng/mL)	Stock solution (Dilution0; 1:5 ⁰)
H1/H2	Leave empty	
A3-H12	Add unknown samples	

Figure 2.2: Setup of the human assay with seven standard samples in duplicates using a 96-well plate.



Well position	Sample	Dilution
A1/A2	Blank control	
B1/B2	MACSPlex Cytokine 10 Standard (0.6 pg/mL or 3.2 pg/mL)	1:15,625 (Dilution6; 1:5 ⁶)
C1/C2	MACSPlex Cytokine 10 Standard (3.2 pg/mL or 16 pg/mL)	1:3125 (Dilution5; 1:5 ⁵)
D1/D2	MACSPlex Cytokine 10 Standard (16 pg/mL or 80 pg/mL)	1:625 (Dilution4; 1:5 ⁴)
E1/E2	MACSPlex Cytokine 10 Standard (80 pg/mL or 400 pg/mL)	1:125 (Dilution3; 1:5 ³)
F1/F2	MACSPlex Cytokine 10 Standard (400 pg/mL or 2 ng/mL)	1:25 (Dilution2; 1:5 ²)
G1/G2	MACSPlex Cytokine 10 Standard (2 ng/mL or 10 ng/mL)	1:5 (Dilution1; 1:5 ¹)
H1/H2	MACSPlex Cytokine 10 Standard (10 ng/mL or 50 ng/mL)	Stock solution (Dilution0; 1:5 ⁰)
A3-H12	Add unknown samples	

Figure 2.3: Setup of the mouse assay with eight standard samples in duplicates using a 96-well plate.

▲ Place the MACSPlex Filter Plate on a non-absorbent surface during loading steps and incubation, i.e. remove any tissues from the surface, to prevent the wells from running dry. Ensure that residual drops under the plate are completely removed to prevent liquid transfer, by placing the plate briefly on a tissue.

▲ Cover unused wells of the filter plate for later use with the adhesive foil provided with the kit.

▲ Washing steps are described for the use of a vacuum manifold. Alternatively, a centrifuge with an adapter for microtiter plates can be used: Put the MACSPlex Filter Plate on top of a conventional 96-flat-bottom microtiter plate with an universal lid and place both into the adapter. Centrifuge at 300xg for 3 minutes at room temperature.

1. Pre-wet required wells of the MACSPlex Filter Plate with 200 μ L of MACSPlex Buffer per well and aspirate off using a vacuum manifold designed to accommodate the filter plate (max. -300 mbar) until the wells are drained.
2. Place the filter plate briefly on a tissue to remove any residual liquid.
3. For generation of standard curve pipette 50 μ L of MACSPlex Buffer or media as a blank control, 50 μ L of each dilution, and the stock solution of the respective MACSPlex Cytokine Standard to the corresponding wells of the filter plate.
4. Add 50 μ L of each unknown sample per well.
5. Resuspend each MACSPlex Cytokine Capture Beads by vortexing for at least 30 seconds.
6. Calculate the amount of required MACSPlex Cytokine Capture Beads by determining the number of analytes and tests (standard

samples and test samples) in the experiment. Add additional 2 tests to the determined test number to ensure a sufficient volume of MACSPlex Cytokine Capture Beads.

7. Use 2 μL of each MACSPlex Cytokine Capture Beads per test and fill up to 20 μL with MACSPlex Buffer per test. For example:

	MACSPlex Cytokine Capture Beads	MACSPlex Buffer
1 analyte		
1 test	1 \times 2 μL	18 μL
20 tests	1 \times 40 μL	360 μL
4 analytes		
1 test	4 \times 2 μL	12 μL
20 tests	4 \times 40 μL	240 μL

8. Transfer 20 μL of diluted MACSPlex Cytokine Capture Beads to each well of standard samples and tests samples.
9. Incubate filter plate for 2 hours protect from light on an orbital shaker (450 rpm).
10. Apply the filter plate to the vacuum manifold and aspirate until wells are drained. Place the filter plate briefly on a tissue to remove any residual liquid.
11. Add 200 μL MACSPlex Buffer to each well and apply the filter plate to the vacuum manifold and aspirate off until wells are drained. Place the filter plate briefly on a tissue to remove residual liquid.
12. Repeat step 11.
13. Calculate the amount of required MACSPlex Detection Reagent by determining the number of analytes and tests (standard samples and test samples) in the experiment. Add additional 2

tests to the determined test number to ensure a sufficient volume of MACSPlex Detection Reagent.

14. Use 2 μL of each MACSPlex Detection Reagent per test and fill up to 100 μL with MACSPlex Buffer per test. For example:

	MACSPlex Cytokine Detection Reagent	MACSPlex Buffer
1 analyte		
1 test	1 \times 2 μL	98 μL
20 tests	1 \times 40 μL	1960 μL
4 analytes		
1 test	4 \times 2 μL	92 μL
20 tests	4 \times 40 μL	1840 μL

15. Transfer 100 μL of diluted MACSPlex Detection Reagent to each well of standard samples and tests samples.
16. Incubate filter plate for 1 hour protect from light on an orbital shaker (450 rpm).
17. Repeat wash steps 10 and 11.
18. Add 200 μL of MACSPlex Buffer to each well.
19. For sample acquisition using MACSQuant® Instruments and the Express Mode place the filter plate onto the Chill 96 Rack. To prevent liquid transfer from the wells, ensure that residual drops under the plate are completely removed by placing the plate briefly on a tissue.

▲ **Note:** Perform the flow cytometric acquisition on the same day, as prolonged storage of samples can result in increased background and reduced sensitivity.

▲ **Note:** Keep samples protected from light by using the protection lid during the flow cytometric acquisition with the MACSQuant Instrument.

2.3.2 Protocol for the assay using 1.5 mL reagent tubes

- ▲ Use polypropylene or polystyrene reagent tubes. Do not use glass vials.
- ▲ Standards should be run as duplicates. The order starts from the blank control (0 pg/mL) moving to the highest concentration (stock solution: 10,000 pg/mL for human or 10,000–50,000 pg/mL for mouse standards).

1. Label reagent tubes for the blank control, each dilution and the stock solution of the MACSPlex Cytokine Standard, and unknown samples.
2. Pipette 50 μL of MACSPlex Buffer or media as blank control, 50 μL of each dilution and the stock solution of the MACSPlex Cytokine Standard into the corresponding reagent tubes. Pipette 50 μL of each unknown sample into the corresponding reagent tube.
3. Resuspend each MACSPlex Cytokine Capture Beads by vortexing for at least 30 seconds.
4. Calculate the amount of required MACSPlex Cytokine Capture Beads by determining the number of analytes and tests (standard samples and test samples) in the experiment. Add additional 2 tests to the determined test number to ensure a sufficient volume of MACSPlex Cytokine Capture Beads.
5. Use 2 μL of each MACSPlex Cytokine Capture Beads per test and fill up to 20 μL with MACSPlex Buffer per test. For example:

	MACSPlex Cytokine Capture Beads	MACSPlex Buffer
1 analyte		
1 test	1 \times 2 μL	18 μL

20 tests	1 \times 40 μL	360 μL
4 analytes		
1 test	4 \times 2 μL	12 μL
20 tests	4 \times 40 μL	240 μL

6. Transfer 20 μL of diluted MACSPlex Cytokine Capture Beads to each tube of standard samples and tests samples.
7. Incubate for 2 hours protect from light on an orbital shaker (1400 rpm).
8. Add 0.5 mL of MACSPlex Buffer to each tube.
9. Centrifuge at 3000 \times g for 5 minutes.
10. Carefully aspirate off the supernatant using fresh pipette tips to avoid carry-over of samples. Leave 20 μL in the tube.
11. Resuspend the MACSPlex Cytokine Capture Bead pellet in each tube by adding 0.5 mL of MACSPlex Buffer and pipetting up and down.
12. Repeat steps 9 and 10.
13. Calculate the amount of required MACSPlex Detection Reagent by determining the number of analytes and tests (standard samples and test samples) in the experiment. Add additional 2 tests to the determined test number to ensure a sufficient volume of MACSPlex Detection Reagent.
14. Use 2 μL of each MACSPlex Detection Reagent per test and add to the remaining 20 μL of supernatant of step 10, fill up to 100 μL with MACSPlex Buffer per test.

15. For example:

	MACSplex Cytokine Detection Reagent	MACSplex Buffer
1 analyte		
1 test	1x2 µL	78 µL
20 tests	1x40 µL	1560 µL
4 analytes		
1 test	4x2 µL	72 µL
20 tests	4x40 µL	1440 µL

16. Transfer 80 µL of diluted MACSplex Detection Reagent to the remaining 20 µL of supernatant in each tube of standard samples and tests samples.
17. Incubate for 1 hour protect from light on an orbital shaker (1400 rpm).
18. Add 0.5 mL of MACSplex Buffer to each tube.
19. Centrifuge at 3000xg for 5 minutes.
20. Carefully aspirate off the supernatant using fresh pipette tips to avoid carry-over of samples. Leave 20 µL in the tube.
21. Resuspend each pellet in 0.5 mL of MACSplex Buffer by pipetting up and down.
22. Repeat steps 18 and 19.
23. Resuspend the pellet in each tube with 200 µL of MACSplex Buffer.
24. For sample acquisition with the MACSQuant® Express Mode, transfer samples to the MACSplex Filter Plate. Place the MACSplex Filter on a non-absorbent surface during loading steps, i.e., remove any tissues from the surface, to prevent the wells from running

dry. Place the filter plate onto a Chill 96 Rack and measure.

▲ **Note:** Acquire cytokine standards first, beginning with the standard samples of the first dilution series in order from the blank control to the highest concentration. Then process the standard samples of the second dilution series in the same order (see figure 2.2 for human or figure 2.3 for mouse assay). Afterwards acquire the unknown samples.

▲ **Note:** Perform the flow cytometric acquisition on the same day, as prolonged storage of samples can result in increased background and reduced sensitivity.

▲ **Note:** Keep samples protected from light by using the protection lid during the flow cytometric acquisition with the MACSQuant Instrument.

3. Flow cytometer set up

The kit includes MACSplex Setup Beads for flow cytometer set up. MACSplex Setup Beads are not required when using the MACSQuant Analyzer or MACSQuant Analyzer 10 but for all other instruments.

The kits is not suitable for use with the MACSQuant VYB.

3.1 Setup of the MACSQuant® Instrument

Calibrate the MACSQuant® Instrument using MACSQuant Calibration Beads (# 130-093-607). For details, refer to the data sheet of the MACSQuant Calibration Beads.

After successful completing the calibration, the MACSQuant Instrument is ready for measurement. No further steps are required as all necessary setup steps are performed automatically during calibration.

3.2 Setup of other flow cytometers

The analysis of MACSplex Cytokine Kit assays require a flow cytometer with a blue (e.g. 488 nm) and a red (e.g. 635 nm) laser, which are capable of detecting FITC, PE, and APC. For the purpose of setting up

these cytometers, MACSplex Setup Beads are included in the kit. For instructions on the setup procedures of other flow cytometers, please refer to the application note "General instructions for MACSplex Cytokine Kits" available on the product's webpage at www.miltenyibiotec.com/130-109-701.

4. Flow cytometric acquisition and data analysis using the MACSQuant® Express Mode

To perform the acquisition and data analysis of the MACSplex Custom Cytokine Assay with the MACSQuant® Instrument it is recommended to use the Express Modes "MACSplex_Standard" and "MACSplex_Sample" to achieve automated measurement and data analysis.

For details refer to the special protocol "Data acquisition and analysis of MACSplex Cytokine Kits using the MACSQuant Analyzer Express Modes" available at www.miltenyibiotec.com/130-109-701 under the Library tab. The minimum version number of the Express Mode package needed to run the assay on the MACSQuant Instrument is at least installExpressModes_28.1.15443 with MACSQuantify™ Software version 2.8 or higher.

To check the version number of your Express Mode package available on your MACSQuant Instrument please select **Help > Info** within the MACSQuantify Software (refer to figure 4.1). The version number of the Express Mode package is increasing with each Express Mode update. Make sure the MACSQuant Instrument contains an Express Mode package with at least the same or higher version number than the special protocol is marked with.

Overview

Relevant parameter for Express Mode experiments	human	mouse
Maximum number of standard positions (including blank) within the serial standard dilution	7	8
Dilution factor of serial standard dilution	5	5




Figure 4.1: Identification of version number of Express Mode package.

5. Troubleshooting

The following section offers solutions for problems that might be encountered when using MACSPlex Assays.

- **Variation between replicate samples:**
MACSPlex Cytokine Capture Beads can aggregate. Mix MACSPlex Cytokine Capture Beads at least 30 seconds before pipetting.
- **Low bead number in samples:**
Mix MACSPlex Cytokine Capture Beads sufficiently before pipetting.
Ensure that the instrument is calibrated for the relevant 96-well plate to avoid aspiration of air.
Avoid aspiration of beads during washing steps. Do not wash or resuspend beads in volumes higher than recommended.
- **High background:**
Avoid cross-well contamination. Ensure pipetting with multichannel pipettes and avoid touching reagent in the plate. The sample may be too concentrated. Test various sample dilutions.
The background may be due to non-specific binding. Increase number of washes to remove excess of MACSPlex Cytokine Detection Reagent.
- **Little or no detection of cytokine in sample:**
Sample may be too dilute. Test various sample dilutions. Use positive and negative control samples.

- **MACSPlex Cytokine Standard Samples show low fluorescence or result in poor standard curves:**
Check that all kit components are prepared and stored properly. Use a new vial of MACSPlex Cytokine Standard for each experiment. Adhere to incubation times as indicated in the protocol.
Assay was exposed to light. Keep the plate and beads covered with aluminum foil or a dark lid during all incubation steps. Use the protection lid during the sample acquisition with the MACSQuant® Instrument.
- **All samples are positive or above the high standard mean fluorescence value:**
The samples may be too concentrated. Dilute the samples further.
- **Beads not in region or gate:**
Ensure proper calibration of the MACSQuant Instrument. Samples containing organic solvents or samples of high viscosity should be diluted or dialyzed, respectively
- **Signal for whole plate is same as background:**
Incorrect amount or no MACSPlex Cytokine Detection Reagent was added.
- **Low signal for standard curve:**
Incubation steps were done at inappropriate temperatures, timings, or agitation.

- **Signals too high, standard curves are saturated:**
Incubation time was too long. Adhere to incubation times as indicated in the protocol.
- **Sample readings are out of range – if below detection limit:**
Samples may be too diluted. Adjust dilutions accordingly.
Samples may contain no analyte. Analyte concentration may be below the detection limit. It may be possible to use higher sample volumes. Please contact Milteny Biotec's Technical Support for appropriate protocol modifications.

If above detection limit:
Samples contain analyte concentrations above 2 ng/mL (mouse) or 10 ng/mL (human). Dilute the samples further.
- **High variation in samples and/or standards:**
Multichannel pipette may not be calibrated.
Plate washing was not uniform.
Samples may have contained high particulate matter or other interfering substances.
Plate agitation was insufficient.
Cross-well contamination. Ensure pipetting with multichannel pipettes and avoid touching reagent in the plate. Change pipette tips for each well when touching the reagent.
- **Filter plate will not vacuum:**
Vacuum pressure is insufficient. Increase vacuum pressure. Samples contain insoluble particles. Centrifuge samples just prior to performing the assay and use supernatant.
High lipid concentration. Centrifuge samples, remove lipid layer, and use supernatant to perform the assay.
- **Plate leakage:**
Vacuum pressure is too high. Adjust vacuum pressure to maximal –300 mbar.
Place the MACSPlex Filter Plate on a non-absorbent surface during filling steps and incubation, i.e., remove any tissues from the surface, to prevent the wells from running dry.
Ensure that residual drops under the plate are completely removed to prevent leakiness of the wells, by placing the plate briefly on a tissue, after each washing step.
Avoid touching the plate filter with the tip of the pipette when adding reagents to the wells. Refer to www.miltenyibiotec.com for all data sheets and protocols. Milteny Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Milteny Biotec contact.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Warnings

Reagents contain sodium azide. Under acidic conditions sodium azide yields hydrazoic acid, which is extremely toxic. Azide compounds should be diluted with running water before discarding. These precautions are recommended to avoid deposits in plumbing where explosive conditions may develop.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein. All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

DendriMACS, MACS, MACSQuant, MACSQuantify, the Miltenyi Biotec logo, and TexMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2021 Miltenyi Biotec and/or its affiliates. All rights reserved.