

Contents

1. Description
 - 1.1 Background information
 - 1.2 Applications
2. Protocol

1. Description

This product is for research use only.

Products StemMACS mRNA Transfection Kits

Order no.	Components
130-132-978	125 µL StemMACS Transfection Reagent 25 mL StemMACS Transfection Buffer
130-132-949	250 µL StemMACS Transfection Reagent 25 mL StemMACS Transfection Buffer
130-104-463	0.5 mL StemMACS Transfection Reagent 25 mL StemMACS Transfection Buffer

Quality control Functionally tested for transfection of human fibroblasts.

Storage Store all components at 2–8°C. Do not freeze. The expiration date is indicated on the vial label.

1.1 Background information

The transient, non-integrative expression of key developmental regulators, recombinases, or markers via mRNA transfection is a powerful tool for modulating cell fate. The StemMACS mRNA Transfection Kit is a novel lipid-based transfection system that has been designed for efficient mRNA delivery into various cell types. StemMACS mRNA Transfection Reagent has been selected for its minimal cytotoxicity and high transfection efficiency. It is therefore particularly suited for sensitive cell lines such as pluripotent stem cells and complex transfection schedules that involve repeated mRNA delivery over several days. The StemMACS mRNA Transfection Kit has been successfully used for transfection of primary human fibroblasts and the generation of iPS cell lines by mRNA reprogramming.

For successful transfection, use a protocol that is optimized for the cell type of interest. StemMACS eGFP mRNA or StemMACS Nuclear eGFP mRNA allow an easy read-out for transfection efficiency and are recommended as positive controls when establishing new transfection protocols. Critical parameters for successful mRNA transfection include mRNA quality and effective suppression of the innate antiviral response to single-stranded RNA. It is recommended to use only high quality mRNA

preparations that have been validated for mRNA transfection and incorporate modified nucleotides such as pseudouridine and 5-methylcytidine (e.g. StemMACS mRNA products). Depending on the application and cell line, supplementation of the cell culture medium with B18R protein may be necessary to ensure complete suppression of the innate antiviral response during transfection.

1.2 Applications

- Transient, non-integrative delivery of mRNA encoded factors into a broad range of cell types, including primary fibroblasts and human iPS cells
- mRNA reprogramming of human fibroblasts
- mRNA-induced differentiation of stem and progenitor cells
- mRNA-induced transdifferentiation of differentiated cells
- mRNA-induced recombination
- Transient labeling with fluorescent proteins

2. Protocol

▲ RNA is susceptible to degradation by exogenous ribonucleases. Wear gloves and use RNase-free tubes, reagents, and pipettes. It is recommended to use low RNA/DNA-binding tubes and barrier pipette tips when handling mRNA.

1. The day before transfection, seed cells in 24-well or 6-well tissue culture plates. For recommended seeding conditions refer to table 1.

▲ **Note:** The transfection can be carried out in the presence of antibiotics (e.g. penicillin, streptomycin) and serum. A medium change prior to transfection is not required.

Table 1: Recommended seeding conditions

	per 24-well plate	per 6-well plate
Recommended seeding density	1–5×10 ⁴ cells	10 ⁴ –10 ⁵ cells
Cell culture medium	475 µL	1400 µL
mRNA transfection complex (added the following day)	25 µL	100 µL
Total volume	500 µL	1500 µL

2. On the day of transfection, thaw mRNA aliquots slowly on ice. Prepare two sterile, RNase-free microfuge tubes for each well to be transfected.
3. Dilute the mRNA in StemMACS Transfection Buffer. Mix by pipetting up and down.

▲ **Note:** The amount of mRNA needed may vary with the cell type, the protein encoded by the mRNA, and the specific requirements of the downstream application. Initially, it is recommended to test at least three different mRNA concentrations in the range of 0.1 to 3 µg mRNA to find the optimal amount for transfection. Refer to table 2 for guidelines on setting up the first experiment.

- Dilute StemMACS Transfection Reagent in StemMACS Transfection Buffer. Mix by pipetting up and down.
Note: In most cases, using a ratio of 1:3 (w/v) of mRNA to StemMACS Transfection Reagent will be appropriate, i.e., use 3 µL of StemMACS Transfection Reagent for each 1 µg of mRNA to be transfected. The volumes given in table 2 are based on this recommended ratio. In case further optimization is needed this ratio may be varied in the range between 1:2 to 1:4 (w/v).

Table 2: Guidelines for preparation of the mRNA transfection complex

24-well format

Tube 1		Tube 2	
mRNA	StemMACS Transfection Buffer	StemMACS Transfection Reagent	StemMACS Transfection Buffer
150 ng	Adjust to 12.5 µL with buffer	0.45 µL	12.15 µL
250 ng	Adjust to 12.5 µL with buffer	0.75 µL	11.75 µL
500 ng	Adjust to 12.5 µL with buffer	1.50 µL	11.50 µL
Total volume: 12.5 µL		Total volume: 12.5 µL	
25 µL mRNA transfection complex			

6-well format

Tube 1		Tube 2	
mRNA	StemMACS Transfection Buffer	StemMACS Transfection Reagent	StemMACS Transfection Buffer
250 ng	Adjust to 50 µL with buffer	0.75 µL	49.25 µL
500 ng	Adjust to 50 µL with buffer	1.50 µL	48.50 µL
1000 ng	Adjust to 50 µL with buffer	3.00 µL	47.00 µL
Total volume: 50 µL		Total volume: 50 µL	
100 µL mRNA transfection complex			

- Combine both dilutions by adding the diluted StemMACS Transfection Reagent to the diluted mRNA.
- Mix by pipetting up and down. Do not vortex.
- Incubate at room temperature for 20 minutes.
- Add the mRNA transfection complex dropwise to the pre-plated cells (i.e. 25 µL of transfection complex per 24-well format, 100 µL per 6-well format). Keep the pipet tip above the surface of the media.
- Gently rock the plates to ensure even distribution of the mRNA transfection complex.
- Incubate cells under standard cell culture conditions.

- After 4 hours, replace the cell culture media with fresh, pre-warmed medium.
Note: This step is highly recommended, when working with sensitive cell lines or performing repeated daily transfections to minimize potential toxicity. The media replacement may be omitted when working with robust cell lines or single transfections.
- Expression of the transfected factor will be detectable 6–24 hours after transfection depending on mRNA and protein stability.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

The Miltenyi Biotec logo and StemMACS are registered trademarks or trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2023 Miltenyi Biotec and/or its affiliates. All rights reserved.