Miltényi Biotec

## CliniMACS® Accessories

Transfer Set Coupler/Coupler  
Luer/Spike Interconnector

**REF** 200-073-902

**QTY** 5

**REF** 200-073-903

**QTY** 5

**Miltényi Biotec B.V. & Co. KG**

Friedrich-Ebert-Straße 68

51429 Bergisch Gladbach | Germany

☎ +49 2204 8306-3803

CE 0123

## Table of contents

|    |                                       |    |
|----|---------------------------------------|----|
| en | Instructions for use . . . . .        | 2  |
| ar | تعليمات الاستخدام . . . . .           | 4  |
| bg | Инструкции за употреба . . . . .      | 6  |
| cs | Návod k použití . . . . .             | 8  |
| da | Brugsanvisning. . . . .               | 10 |
| de | Gebrauchsanweisung . . . . .          | 12 |
| el | Οδηγίες χρήσης . . . . .              | 14 |
| es | Instrucciones de uso . . . . .        | 16 |
| et | Kasutusjuhend . . . . .               | 18 |
| fi | Käyttöohje . . . . .                  | 20 |
| fr | Notice d'utilisation . . . . .        | 22 |
| he | הוראות שימוש . . . . .                | 24 |
| hr | Upute za upotrebu . . . . .           | 26 |
| hu | Használati útmutató . . . . .         | 28 |
| it | Istruzioni per l'uso. . . . .         | 30 |
| lt | Naudojimo instrukcija . . . . .       | 32 |
| lv | Lietošanas pamācība . . . . .         | 34 |
| nl | Gebruiksaanwijzing. . . . .           | 36 |
| no | Bruksanvisning. . . . .               | 38 |
| pl | Instrukcja używania. . . . .          | 40 |
| pt | Instruções de utilização . . . . .    | 42 |
| ro | Instrucțiuni de utilizare. . . . .    | 44 |
| ru | Инструкции по использованию . . . . . | 46 |
| sk | Návod na použitie. . . . .            | 48 |
| sl | Navodila za uporabo . . . . .         | 50 |
| sr | Uputstvo za upotrebu . . . . .        | 52 |
| sv | Bruksanvisning. . . . .               | 54 |
| tr | Kullanma talimatları. . . . .         | 56 |

# Instructions for use

## Contents

**QTY** Five Transfer Set Coupler/Coupler or  
Five Luer/Spike Interconnector

## Intended purpose

The CliniMACS® Accessories are intended for use in combination with CliniMACS Tubing Sets and CliniMACS Prodigy® Tubing Sets:

The Transfer Set Coupler/Coupler is intended for interconnecting different containers or for connecting containers with the CliniMACS Tubing Sets or the CliniMACS Prodigy Tubing Sets to channel solutions or human cell suspensions to assist the enrichment or depletion of human cells with a CliniMACS System only.

The Luer/Spike Interconnector is intended for interconnecting different containers or for connecting containers with the CliniMACS Tubing Sets or the CliniMACS Prodigy Tubing Sets to channel solutions or human cell suspensions or for sample collection to assist the enrichment or depletion of human cells with a CliniMACS System only.

## Limitation

Miltenyi Biotec as the manufacturer of the CliniMACS System does not give any recommendations regarding the use of separated cells for therapeutic purposes and does not make any claims regarding a clinical benefit.

## Side-effects

Side-effects have not been reported.

## Precautions

The corresponding CliniMACS User Manual must be followed.

Aseptic working procedures must be applied for the unpacking, assembly, and use of the product.

After installation of the product check all luer connections for tightness. If necessary retighten the luer connections.

For the manufacturing and use of target cells in humans the national legislation and regulations – e.g. for the EU the Directive 2004/23/EC (human tissues and cells) or the Directive 2002/98/EC (human blood and blood components) – must be followed. Thus, any clinical application of the target cells is exclusively within the responsibility of the user of a CliniMACS System.

The separation of cells using the product must be performed by trained operators only.

All materials that have been in contact with blood and blood components must be treated and disposed as biohazardous materials. Therefore, standard hospital or institutional requirements and national legislations must be observed.

## Warnings

Do not use after the use-by date printed on the product label.

Do not use if package is damaged. Before opening, inspect the packaging for damage, punctures, or tears. Use the product only if package is undamaged and sealed.

Do not use if the product is damaged. Do not use if any leakages of the product are observed during priming or separation.

Do not connect the product directly to the patient.

Do not store blood, blood fractions or cell fractions in the product.

Do not re-use. The re-use of the product or its parts leads to the endangerment of the patient due to biological contamination and to inefficient cell separation.

## Storage

The product must be stored at room temperature (+15 °C to +30 °C [+59 °F to +86 °F]).

## Further information

The product is sterilized using irradiation.

The depicted symbol indicates the single sterile barrier of the product with protective packaging outside.

Any serious incident that has occurred in relation to this product should be reported to Miltenyi Biotec B.V. & Co.KG - using the contact information provided - and the competent authority of the member state in which the user of this product is established.

The product contains a non-pyrogenic, sterile fluid path.

## Performance

The performance of the product depends on the chosen application. For information on the different applications, refer to the corresponding CliniMACS User Manual.

## Handling instructions

Refer to the CliniMACS User Manual of the corresponding application for the validated combinations of instrument, separation programs, reagents, tubing sets, and buffer.

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Miltenyi Biotec Swiss AG
Miltenyi Biotec Ltd.