Miltenyi Biotec

CryoMACS® Freezing Bags

Instructions for use

CryoMACS® Freezing Bags

Instructions for use

Issued: 2021-04
35210/06

CE 0123

For U.S. CAUTION:
Federal (USA) law restricts this device
to sale by or on the order of a physician.

Miltenyi Biotec B.V. & Co. KG

Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany

Miltenyi Biotec Technical Support:

☎ +49 2204 8306-8484
☎ +49 2204 8306-80 (hotline)
☎ +49 2204 8306-89
✉ macsde@miltenyi.com
🏠 www.miltenyibiotec.com

Table of contents

en	Instructions for use	4
ar	إرشادات الاستخدام	7
cs	Návod k použití	8
da	Brugsanvisning	10
de	Gebrauchsanweisung	12
el	Οδηγίες χρήσης	14
es	Instrucciones de uso	16
fi	Käyttöohje	18
fr	Notice d'utilisation	20
he	הוראות שימוש	23
hr	Upute za upotrebu	24
hu	Használati útmutató	26
it	Istruzioni per l'uso	28
nl	Gebruiksaanwijzing	30
no	Bruksanvisning	32
pl	Instrukcja używania	34
pt	Instruções de utilização	36
ro	Instrucțiuni de utilizare	38
ru	Инструкции по использованию	40
sk	Návod na použitie	42
sl	Navodila za uporabo	44
sv	Bruksanvisning	46
tr	Kullanma talimatları	48

Contents

Each CryoMACS® Freezing Bag consists of one freezing bag with access ports as the primary containment for hematopoietic progenitor cells and one overwrap bag as a secondary containment.

Part number (P/N)	Catalogue number (REF)	Product	Nominal volume	Recommended fill volume
74400	200-074-400	CryoMACS Freezing Bag 50	50 mL	10–20 mL
74401	200-074-401	CryoMACS Freezing Bag 250	250 mL	30–70 mL
74402	200-074-402	CryoMACS Freezing Bag 500	500 mL	55–100 mL
74403	200-074-403	CryoMACS Freezing Bag 750	750 mL	80–190 mL
74404	200-074-404	CryoMACS Freezing Bag 1000	1000 mL	125–270 mL

The fill volume recommendations are based on freezing in a horizontal position with the use of a protective metal cassette with an internal thickness of 4 mm. If any other storage form is used, the optimal fill volume and freezing conditions must be validated by the user.

Intended purpose

The CryoMACS Freezing Bags are intended for a single cycle of freezing, storage (down to –196 °C [–321 °F]), and subsequent thawing (at +37 °C [+99 °F]) of hematopoietic progenitor cells.

Precautions

All materials that have been in contact with blood or blood components, including hematopoietic progenitor cells, must be treated and disposed as biohazardous materials. Therefore, standard hospital or institutional requirements and national legislations must be observed.

Before application to humans, the suitability of the stored hematopoietic progenitor cells must be demonstrated regarding the indication, quality, and quantity. For the use of the stored hematopoietic progenitor cells in humans the national legislation and regulations must be followed.

Any use of disinfectants on the freezing bag and the overwrap bag of the CryoMACS Freezing Bag is exclusively within the responsibility of the user and risks must be assessed within the user's own risk management.

Full performance of this product requires correct operation by the user after filling of the product according to the instructions for use. Pay special attention to use of the overwrap bag and metal cassettes as well as drying of the product and removal of residual air.

The user of the CryoMACS Freezing Bag is responsible for the selection of adequate cryopreservation protocols for the target components. Cryo-preservation using these bags should be performed by trained operators only, such as healthcare professionals in a clinical environment. Thus,

CryoMACS Freezing Bags are not intended to be used by lay-persons. Handle frozen bags with care.

Warnings

Do not use after the use-by date indicated on the product.

Do not use if package is damaged. Use product only if package is undamaged and sealed and if port and tubing closures are in place. Use only intact bags without any signs of leakage. Use only the provided labeling pouch for identification of contents. Use only permanent sealing for tubing closure. Administer DMSO only at a concentration up to 20% (v/v) to the bag. DMSO can cause degradation of plastic materials.

Do not re-use. Re-use causes contamination of the interior fluid path. The product cannot be sterilized after use. The fluid path can be compromised by the thermal stress of re-use and the irreversible breach of the bag closures.

Storage

Store the product at +5 °C to +35 °C (+41 °F to +95 °F) until usage.

Keep away from sunlight.

Further information

Sterilized using irradiation.

The sterile barrier (SB) symbol indicates the sterile packaging of the product.

Any serious incident that has occurred in relation to this product should be reported to Miltenyi Biotec B.V. & Co. KG – using the contact information provided – and the competent authority of the member state in which the user of this product is established.

Handling instructions

Follow these handling instructions for the product and institutional protocols for filling, freezing (including freezing parameters), storage, and thawing. Follow aseptic working procedures for fluid transfer.

Labeling: Place an identification tag into the label pocket of the freezing bag and spot seal intermittently along the edge. Do not seal the pocket completely.

Fluid transfer (in): Use only the filling assembly for transfer of fluids into the bag. Do not use the spike ports. Use only the injection port for transfer of fluids containing DMSO. Attach transfer devices to either of the luer connectors, the injection port or by sterile docking to the PVC tubing (inner/outer diameter 2.6/4.1 mm) distal to the injection port. Observe the instructions for use supplied by the manufacturer of the sterile docking device. Carefully remove any residual air from the freezing bag after filling.

WARNING

The following sealing procedures are performed on ethyl vinyl acetate (EVA) material, adjust dielectric or radio frequency (RF) sealing procedures accordingly.

After filling, close the bag by sealing the EVA tubing. Test the seal integrity and remove the tubing distal to the seal.

Freezing: Ensure during this procedure that all surfaces of the bags and the cassette are completely dry before freezing, because ice particles on the surface may damage the bags. Place the freezing bag into the overwrap bag and vacuum seal. Place the doubled bag into a protective metal cassette. Avoid bending, folding, or kinks of the bag when closing the cassette. Seal off the tubing proximal to the freezing bag before freezing. Use of a controlled rate freezer is recommended.

Storage: Storage should be performed in the vapor phase of liquid nitrogen. Temperature fluctuations during storage should be avoided.

Thawing: If the bag is stored in the liquid phase of nitrogen, place the bag into the vapor phase of nitrogen for a minimum of 4 hours before removal from the nitrogen tank and thawing. Remove the bag with its overwrap from the cassette and inspect for leaks or ruptures before thawing in a water bath at +37 °C (+99 °F).

WARNING

Entry of air or nitrogen into the bag during storage may cause rapid expansion and subsequent bursting of the bag during rewarming.

Fluid transfer (out): Completely dry the bag surface, remove the overwrap bag, twist off the top part of one of the spike ports and carefully access the bag contents by puncturing the port seal with the spike of a sterile transfusion assembly.

Description of components

The components of a CryoMACS Freezing Bag and the filling assembly are described in the drawing shown on the fold-out page.

No.	CryoMACS Freezing Bag components
1	Spike ports
2	Sealing area
3	EVA tubing
4	Label pouch
5	Overwrap bag
No.	Filling assembly components
6	Injection port
7	PVC tubing
8	Roller clamp
9	Male luer lock connector
10	Female luer lock connectors

Glossary of symbols

An explanation of the symbols used for the labeling of the products can be found in the glossary on page 50.

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Miltenyi Biotec Swiss AG
Miltenyi Biotec Ltd.