



Special protocol

Index

1. Description
 - 1.1 Background
 - 1.2 Reagent and instrument requirements
 - 1.3 Reagent preparation
2. Protocol for immunohistochemical staining of paraffin-embedded tissue sections
3. Example of immunohistochemical staining for CD133
4. References

1. Description

1.1 Background

CD133 (also known as AC133) is a five-transmembrane domain glycoprotein that is selectively expressed on populations of hematopoietic stem and progenitor cells deriving from adult and fetal bone marrow, cord blood, and peripheral blood.¹ In addition, CD133 is known to be a marker of stem cells of a variety of nonhematopoietic tissues, including neural and glial cells in the fetal brain² and neural progenitor cells in human skin³. Furthermore, CD133 was found on stem cells of prostatic epithelia⁴, muscle tissue⁵, kidney⁶, liver⁷, and corneal stroma⁸. CD133 has also been found to be expressed on cancer stem cells of a variety of cellular and tissue origin.⁹ In human fetal tissue, CD133 has been found to be expressed on the neural tube, gut, and kidney¹⁰ and highly in fetal liver – the site of embryonic hematopoiesis.

1.2 Reagent and instrument requirements

- 100% ethanol
- 100% methanol
- ChemMate™ Target Retrieval Solution (10×) (DakoCytomation, # S 2031)
- 30% (w/w) Hydrogen peroxide solution (H₂O₂)
- CD133/1 (AC133) pure, human (Miltenyi Biotec, # 130-090-422)
- VECTASTAIN® ABC kit (Mouse IgG) (Vector Laboratories, # PK-4002)
 - ▲ **Note:** The VECTASTAIN ABC kit contains the following reagents: blocking serum, biotinylated affinity-purified anti-mouse immunoglobulin antibody, Reagent A (Avidin DH solution), and Reagent B (biotinylated horseradish peroxidase A).
- 3,3'-Diaminobenzidine tetrahydrochloride (DAB) (Sigma, # D5637)
- Nickel(II) chloride (NiCl₂)
- Tris base
- Sodium chloride (NaCl)
- Hydrochloric acid (HCl)
- Methyl green solution or Meyer's hemalum
- Phosphate-buffered saline (PBS)
- Deionized water

Immunohistochemical detection of CD133 human

- Tap water
- H₂O₂
- Water bath
- Microwave
- Fluorescence microscope
- Hydrophobic pen
- Cover slides of appropriate length
- Hellendahl jars
- Plastic Hellendahl jars
- Humidified chamber
- Roti®-Histol (Roth, # 6640.1)
- Roti®-Histokitt (Roth, # 6638.1)

1.3 Reagent preparation

Prepare the following stock solutions:

- Ethanol dilution series: 100%, 96%, 80%, 70% ethanol in deionized water. Store at room temperature.
- 10× Tris buffer: Dissolve 61 g of Tris base and 116.9 g of NaCl in 1 L of deionized water. Adjust pH to 7.6 (±0.1) with HCl. Store at room temperature.
- DAB stock solution: Dissolve 5 g of DAB tetrahydrochloride in 132 mL of 1× Tris buffer. Prepare 4 mL aliquots and store at -20 °C.

Prepare the following working solutions:

▲ **Note:** Always prepare fresh working solutions.

- Target Retrieval Solution: Dilute 10× Target Retrieval Solution 1:10 with deionized water.
- 3% H₂O₂: Dilute 30% H₂O₂ solution 1:10 with 100% methanol.
- VECTASTAIN ABC kit: Prepare reagents according to the manufacturer's instructions.
- NiCl₂ working solution: Dissolve 8 g of NiCl₂ in 100 mL of deionized water. Store at 4 °C.
 - ▲ **Note:** Caution! NiCl₂ is toxic.
- DAB tetrahydrochloride working solution: Dilute 4 mL of the DAB stock solution with 175 mL of 1× Tris buffer. Warm to 37 °C in a water bath.

2. Protocol for immunohistochemical staining of paraffin-embedded tissue sections


1. Prepare sections of Paraffin-embedded tissues of no more than 2 µm in thickness.



2. Deparaffinize and rehydrate the tissue sections by serial immersion in Hellendahl jars containing the following:
 - 2×5 min Roti-Histol
 - 2×2 min 100% ethanol
 - 2×2 min 96% ethanol
 - 2×2 min 80% ethanol
 - 2×2 min 70% ethanol
 - 1×5 min PBS
 - 1×brief rinse with deionized water
3. For antigen retrieval, fill a plastic Hellendahl jar with 1× Target Retrieval Solution and submerge the tissue sections completely in the solution. Heat uncovered at 250–300 W in a microwave for 15 minutes. Replace any evaporated Target Retrieval Solution and leave to cool to room temperature.
4. Wash 1×5 minutes with deionized water, followed by 1×5 minutes with PBS in a Hellendahl jar.
5. Block endogenous peroxidase activity by immersing slide in 3% H₂O₂ solution for 30 min at room temperature.
 - ▲ **Note:** For tissue sections with low endogenous peroxide activity, e.g. organ tissue, the incubation time can be shortened to 10 minutes.
6. Wash 2×3 minutes with fresh PBS in a Hellendahl jar.
7. Dab slide dry around the section and encircle it with a hydrophobic pen.
8. Block non-specific binding of CD133/1 (AC133) pure antibody by incubating slide in a humidified chamber at room temperature for 20 minutes in the VECTASTAIN ABC kit blocking serum.
9. Drain off supernatant by tilting the slide and then blot dry. Do not wash the slide.
10. Pipette 50–250 µL of CD133/1 (AC133) pure antibody, diluted in the VECTASTAIN ABC kit blocking serum, to each section and incubate slide for 1 hour at room temperature in a humidified chamber.
 - ▲ **Note:** The optimal working concentration of antibody should be pre-determined as it will vary according to tissue type and thickness of sections.
 - ▲ **Note:** For a negative control, use VECTASTAIN ABC kit blocking serum only.
11. Wash 2×3 minutes with fresh PBS in a Hellendahl jar.
12. Pipette carefully 50–250 µL of the VECTASTAIN ABC kit biotinylated secondary antibody, diluted in the blocking serum, to each section and incubate slide for 30 minutes at room temperature in a humidified chamber.
 - ▲ **Note:** During this incubation step, prepare the VECTASTAIN ABC reagent according to the manufacturer's instructions. This reagent must be incubated at room temperature for 30 minutes before use.
13. Wash 2×3 minutes with fresh PBS in a Hellendahl jar.
14. Pipette carefully 50–250 µL of freshly prepared VECTASTAIN ABC reagent to each section. Incubate for 30 minutes at room temperature in a humidified chamber.
15. Wash 2×3 minutes with fresh PBS in a Hellendahl jar.
16. Add 0.1 mL of the 3% H₂O₂ working solution and 1 mL of the NiCl₂ working solution to the pre-warmed DAB tetrahydrochloride working solution.
 - ▲ **Note:** Caution! NiCl₂ is toxic.
17. Incubate tissue sections in the DAB working solution in a Hellendahl jar for 10 minutes at 37 °C in the water bath.
18. Wash 2×3 minutes with tap water.
19. Counterstain sections by incubating for 1–2 minutes at room temperature in Methyl Green solution or Meyer's hemalum. Shake occasionally.
20. Dehydrate the tissue sections by serial immersion in Hellendahl jars containing the following:
 - 2×brief immersion in 96% ethanol
 - 3×brief immersion in 100% ethanol
 - ▲ **Note:** Last ethanol step must be clear after immersion of the tissue sections.
 - 2×5 min Roti-Histol
21. Mount sections directly using the Roti-Histokitt and apply coverslip.


3. Example immunohistochemical staining for CD133

A. Paraffin-embedded tissue sections of a human Pineoblastoma were stained with CD133/1 (AC133) pure and developed using DAB as the substrate (brown color). Hemalum was used for counterstaining.



Courtesy of Dr. Michel Mittelbronn, University of Tübingen, Germany.

B. Paraffin-embedded tissue sections of a human Pineoblastoma were stained with CD133/1 (AC133) pure and developed using the VECTASTAIN ABC kit with DAB as the substrate (brown color). Methyl Green was used for counterstaining.



4. References

1. Yin, AH. *et al.* (1997) AC133, a novel marker for human hematopoietic stem cells and progenitor cells. *Blood* 90: 5002–5012.
2. Yu, S. *et al.* (2004) Isolation and characterization of the CD133⁺ precursors from the ventricular zone of human fetal brain by magnetic affinity cell sorting. *Biotechnol. Lett.* 26: 1131–1136.
3. Belicchi, M. *et al.* (2004) Human skin-derived stem cells migrate throughout forebrain and differentiate into astrocytes after injection into adult mouse brain. *J. Neurosci. Res.* 77: 475–486.
4. Richardson, GD. *et al.* (2004) CD133, a novel marker for human prostatic epithelial stem cells. *J. Cell Sci.* 117: 3539.
5. Alessandri, G. *et al.* (2004) Isolation and culture of human muscle-derived stem cells able to differentiate into myogenic and neurogenic cell lineages. *Lancet* 364:1872–1883.
6. Bussolati, B. *et al.* (2005) Isolation of renal progenitor cells from adult human kidney. *Am. J. Pathol.* 166: 545–555.
7. Laurson, J. *et al.* (2005) Hepatocyte progenitors in man and in rodents – multiple pathways, multiple candidates. *Int. J. Exp. Pathol.* 86: 1–18.
8. Thill, M. *et al.* (2004) Identification of a population of CD133⁺ precursor cells in the stroma of human cornea. MEETING ABSTRACT
9. Collins, A. *et al.* (2005) Prospective identification of tumorigenic prostate cancer stem cells. *Cancer Res.* 65: 10946–10951.
10. Corbeil, D. *et al.* (2000) The human AC133 hematopoietic stem cell antigen is also expressed in epithelial cells and targeted to plasma membrane protrusions. *J. Biol. Chem.* 275: 5512–5520.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS is a registered trademark of Miltenyi Biotec GmbH.

Roti is a registered trademark of Carl Roth GmbH.

ChemMate is a trademark of Ventana Medical Systems, Inc.

VECTASTAIN is a registered trademark of Vector Laboratories.

© 2006 Miltenyi Biotec GmbH. Printed in Germany.