



Neural Tissue Dissociation Kit Enzyme Mix 2

Order no. 130-094-160

This product can be used as supplement with the Neural Tissue Dissociation Kit (P) (# 130-092-628) or the Neural Tissue Dissociation Kit (T) (# 130-093-231) for the repetition of incubation with Enzyme Mix 2.

1. Description

Components 3 vials containing:
1.5 mL of Solution 3 (sterile)
5–15 mg* of Solution 4
1 mL of Storage Buffer for Solution 4

Storage Store all components at 2–8 °C. Prior to first usage, Storage Buffer for Solution 4 should be added to Solution 4, sterile filtered, aliquoted and stored at –20 °C. Aliquots are stable for 4 weeks at 2–8 °C. The expiration date is indicated on the kit box label.

*A variability in ranges of mg of powder contained within the vial is possible—according to variability of powder activity.

2. Preparation

For details on usage, please refer to the data sheet of the respective Neural Tissue Dissociation Kit.

All protocols and data sheets are available at www.miltenyibiotec.com.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

MACS is a registered trademark of Miltenyi Biotec GmbH.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

Copyright © 2009 Miltenyi Biotec GmbH. All rights reserved.

140-002-477-02