



Column Extensions

Order no. 130-041-419

Index

1. Description
 - 1.1 Example applications
2. Instructions for Assembly

1. Description

Components 100 Column Extensions for optimal positioning of the CS Columns in the magnetic field of a SuperMACS™ Separator or SuperMACS II Separator, sterile packed.

Storage Store dry.

1.1 Example application

- Magnetic separation with CS Columns in SuperMACS Separator or SuperMACS II Separator.

2. Instructions for Assembly

1. Remove yellow cap from the bottom of the CS Column and replace by the column extension. Attach tightly to the bottom fitting of the CS Column.
2. Attach CS Column with extension to port A of the 3-way-stop cock.
3. For further assembly and filling of the column refer to CS Column data sheet.

Warranty

The products sold hereunder are warranted only to be free from defects in workmanship and material at the time of delivery to the customer. Miltenyi Biotec GmbH makes no warranty or representation, either expressed or implied, with respect to the fitness of a product for a particular purpose. There are no warranties, expressed or implied, which extend beyond the technical specifications of the products. Miltenyi Biotec GmbH's liability is limited to either replacement of the products or refund of the purchase price. Miltenyi Biotec GmbH is not liable for any property damage, personal injury or economic loss caused by the product.

Unless otherwise specifically indicated, all Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.

MACS is a registered trademark of Miltenyi Biotec GmbH.
SuperMACS is a trademark of Miltenyi Biotec GmbH.

© 2006 Miltenyi Biotec GmbH. Printed in Germany.

Assembly of the CS Column for use in the SuperMACS™ Separator or SuperMACS II Separator.

