



Miltenyi Biotec

Antibody cross-reactivity information

Antibody	Clone	Rhesus monkey	Cynomolgus monkey	Pigtail monkey	Chimpanzee	Baboon	Cottontop tamarin	African green monkey	Common marmoset	Other species
		<i>Macaca mulatta</i>	<i>Macaca fascicularis</i>	<i>Macaca nemestrina</i>	<i>Pan troglodytes</i>	<i>Papio anubis</i>	<i>Sanguinus oedipus</i>	<i>Chlorocebus aethiops</i>	<i>Callithrix jacchus</i>	
Anti-FoxP3	3G3	√	√*					no		
Anti-IFN-α	LT27:295							√ ⁷		
Anti-IFN-γ	45-15	√	√							
Anti-Igλ Light Chain	IS7-24C7	√	√							
Anti-IFN-γ	45-15	√	√							
Anti-IL-2	N7.48A	√	√					√ ⁷		
Anti-IL-4	7A3-3	√	√							
Anti-IL-5	JES1-39D10	√	√							
Anti-IL-10	B-T10	√	√					no		
Anti-IL-12	8.6	√	√							
Anti-IL-13	JES 10-5A2.2	√								
Anti-IL-17	CZ8-23G1	√	√							
Anti-KIR2D	NKVFS1	√	no							
Anti-MSCA-1	W8B2									pig, goat
CD1c	AD5-8E7	√	√	√				√ ⁷		
CD2	LT2	√	√		√ ¹			√ ⁷		squirrel monkey ¹
CD3	OKT3	no								horse ⁶
CD3	10D12	√	√	√				√ ⁷		
CD4	M-T466	√	√					√ ⁷		horse ⁶
CD6	M-T411	√	√							
CD8	BW135/80	√	√					√ ⁷		
CD11b	M1/70.15.11.5	√	√					√ ⁷		horse ⁵
CD14	Tük4	√	√	√			√ ²	√ ⁷		pig, cow, sheep, goat, dog, mink, rabbit ³ , cat
CD16	VEP13	√	√					√ ⁷		
CD20	LT20	√	√					√ ⁷		
CD25	4E3	√	√					√ ⁷		
CD27	M-T271	√	√			√ ⁴				
CD28	15E8	√	√					√ ⁷		
CD31	AC128	√	no							
CD33	AC104.3E3	√	√							

Antibody	Clone	Rhesus monkey <i>Macaca mulatta</i>	Cynomolgus monkey <i>Macaca fascicularis</i>	Pigtail monkey <i>Macaca nemestrina</i>	Chimpanzee <i>Pan troglodytes</i>	Baboon <i>Papio anubis</i>	Cottontop tamarin <i>Sanguinus oedipus</i>	Africangreen monkey <i>Chlorocebus aethiops</i>	Common marmoset <i>Callithrix jacchus</i>	Other species
CD34	AC136	no								horse ^{5,6}
CD39	MZ18-23C8	√								
CD40	HB14	√	√							
CD45	MB4-6D6	√	√	√						human
CD45R (B220)	RA3-6B2									horse ⁵
CD45RA	T6D11	√	√							
CD49d	MZ-24A9	√	√							
CD56	AF12-7H3	√	√							
CD61	Y2/51	√								
CD66abce	TET2	√	√							
CD69	FN50	√	√			√ ⁴		√ ⁷		
CD95	DX2	√	√			√ ⁴		√ ⁷		
CD123	AC145							no		
CD137	4B4-1	√	√			√ ⁴				
CD154	5C8	√	√					√ ⁷		
CD226	DX11	√	√							
CD271 (LNGFR)	ME20.4-1.H4									dog, cat, rabbit, pig, raccoon ⁸
CD303 (BDCA-2)	AC144	no	√	√	√ [*]	no		√ ⁷	no	
CD314 (NKG2D)	BAT221	√	√	√ ⁹						
CD337 (NKp30)	AF29-4D12	√	√							

* weak staining
√ tested to be reactive
no tested to be non-reactive

References:

- Ozwarra *et al.* (1997) *J. Med. Primatol.* 26: 164–171.
- Wilson *et al.* (1995) *J. Immunol. Meth.* 178: 195–200.
- Jacobsen *et al.* (1993) *Vet. Immunol. Immunopathol.* 39: 461–466.
- Lin *et al.* (2002) *Leucocyte Typing VII*, 519–524.
- Steinbach *et al.* (2007) *Vet. Immunol. Immunopathol.* 119: 100–105.
- Ibrahim and Steinbach (2007) *Vet. Immunol. Immunopathol.* 119: 81–91.
- Information provided by Dr. D. Kunkel, Institut Pasteur, Paris.
- Tremere *et al.* (2000) *Neuroreport* 11: 2177–2183.
- Mavilio *et al.* (2005) *Blood* 106: 1718–1725.



Miltenyi Biotec

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

www.miltenyibiotec.com

Miltenyi Biotec Inc.
12740 Earhart Avenue
Auburn, CA 95602, USA
Phone 800 FOR MACS,
+1 530 888 8871
Fax +1 530 888 8925
macs@miltenyibiotec.com

Miltenyi Biotec Australia Pty. Ltd.
Phone +61 02 8877 7400
macs@miltenyibiotec.com.au

Miltenyi Biotec B.V. (Benelux)
macs@miltenyibiotec.nl
Customer service Netherlands
Phone 0800 4020120
Customer service Belgium
Phone 0800 94016
Customer service Luxembourg
Phone 800 24971

Miltenyi Biotec Trading (Shanghai) Co., Ltd. (P.R. China)
Phone +86 21 6235 1005
macs@miltenyibiotec.com.cn

Miltenyi Biotec SAS (France)
Phone +33 1 56 98 16 16
macs@miltenyibiotec.fr

Miltenyi Biotec S.r.l. (Italy)
Phone +39 051 646 0411
macs@miltenyibiotec.it

Miltenyi Biotec K.K. (Japan)
Phone +81 3 5646 8910
macs@miltenyibiotec.jp

Miltenyi Biotec Asia Pacific Pte. Ltd. (Singapore)
Phone +65 6238 8183
macs@miltenyibiotec.com.sg

Miltenyi Biotec S.L. (Spain)
Phone +34 91 512 12 90
macs@miltenyibiotec.es

Miltenyi Biotec Ltd. (UK)
Phone +44 1483 799 800
macs@miltenyibiotec.co.uk

MACS is a registered trademark of Miltenyi Biotec GmbH.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for diagnostic or therapeutic use.