




Excerpt from

Miltenyi Biotec Flyer 2005

autoMACS™ Separator



Technical specifications

| | |
|-------------------------|-----------------------------------|
| Programs | 9 preset |
| Protocol | single or double enrichment |
| Speed (1 mL) | 10 ⁷ cells/sec |
| Program duration (1 mL) | 2.0–15 min* |
| Column capacity | 2 × 10 ⁸ labeled cells |
| Input volume | 0.25 mL–50 mL |
| Output volume | 0.5 mL–50 mL |
| Size | 500 × 370 × 370 mm |
| Weight | 21 kg |
| Input voltage | 100–250 V AC 50–60 Hz |
| Power consumption | 50 VA |

* dependent on the separation program chosen.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use.

Miltenyi Biotec

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

www.miltenyibiotec.com

Miltenyi Biotec Inc.
12740 Earhart Avenue
Auburn CA 95602, USA
Phone 800 FOR MACS,
+1 530 888 8871
Fax +1 530 888 8925
macs@miltenyibiotec.com

**Miltenyi Biotec Pty. Ltd.
(Australia)**
Phone +61 02 8877 7400
macs@miltenyibiotec.com.au

Miltenyi Biotec B. V. (Benelux)
macs@miltenyibiotec.nl
Customer service, Netherlands
Phone 0800 4020120
Customer service, Belgium
Phone 0800 94016
Customer service, Luxembourg
Phone 800 24971

Miltenyi Biotec Shanghai Office
Phone +86 21 6235 1005
macs@miltenyibiotec.com.cn

Miltenyi Biotec (France)
Phone +33 1 56 98 16 16
macs@miltenyibiotec.fr

Miltenyi Biotec S.r.l. (Italy)
Phone +39 051 646 0411
macs@miltenyibiotec.it

Miltenyi Biotec K.K. (Japan)
Phone +81 3 56 46 8910
macs@miltenyibiotec.jp

**Miltenyi Biotec Asia Pacific
Pte. Ltd. (Singapore)**
Phone +65 6238 8183
macs@miltenyibiotec.com.sg

Miltenyi Biotec S.L. (Spain)
Phone +34 91 512 12 90
macs@miltenyibiotec.es

Miltenyi Biotec Ltd. (UK)
Phone +44 1483 799 800
macs@miltenyibiotec.co.uk

