

Cell type-specific Total RNA

— ready-to-use!

Gene expression analysis from distinct hematopoietic cell populations

There is no need to collect large pools of blood and to separate specific cell populations. The high-quality RNA is ready to use—as sample or reference for gene expression profiling such as quantitative PCR or microarray analysis. Using the total RNA as reference allows comparison of treated cell culture samples or patient samples to healthy, untreated conditions. Alternatively, the RNA can be used to clone target genes following cDNA synthesis or RT-PCR.


Reliable cell-type specific RNA

Using healthy donor blood, specific cell fractions are separated with the standard method for cell isolation, MACS® Technology. This technology is particularly advantageous for the isolation of rare and immunologically relevant cell populations such as hematopoietic progenitor cells (CD34⁺) and regulatory T cells (CD4⁺ CD25⁺). The extremely small MicroBeads and gentle MACS Column Technology enable cell isolation with generally no effect on cell function or activity.

High-quality RNA

Separated cell fractions from multiple donors are pooled to contain cells of different donor ages and equal gender distribution (50% female, 50% male). RNA is then isolated and extensively quality controlled to provide highly pure cell-type specific total RNA.

For more information please refer to www.MACSmolecular.com.



Real-time PCR analysis of Regulatory T cell (CD4⁺ CD25⁺) Total RNA. Quantitative RT-PCR with intron-spanning primers for a 88 bp-fragment of Hypoxanthineguanine-phosphoribosyltransferase (HPRT) low-copy housekeeping gene (LightCycler®, Roche).

Product	RNA	Order no.
Hematopoietic Progenitor Cell (CD34 ⁺) Total RNA, human	0.2 µg	130-093-167
Regulatory T Cell (CD4 ⁺ CD25 ⁺) Total RNA, human	0.5 µg	130-093-165
T Helper Cell (CD4 ⁺) Total RNA, human	5 µg	130-093-163
Cytotoxic T Cell (CD8 ⁺) Total RNA, human	5 µg	130-093-168
T Cell (CD3 ⁺) Total RNA, human	5 µg	130-093-164
B Cell (CD19 ⁺) Total RNA, human	5 µg	130-093-169
Monocyte (CD14 ⁺) Total RNA, human	5 µg	130-093-166



Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

www.miltenyibiotec.com

Miltenyi Biotec Inc.
12740 Earhart Avenue
Auburn CA 95602, USA
Phone 800 FOR MACS,
+1 530 888 8871
Fax +1 530 888 8925
macs@miltenyibiotec.com

Miltenyi Biotec Pty. Ltd. (Australia)
Phone +61 02 8877 7400
macs@miltenyibiotec.com.au

Miltenyi Biotec B. V. (Benelux)
macs@miltenyibiotec.nl
Customer service, Netherlands
Phone 0800 4020120
Customer service, Belgium
Phone 0800 94016
Customer service, Luxembourg
Phone 800 24971

Miltenyi Biotec Shanghai Office
Phone +86 21 6235 1005
macs@miltenyibiotec.com.cn

Miltenyi Biotec SAS (France)
Phone +33 1 56 98 16 16
macs@miltenyibiotec.fr

Miltenyi Biotec S.r.l. (Italy)
Phone +39 051 646 0411
macs@miltenyibiotec.it

Miltenyi Biotec K.K. (Japan)
Phone +81 3 56 46 8910
macs@miltenyibiotec.jp

Miltenyi Biotec Asia Pacific Pte. Ltd. (Singapore)
Phone +65 6238 8183
macs@miltenyibiotec.com.sg

Miltenyi Biotec S.L. (Spain)
Phone +34 91 512 12 90
macs@miltenyibiotec.es

Miltenyi Biotec Ltd. (UK)
Phone +44 1483 799 800
macs@miltenyibiotec.co.uk