



GMP quality

TexMACS™ GMP Medium

High-performance complete T cell expansion medium

Reliable

High lot-to-lot consistency
and extensive documentation

Optimized

For the demands of T cells and Treg cells

Trusted

Benefit from over 14 years of experience
in manufacturing of GMP products

► miltenyibiotec.com/gmp

TexMACS™ GMP Medium

This medium has been specifically developed for the demands of T cells and Treg cells. TexMACS™ GMP Medium has been designed for high-performance T cell growth, high cell viability, and consistency under serum-free conditions. TexMACS GMP Medium is serum and animal-component free.

- Contains pharmaceutical-grade human serum albumin
- Optimized formulation with or without phenol red
- Complete medium with glucose and stable glutamine (L-alanyl-L-glutamine)
- QC functionality test on every batch
- Filled in flexible bags (connection via Luer fitting or sterile welding)
- Customized formats on request

TexMACS™ Medium is also available as research grade

- Optimized formulation with phenol red
- Filled in bottles

MACS® GMP Products – GMP-compliant manufacturing of T cells

The quality of starting materials dramatically influences the properties of the final cellular product. For the development of applications involving cellular therapies, strict procedures according to GMP guidelines are mandatory. All materials used, such as antibodies, cell culture bags, media and cytokines, must comply with GMP specifications.

GMP manufacturing

- Manufactured in our state-of-the-art facility in Teterow, Germany, in accordance with Good Manufacturing Practice standards
- Designed according to the recommendations of the United States Pharmacopeia, chapter <1043> on ancillary materials
- Manufactured and tested under a certified ISO 13485 quality system

Safety

- Intended for research use and *ex vivo* cell culture processing
- Free of animal-derived contaminants

Extensive documentation

- Batch-specific certificate of analysis (CoA)
- Product data sheet
- Regulatory support document cross reference

Product	Format	Capacity	Order no.
TexMACS GMP Medium	Flexible bag	2000 mL	170-076-306
TexMACS GMP Medium	Bottle	1000 mL	170-076-307
TexMACS GMP Medium (Phenol Red)	Bottle	1000 mL	170-076-309
TexMACS Medium – research grade	Bottle	500 mL	130-097-196

For an extensive overview of MACS® GMP Products please visit www.miltenyibiotec.com/gmp



Miltenyi Biotec

Miltenyi Biotec GmbH
Friedrich-Ebert-Straße 68
51429 Bergisch Gladbach
Germany
Phone +49 2204 8306-0
Fax +49 2204 85197
macs@miltenyibiotec.de

Miltenyi Biotec Inc.
2303 Lindbergh Street
Auburn, CA 95602-9562, USA
Phone 800 FOR MACS
Phone +1 530 888 8871
Fax +1 877 591 1060
macs@miltenyibiotec.com

Miltenyi Biotec Australia Pty. Ltd.
Phone +61 2 8877 7400
macs@miltenyibiotec.com.au

Miltenyi Biotec B.V. (Benelux)
macs@miltenyibiotec.nl
Customer service The Netherlands
Phone 0800 4020120
Customer service Belgium
Phone 0800 94016
Customer service Luxembourg
Phone 800 24971

Miltenyi Biotec Technology & Trading (Shanghai) Co., Ltd. (P.R. China)
Phone +86 21 62351005
macs@miltenyibiotec.com.cn

Miltenyi Biotec SAS (France)
Phone +33 1 56 98 16 16
macs@miltenyibiotec.fr

Miltenyi Biotec S.r.l. (Italy)
Phone +39 051 646 0411
macs@miltenyibiotec.it

Miltenyi Biotec K.K. (Japan)
Phone +81 3 5646 8910
macs@miltenyibiotec.jp

Miltenyi Biotec Norden AB (Nordics and Baltics)
macs@miltenyibiotec.se
Customer service Sweden
Phone 0200-111 800
Customer service Denmark
Phone 80 20 30 10
Customer service Norway, Finland, Iceland, and Baltic countries
Phone +46 46 280 72 80

Miltenyi Biotec Asia Pacific Pte. Ltd. (Singapore)
Phone +65 6238 8183
macs@miltenyibiotec.com.sg

Miltenyi Biotec Korea Co., Ltd (South Korea)
Phone +82 2 555 1988
macs@miltenyibiotec.co.kr

Miltenyi Biotec S.L. (Spain)
Phone +34 91 512 12 90
macs@miltenyibiotec.es

Miltenyi Biotec Ltd. (UK)
Phone +44 1483 799 800
macs@miltenyibiotec.co.uk

www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

MACS GMP Products are for research use and *ex vivo* cell culture processing only, and are not intended for human *in vivo* applications. For regulatory status in the USA, please contact your local representative. MACS and TexMACS are registered trademarks or trademarks of Miltenyi Biotec GmbH. Copyright © 2018 Miltenyi Biotec GmbH. All rights reserved.