

Contents

1. Description

1.1 Background information

1.2 Applications

2. References

1. Description

Products	Mouse IFN- γ , research grade. Recombinant mouse interferon γ .						
	<table border="1"> <thead> <tr> <th>Content in μg</th><th>Order no.</th></tr> </thead> <tbody> <tr> <td>10</td><td>130-105-790</td></tr> <tr> <td>25</td><td>130-105-785</td></tr> </tbody> </table>	Content in μ g	Order no.	10	130-105-790	25	130-105-785
Content in μ g	Order no.						
10	130-105-790						
25	130-105-785						
Biological activity	The ED ₅₀ is ≤ 1.0 ng/mL corresponding to an activity of $\geq 1 \times 10^6$ IU/mg. Note: The ED ₅₀ is determined by inhibition assay using murine WEHI-279 cells according to Wong <i>et al.</i> ¹ The inhibition assay was calibrated with the international standard for mouse IFN- γ (NR-3081).						
Primary structure	Single non-glycosylated polypeptide chain without N-terminal methionine (133 amino acid residues).						
Molecular mass	15.5 kDa.						
Source	Produced in <i>E. coli</i> .						
Product format	Lyophilized from a filtered (0.2 μ m) buffer solution.						
Stabilizer	Mannitol and trehalose.						
Purity	>90% as determined by SDS-PAGE analysis.						
Endotoxin level	Low endotoxin (<1.0 EU/ μ g cytokine) as determined by Limulus Amebocyte Lysate (LAL) assay.						
Storage	Lyophilized Mouse IFN- γ , research grade should be stored at -20 °C. The expiration date is indicated on the vial label. Upon reconstitution aliquots should be stored at -20 °C or below. Avoid repeated freeze-thaw cycles.						
Reconstitution	It is recommended to reconstitute lyophilized Mouse IFN- γ , research grade with deionized sterile-filtered water to a final concentration of 0.1–1.0 mg/mL in a minimal volume of 100 μ L. Further dilutions should be prepared with 0.1% bovine serum albumin (BSA) or human serum albumin (HSA) in phosphate-buffered saline.						

1.1 Background information

Interferon γ (IFN- γ) is a homodimer produced by a variety of immune cells, mainly activated T cells and NK cells. In addition to its antiviral activity², IFN- γ is a potent activator of macrophages, has antiproliferative effects on transformed cells, and can potentiate the antiviral and antitumor effects of type I interferons. Transgenic mice expressing IFN- γ in the pancreas have been shown to develop autoimmune diabetes, caused by infiltration of lymphocytes. IFN- γ knockout mice develop normally, but are highly susceptible to infection with low numbers of certain intracellular bacterial pathogens and fungi. Mouse IFN- γ shows approximately 40% amino acid sequence homologous with human IFN- γ . It is highly specific without showing cross-reactivity between species.

1.2 Applications

Mouse IFN- γ may be used for a variety of applications, including:

- Macrophage stimulation to study antimicrobial and anti-tumor effects.
- *In vitro* maturation of DCs, e.g. of MDCs to study the influence of HIV-1 infection.

Optimal concentration for a specific application should be determined by a dose-response experiment.

2. References

1. Wong, G. H. *et al.* (1983) Interferon-gamma induces enhanced expression of Ia and H-2 antigens on B lymphoid, macrophage, and myeloid cell lines. *J Immunol* 131: 788–793.
2. Groot, F. *et al.* (2006) Opposing roles of blood myeloid and plasmacytoid dendritic cells in HIV-1 infection of T cells: transmission facilitation versus replication inhibition. *Blood* 108: 1957–1964.

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com for local Miltenyi Biotec Technical Support contact information.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE

TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

The Miltenyi Biotec logo is a registered trademark or trademark of Miltenyi Biotec and/or its affiliates in various countries worldwide.

Copyright © 2023 Miltenyi Biotec and/or its affiliates. All rights reserved.