




MACS® Marker Screen Plate allocations

MACS Marker Screen Plate 1

	1	2	3	4	5	6	7	8	9	10	11	12
A	Empty	CD1a	CD1b	CD1c (BDCA-1)	CD1d	CD2	CD3	CD4 (VIT4)	CD5	CD6	CD7	CD8
B	CD9	CD10	CD11a	CD11b	CD11c	CD13	CD14	CD15	CD16	CD18	CD19	CD20
C	CD21	CD22	CD23	CD24	CD25	CD26	CD27	CD28	CD29	CD30	CD31	CD32
D	CD33	CD34	CD35	CD36	CD37	CD38	CD39	CD40	CD41a	CD41b	CD42a	CD42b
E	CD43	CD44	CD45	CD45RA	CD45RB	CD45RO	CD46	CD47	CD48	CD49a	CD49b	CD49c
F	CD49d	CD49e	CD49f	CD51	CD51/CD61	CD52	CD53	CD54 (ICAM-1)	CD55 (DAF)	CD56	CD57	CD58 (LFA-3)
G	CD61	CD62E	CD62L	CD62P	CD63	CD64	CD66acde	CD66b	CD66c	CD68	CD69	CD70
H	CD71	CD72	CD73	CD74	CD79a	CD79b	CD80	CD81	CD82	CD83	CD84	CD85a (ILT5)

Antigen	Clone	Isotype	Order no.
CD1b	SN13	mouse IgG1k	130-098-337
CD1c (BDCA-1)	AD5-8E7	mouse IgG2a	130-098-005
CD1d	51.1	mouse IgG2bk	130-099-980
CD2	LT2	mouse IgG2b	130-098-576
CD3	BW264/56	mouse IgG2a	130-098-156
CD4 (VIT4)	VIT4	mouse IgG2a	130-098-158
CD5	UCHT2	mouse IgG1k	130-100-145
CD6	M-T411	mouse IgG1	130-100-012

Antigen	Clone	Isotype	Order no.
CD7	CD7-6B7	mouse IgG2ak	130-105-900
CD8	BW135/80	mouse IgG2a	130-098-078
CD9	SN4 C3-3A2	mouse IgG1	130-103-989
CD10	97C5	mouse IgG1	130-100-082
CD11a	REA378	recombinant human IgG1	130-105-480
CD11b	M1/70.15.11.5	rat IgG2bk	130-098-088
CD11c	MJ4-27G12	mouse IgG2b	130-092-412
CD13	REA263	recombinant human IgG1	130-103-734
CD14	TÜK4	mouse IgG2a	130-098-070
CD15	VIMC6	mouse IgM	130-098-008
CD16	REA423	recombinant human IgG1	130-106-763
CD18	TS1/18	mouse IgG1k	130-101-226
CD19	LT19	mouse IgG1	130-098-069
CD20	LT20	mouse IgG1	130-097-619
CD21	HB5	mouse IgG2ak	130-101-661
CD22	REA340	recombinant human IgG1	130-105-085
CD23	M-L23.4	mouse IgG1	130-100-302
CD24	32D12	mouse IgG1	130-099-399
CD25	4E3	mouse IgG2b	130-098-213
CD26	FR10-11G9	mouse IgG2a	130-099-737
CD27	REA499	recombinant human IgG1	130-108-253
CD28	15E8	mouse IgG1	130-099-191
CD29	TS2/16	mouse IgG1k	130-101-271
CD30	Ki-2	mouse IgG1	130-098-615
CD31	AC128	mouse IgG1	130-098-174
CD32	2E1	mouse IgG2a	130-097-525
CD33	AC104.3E3	mouse IgG1	130-098-864
CD34	AC136	mouse IgG2a	130-098-139
CD35	E11	mouse IgG1k	130-099-923
CD36	AC106	mouse IgG2a	130-100-307
CD37	REA366	recombinant human IgG1	130-105-287
CD38	REA572	recombinant human IgG1	130-108-863
CD39	MZ18-23C8	mouse IgG1	130-100-459
CD40	HB14	mouse IgG1k	130-099-186
CD41a	REA386	recombinant human IgG1	130-105-613
CD41b	REA336	recombinant human IgG1	130-105-072
CD42a	REA209	recombinant human IgG1	130-100-934
CD42b	REA185	recombinant human IgG1	130-100-208
CD43	DF-T1	mouse IgG1	130-097-367
CD44	DB105	mouse IgG1	130-098-110
CD45	5B1	mouse IgG2a	130-098-143
CD45RA	REA562	recombinant human IgG1	130-108-782
CD45RB	REA119	recombinant human IgG1	130-100-393
CD45RO	UCHL1	mouse IgG2ak	130-099-395
CD46	REA312	recombinant human IgG1	130-104-558
CD47	REA220	recombinant human IgG1	130-101-337
CD48	REA426	recombinant human IgG1	130-106-517
CD49a	TS2/7	mouse IgG1k	130-101-400
CD49b	REA188	recombinant human IgG1	130-100-348
CD49c	REA360	recombinant human IgG1	130-105-407

Antigen	Clone	Isotype	Order no.
CD49d	REA545	recombinant human IgG1	130-108-258
CD49e	NKI-SAM1	mouse IgG2bk	130-097-221
CD49f	GoH3	rat IgG2ak	130-100-147
CD51	REA181	recombinant human IgG1	130-100-552
CD51/CD61	23C6	mouse IgG1k	130-103-745
CD52	REA164	recombinant human IgG1	130-099-630
CD53	REA259	recombinant human IgG1	130-101-795
CD54 (ICAM-1)	REA266	recombinant human IgG1	130-103-910
CD55 (DAF)	JS11	mouse IgG1k	130-101-802
CD56	REA196	recombinant human IgG1	130-100-697
CD57	TB03	mouse IgM	130-099-748
CD58 (LFA-3)	TS2/9	mouse IgG1k	130-101-191
CD61	Y2/51	mouse IgG1	130-098-577
CD62E	REA280	recombinant human IgG1	130-104-685
CD62L	145/15	mouse IgG1	130-099-252
CD62P	REA389	recombinant human IgG1	130-105-715
CD63	H5C6	mouse IgG1k	130-100-173
CD64	10.1.1	mouse IgG1k	130-100-417
CD66acde	REA428	recombinant human IgG1	130-106-418
CD66b	REA306	recombinant human IgG1	130-104-415
CD66c	REA414	recombinant human IgG1	130-106-379
CD68	Y1/82A	mouse IgG2bk	130-099-903
CD69	FN50	mouse IgG1k	130-098-893
CD70	REA292	recombinant human IgG1	130-104-358
CD71	AC102	mouse IgG2a	130-099-239
CD72	REA231	recombinant human IgG1	130-101-334
CD73	AD2	mouse IgG1k	130-097-945
CD74	5-329	mouse IgG1	130-101-543
CD79a	HM47	mouse IgG1k	130-098-693
CD79b	REA120	recombinant human IgG1	130-099-103
CD80	2D10	mouse IgG1k	130-099-710
CD81	REA513	recombinant human IgG1	130-107-983
CD82	REA221	recombinant human IgG1	130-101-311
CD83	HB15	mouse IgG1k	130-098-889
CD84	MZ18-21F6	mouse IgG1	130-093-276
CD85a (ILT5)	REA207	recombinant human IgG1	130-101-665

MACS Marker Screen Plate 2

	1	2	3	4	5	6	7	8	9	10	11	12
A	Empty	CD85d (ILT4)	CD85g (ILT7)	CD85h (ILT1)	CD85j (ILT2)	CD85k (ILT3)	CD86	CD87	CD88 (C5AR)	CD89	CD90	CD93
B	CD94	CD95	CD96 (TACTILE)	CD97	CD98	CD99	CD100	CD101	CD103	CD104 (Integrin β 4)	CD105	CD106
C	CD107a	CD107b	CD110	CD111	CD114	CD116	CD117 (A3C6E2)	CD119	CD120a	CD122 (IL-2R β)	CD123	CD126 (IL-6RA)
D	CD127	CD131	CD132	CD133/1 (AC133)	CD133/2 (293C3)	CD134 (OX40)	CD135	CD137	CD137L (4-1BBL)	CD138	CD140b	CD141 (BDCA-3)
E	CD142	CD143 (ACE)	CD144 (VE-Cadherin)	CD146	CD147	CD148	CD150 (SLAM)	CD151	CD152	CD154	CD155	CD156a (ADAM8)
F	CD156c (ADAM10)	CD157 (BST-1)	CD158a (KIR2DL1)	CD158a/h (KIR2DL1/DS1)	CD158b (KIR2DL2/DL3)	CD158b2 (KIR2DL3)	CD158e (KIR3DL1)	CD158e1/e2	CD158f	CD158i (KIR2DS4)	CD159a (NKG2A)	CD159c (NKG2C)
G	CD161	CD162	CD163	CD164	CD165	CD166	CD169 (Siglec-1)	CD170 (Siglec-5)	CD171 (L1CAM)	CD172a (SIRPa)	CD172b (SIRP β)	CD172g (SIRPy)
H	CD177	CD178	CD180 (RP105)	CD181 (CXCR1)	CD182 (CXCR2)	CD183 (CXCR3)	CD184 (CXCR4)	CD185 (CXCR5)	CD191 (CCR1)	CD192 (CCR2)	CD193 (CCR3)	CD194 (CCR4)

Antigen	Clone	Isotype	Order no.
CD85d (ILT4)	REA184	recombinant human IgG1	130-100-561
CD85g (ILT7)	REA100	recombinant human IgG1	130-099-008
CD85h (ILT1)	REA219	recombinant human IgG1	130-100-920
CD85j (ILT2)	GHI/75	mouse IgG2bk	130-098-446
CD85k (ILT3)	REA141	recombinant human IgG1	130-099-799
CD86	FM95	mouse IgG1	130-097-920
CD87	VIM5	mouse IgG1k	130-099-470
CD88 (C5AR)	S5/1	mouse IgG2ak	130-104-335
CD89	REA234	recombinant human IgG1	130-101-384
CD90	DG3	mouse IgG1	130-097-935
CD93	VIMD2	mouse IgG1k	130-098-441
CD94	REA113	recombinant human IgG1	130-098-978
CD95	DX2	mouse IgG1k	130-099-738
CD96 (TACTILE)	REA195	recombinant human IgG1	130-101-031
CD97	VIM3b	mouse IgG1k	130-097-100
CD98	REA387	recombinant human IgG1	130-105-709
CD99	3B2/TA8	mouse IgG2ak	130-104-365
CD100	REA316	recombinant human IgG1	130-104-674
CD101	BB27	mouse IgG1k	130-095-132
CD103	Ber-ACT 8	mouse IgG1k	130-103-710
CD104 (Integrin β 4)	REA236	recombinant human IgG1	130-101-438
CD105	43A4E1	mouse IgG1	130-099-125
CD106	REA269	recombinant human IgG1	130-104-164
CD107a	H4A3	mouse IgG1k	130-100-314

Antigen	Clone	Isotype	Order no.
CD107b	H4B4	mouse IgG1k	130-103-993
CD110	REA250	recombinant human IgG1	130-101-647
CD111	R1.302	mouse IgG1k	130-103-904
CD114	LMM741	mouse IgG1k	130-097-308
CD116	REA211	recombinant human IgG1	130-100-986
CD117 (A3C6E2)	A3C6E2	mouse IgG1k	130-098-207
CD119	REA161	recombinant human IgG1	130-099-920
CD120a	REA252	recombinant human IgG1	130-106-361
CD122 (IL-2Rβ)	REA167	recombinant human IgG1	130-100-248
CD123	AC145	mouse IgG2a	130-098-873
CD126 (IL-6RA)	REA291	recombinant human IgG1	130-104-138
CD127	MB15-18C9	mouse IgG2a	130-098-097
CD131	REA249	recombinant human IgG1	130-104-035
CD132	REA313	recombinant human IgG1	130-104-660
CD133/1 (AC133)	AC133	mouse IgG1k	130-098-829
CD133/2 (293C3)	293C3	mouse IgG2bk	130-098-129
CD134 (OX40)	ACT35	mouse IgG1k	130-095-272
CD135	BV10A4H2	mouse IgG1k	130-106-356
CD137	4B4-1	mouse IgG1k	130-098-912
CD137L (4-1BBL)	REA254	recombinant human IgG1	130-103-721
CD138	44F9	mouse IgG1	130-098-747
CD140b	REA363	recombinant human IgG1	130-105-322
CD141 (BDCA-3)	AD5-14H12	mouse IgG1	130-098-840
CD142	HTF-1	mouse IgG1k	130-098-745
CD143 (ACE)	REA522	recombinant human IgG1	130-108-014
CD144 (VE-Cadherin)	REA199	recombinant human IgG1	130-100-709
CD146	541-10B2	mouse IgG1	130-097-942
CD147	REA282	recombinant human IgG1	130-104-540
CD148	REA204	recombinant human IgG1	130-101-183
CD150 (SLAM)	REA151	recombinant human IgG1	130-099-722
CD151	REA265	recombinant human IgG1	130-103-728
CD152	BNI3	mouse IgG2ak	130-097-682
CD154	5C8	mouse IgG2a	130-098-883
CD155	PV404.19	mouse IgG1	130-105-906
CD156a (ADAM8)	REA331	recombinant human IgG1	130-104-896
CD156c (ADAM10)	REA309	recombinant human IgG1	130-104-426
CD157 (BST-1)	REA465	recombinant human IgG1	130-106-982
CD158a (KIR2DL1)	REA284	recombinant human IgG1	130-103-968
CD158a/h (KIR2DL1/DS1)	11PB6	mouse IgG1	130-099-715
CD158b (KIR2DL2/DL3)	DX27	mouse IgG2ak	130-099-378
CD158b2 (KIR2DL3)	REA147	recombinant human IgG1	130-100-130
CD158e (KIR3DL1)	DX9	mouse IgG1k	130-099-839
CD158e1/e2	REA168	recombinant human IgG1	130-104-534
CD158f	UP-R1	mouse IgG1k	130-098-567
CD158i (KIR2DS4)	JJC11.6	mouse IgG1	130-099-709
CD159a (NKG2A)	REA110	recombinant human IgG1	130-098-809
CD159c (NKG2C)	REA205	recombinant human IgG1	130-103-701
CD161	191B8	mouse IgG2a	130-098-908
CD162	REA319	recombinant human IgG1	130-104-755
CD163	GHI/61.1	mouse IgG1k	130-100-612

Antigen	Clone	Isotype	Order no.
CD164	67D2	mouse IgG1κ	130-108-105
CD165	SN2	mouse IgG1κ	130-098-694
CD166	REA442	recombinant human IgG1	130-106-619
CD169 (Siglec-1)	7-239	mouse IgG1κ	130-098-643
CD170 (Siglec-5)	1A5	mouse IgG1κ	130-101-760
CD171 (L1CAM)	REA163	recombinant human IgG1	130-100-686
CD172a (SIRPα)	REA144	recombinant human IgG1	130-099-785
CD172b (SIRPβ)	B4B6	mouse IgG1κ	130-105-484
CD172g (SIRPγ)	REA479	recombinant human IgG1	130-107-153
CD177	REA258	recombinant human IgG1	130-101-512
CD178	NOK-1	mouse IgG1κ	130-096-458
CD180 (RP105)	MHR73-11	mouse IgG1κ	130-098-490
CD181 (CXCR1)	8F1	mouse IgG2bk	130-105-394
CD182 (CXCR2)	REA208	recombinant human IgG1	130-100-902
CD183 (CXCR3)	REA232	recombinant human IgG1	130-101-378
CD184 (CXCR4)	12G5	mouse IgG2ak	130-098-357
CD185 (CXCR5)	REA103	recombinant human IgG1	130-098-422
CD191 (CCR1)	REA158	recombinant human IgG1	130-100-360
CD192 (CCR2)	REA264	recombinant human IgG1	130-103-900
CD193 (CCR3)	REA574	recombinant human IgG1	130-108-920
CD194 (CCR4)	REA279	recombinant human IgG1	130-103-883

MACS Marker Screen Plate 3

	1	2	3	4	5	6	7	8	9	10	11	12
A	Empty	CD195 (CCR5)	CD196 (CCR6)	CD197 (CCR7)	CD199 (CCR9)	CD200	CD201 (EPCR)	CD202b (TIE-2)	CD203c	CD204	CD205	CD206
B	CD207 (Langerin)	CD208 (DC-LAMP)	CD209 (DC-SIGN)	CD212	CD213a2 (IL-13Ra2)	CD217	CD218 (IL-18Ra)	CD220	CD221 (IGF-1R)	CD222	CD223	CD226 (DNAM-1)
C	CD227 (Muc-1)	CD229 (Ly-9)	CD230 (PrP)	CD231 (TALLA)	CD233	CD234 (DARC)	CD235a (Glycophorin A)	CD238 (KEL)	CD239 (BCAM)	CD240DCE	CD243 (ABCB1)	CD244 (2B4)
D	CD246 (ALK)	CD253	CD256 (APRIL)	CD258 (LIGHT)	CD262	CD263 (TRAIL-R3)	CD266 (FN14)	CD268	CD269 (BCMA)	CD270 (HVEM)	CD271 (LNGFR)	CD273 (PD-L2)
E	CD275 (B7-H2)	CD276	CD277	CD278 (ICOS)	CD279 (PD-1)	CD282 (TLR2)	CD283	CD284	CD286 (TLR6)	CD295 (LEPR)	CD298	CD300e (IREM-2)
F	CD300f (IREM-1)	CD302 (CLEC13A)	CD303 (BDCA-2)	CD304 (BDCA-4/Neuropilin-1)	CD305 (LAIR-1)	CD307a (FcRL1)	CD307b (FcRL2)	CD307e (FcRL5)	CD309 (VEGFR-2/KDR)	CD310 (VEGFR-3)	CD312 (EMR2)	CD314 (NKG2D)
G	CD317 (BST2)	CD318 (CDCP1)	CD319 (CRACC)	CD324 (E-Cadherin)	CD326 (EpCAM)	CD328 (Siglec-7)	CD329 (Siglec-9)	CD335 (NKp46)	CD336 (NKp44)	CD337 (NKp30)	CD338 (ABCG2)	ErbB-2 (CD340)
H	CD344 (Frizzled-4)	CD352 (NTB-A)	CD353 (SLAMF8)	CD354 (TREM-1)	CRTAM	CD360 (IL-21R)	CD362 (Syndecan-2)	A2B5	Podoplanin	C3a Receptor	CCL2 (MCP-1)	CCR10

Antigen	Clone	Isotype	Order no.
CD195 (CCR5)	REA245	recombinant human IgG1	130-106-262
CD196 (CCR6)	REA190	recombinant human IgG1	130-100-374
CD197 (CCR7)	REA546	recombinant human IgG1	130-108-310
CD199 (CCR9)	REA469	recombinant human IgG1	130-107-170
CD200	OX-104	mouse IgG1κ	130-096-816
CD201 (EPCR)	REA337	recombinant human IgG1	130-105-299
CD202b (TIE-2)	REA198	recombinant human IgG1	130-101-632
CD203c	FR 3-16A11	mouse IgG1	130-098-888
CD204	REA460	recombinant human IgG1	130-107-062
CD205	HD30	mouse IgG1κ	130-096-378
CD206	DCN228	mouse IgG1	130-099-731
CD207 (Langerin)	MB22-9F5	mouse IgG1	130-098-364
CD208 (DC-LAMP)	REA295	recombinant human IgG1	130-104-411
CD209 (DC-SIGN)	DCN47.5	mouse IgG1	130-099-727
CD212	REA242	recombinant human IgG1	130-101-539
CD213α2 (IL-13Ra2)	REA308	recombinant human IgG1	130-104-554
CD217	REA290	recombinant human IgG1	130-104-769
CD218 (IL-18Ra)	H44	mouse IgG1κ	130-101-653
CD220	REA260	recombinant human IgG1	130-103-718
CD221 (IGF-1R)	REA271	recombinant human IgG1	130-103-973
CD222	REA187	recombinant human IgG1	130-100-863
CD223	REA351	recombinant human IgG1	130-105-496
CD226 (DNAM-1)	DX11	mouse IgG1	130-100-436
CD227 (Muc-1)	REA448	recombinant human IgG1	130-106-836
CD229 (Ly-9)	Hly9.1.25	mouse IgG1κ	130-100-114
CD230 (PrP)	REA203	recombinant human IgG1	130-101-284
CD231 (TALLA)	SN1a	mouse IgG1κ	130-098-500
CD233	REA368	recombinant human IgG1	130-105-729
CD234 (DARC)	REA376	recombinant human IgG1	130-105-731
CD235a (Glycophorin A)	REA175	recombinant human IgG1	130-100-262
CD238 (KEL)	REA330	recombinant human IgG1	130-105-003
CD239 (BCAM)	REA276	recombinant human IgG1	130-103-913
CD240DCE	REA327	recombinant human IgG1	130-104-864
CD243 (ABCB1)	REA495	recombinant human IgG1	130-107-526
CD244 (2B4)	REA112	recombinant human IgG1	130-099-073
CD246 (ALK)	REA425	recombinant human IgG1	130-106-494
CD253	RIK-2.1	mouse IgG1κ	130-097-314
CD256 (APRIL)	REA347	recombinant human IgG1	130-105-388
CD258 (LIGHT)	REA244	recombinant human IgG1	130-103-725
CD262	DJR2-4	mouse IgG1κ	130-097-249
CD263 (TRAIL-R3)	DJR3	mouse IgG1κ	130-104-788
CD266 (FN14)	ITEM-4	mouse IgG2bκ	130-104-330
CD268	11C1	mouse IgG1κ	130-097-648
CD269 (BCMA)	REA315	recombinant human IgG1	130-104-551
CD270 (HVEM)	REA247	recombinant human IgG1	130-101-616
CD271 (LNGFR)	ME20.4-1.H4	mouse IgG1κ	130-098-112
CD273 (PD-L2)	MIH18	mouse IgG1κ	130-098-652
CD275 (B7-H2)	MIH12	mouse IgG1κ	130-098-738
CD276	FM276	mouse IgG2b	130-095-522
CD277	BT3.1	mouse IgG1κ	130-108-249

Antigen	Clone	Isotype	Order no.
CD278 (ICOS)	REA192	recombinant human IgG1	130-100-719
CD279 (PD-1)	PD1.3.1.3	mouse IgG2b	130-100-440
CD282 (TLR2)	REA109	recombinant human IgG1	130-099-020
CD283	TLR 3.7	mouse IgG1	130-100-383
CD284	HTA125	mouse IgG2ak	130-100-150
CD286 (TLR6)	REA382	recombinant human IgG1	130-106-632
CD295 (LEPR)	REA361	recombinant human IgG1	130-105-211
CD298	REA217	recombinant human IgG1	130-101-263
CD300e (IREM-2)	UP-H2	mouse IgG1k	130-101-788
CD300f (IREM-1)	UP-D2	mouse IgG1k	130-098-471
CD302 (CLEC13A)	REA509	recombinant human IgG1	130-107-850
CD303 (BDCA-2)	AC144	mouse IgG1	130-097-931
CD304 (BDCA-4/Neuropilin-1)	AD5-17F6	mouse IgG1	130-098-874
CD305 (LAIR-1)	REA447	recombinant human IgG1	130-106-850
CD307a (FcRL1)	REA440	recombinant human IgG1	130-106-497
CD307b (FcRL2)	REA474	recombinant human IgG1	130-107-493
CD307e (FcRL5)	REA391	recombinant human IgG1	130-106-045
CD309 (VEGFR-2/KDR)	ES8-20E6	mouse IgG1	130-098-910
CD310 (VEGFR-3)	REA419	recombinant human IgG1	130-106-624
CD312 (EMR2)	REA302	recombinant human IgG1	130-104-656
CD314 (NKG2D)	BAT221	mouse IgG1	130-099-215
CD317 (BST2)	REA202	recombinant human IgG1	130-101-660
CD318 (CDCP1)	REA194	recombinant human IgG1	130-101-208
CD319 (CRACC)	REA150	recombinant human IgG1	130-099-569
CD324 (E-Cadherin)	67A4	mouse IgG1k	130-099-723
CD326 (EpCAM)	HEA-125	mouse IgG1	130-098-118
CD328 (Siglec-7)	REA214	recombinant human IgG1	130-101-013
CD329 (Siglec-9)	REA492	recombinant human IgG1	130-107-649
CD335 (NKP46)	9E2	mouse IgG1k	130-099-092
CD336 (NKP44)	2.29	mouse IgG1	130-099-414
CD337 (NKP30)	AF29-4D12	mouse IgG1	130-099-380
CD338 (ABCG2)	5D3	mouse IgG2bk	130-105-011
ErbB-2 (CD340)	24D2	mouse IgG1k	130-106-754
CD344 (Frizzled-4)	CH3A4A7	mouse IgG1k	130-106-614
CD352 (NTB-A)	REA339	recombinant human IgG1	130-105-595
CD353 (SLAMF8)	REA394	recombinant human IgG1	130-106-734
CD354 (TREM-1)	REA213	recombinant human IgG1	130-101-055
CRTAM	Cr24.1	mouse IgG2ak	130-096-158
CD360 (IL-21R)	REA233	recombinant human IgG1	130-101-478
CD362 (Syndecan-2)	REA468	recombinant human IgG1	130-107-534
A2B5	105HB29	mouse IgM	130-098-039
Podoplanin	REA446	recombinant human IgG1	130-107-016
C3a Receptor	hC3aRZ8	mouse IgG2b	130-108-110
CCL2 (MCP-1)	REA485	recombinant human IgG1	130-107-536
CCR10	REA326	recombinant human IgG1	130-104-867

MACS Marker Screen Plate 4

	1	2	3	4	5	6	7	8	9	10	11	12
A	Empty	ChemR23 (CMKLR1)	CLA	CLEC2A	CLEC12A	CLEC9A	CLIP	CX3CR1	DCIR	DLL1	DLL4	DR3
B	ErbB-3 (Her-3)	FcεRIα	Fibroblast	fMLP receptor	Galectin 3	Galectin-9	GARP (LRRC32)	GITR	GLAST (ACSA-1)	GPR56	HLA Class I B8	HLA Class I Bw4
C	HLA Class I Bw6	HLA-A2	HLA-A2, A28	HLA-A9	HLA-B12	HLA-B7, B27	HLA-DM	HLA-DQ	HLA-DR	HLA-DR, DP, DQ	HLA-E	IFNAR2
D	IL-1RAcP	IL-12R β2	IL-17RC	iNKT	Integrin β7	Jagged2	KIR2D	KLRG1 (MAFA)	LAP (TGF-β1)	LGR5	LT-βR	Melanoma (MCSP)
E	HLA-ABC	MICA/ MICB	MSCA-1 (W8B2)	NKp80	Notch1	Notch2	Notch3	O4	OSCAR	Perforin	Plexin-D1	PSA-NCAM
F	PSMA	PTK7 (CCK-4)	ROR1	Siglec-10	Siglec-5/ Siglec-14	Siglec-8	Slan (M-DC8)	SSEA-1	SSEA-4	SSEA-5	SUSD2	TCR-Vα7.2
G	TCR Vβ11	TCR Vβ14	TCR Vβ16	TCR Vβ23	TCR Vγ9	TCR-Vδ1	TCR-Vδ2	TCRα/β	TCRγ/δ	TdT	TIM-1	TIM-3
H	TLT-2	TRA-1-85 (CD147)	TSPAN8	Mouse IgG1 – isotype control	Mouse IgG2a – isotype control	Mouse IgG2b – isotype control	Mouse IgM – isotype control	Rat IgG1 – isotype control	Rat IgG2a – isotype control	Rat IgG2b – isotype control	Rat IgM – isotype control	REA Control (S)

Antigen	Clone	Isotype	Order no.
ChemR23 (CMKLR1)	REA455	recombinant human IgG1	130-107-643
CLA	HECA-452	rat IgMk	130-098-572
CLEC2A	REA544	recombinant human IgG1	130-108-748
CLEC12A	REA431	recombinant human IgG1	130-106-483
CLEC9A	8F9	mouse IgG2a	130-100-406
CLIP	REA296	recombinant human IgG1	130-104-421
CX3CR1	REA385	recombinant human IgG1	130-105-897
DCIR	REA329	recombinant human IgG1	130-105-029
DLL1	MHD1-314	mouse IgG1k	130-096-961
DLL4	MHD4-46	mouse IgG1k	130-100-407
DR3	JD3	mouse IgG1k	130-105-077
ErbB-3 (Her-3)	REA508	recombinant human IgG1	130-107-962
FcεRIα	CRA1	mouse IgG2b	130-095-980
Fibroblast	REA165	recombinant human IgG1	130-100-133
fMLP receptor	REA169	recombinant human IgG1	130-100-813
Galectin 3	M3/38	rat IgG2ak	130-101-313
Galectin-9	REA435	recombinant human IgG1	130-106-523
GARP (LRRC32)	REA166	recombinant human IgG1	130-103-890
GITR	DT5D3	mouse IgG1	130-098-574
GLAST (ACSA-1)	ACSA-1	mouse IgG2a	130-098-803
GPR56	REA467	recombinant human IgG1	130-106-976
HLA Class I B8	REA145	recombinant human IgG1	130-099-595
HLA Class I Bw4	REA274	recombinant human IgG1	130-103-918
HLA Class I Bw6	REA143	recombinant human IgG1	130-099-857

Antigen	Clone	Isotype	Order no.
HLA-A2	REA517	recombinant human IgG1	130-107-993
HLA-A2, A28	REA142	recombinant human IgG1	130-099-585
HLA-A9	REA127	recombinant human IgG1	130-099-540
HLA-B12	REA138	recombinant human IgG1	130-099-854
HLA-B7, B27	REA176	recombinant human IgG1	130-106-048
HLA-DM	REA406	recombinant human IgG1	130-106-169
HLA-DQ	REA303	recombinant human IgG1	130-104-546
HLA-DR	AC122	mouse IgG2a	130-098-178
HLA-DR, DP, DQ	REA332	recombinant human IgG1	130-104-870
HLA-E	3D12	mouse IgG1k	130-096-846
IFNAR2	REA124	recombinant human IgG1	130-099-560
IL-1RAcP	REA558	recombinant human IgG1	130-108-755
IL-12R β 2	REA333	recombinant human IgG1	130-105-036
IL-17RC	REA571	recombinant human IgG1	130-109-198
iNKT	6B11	mouse IgG1	130-098-130
Integrin β 7	REA441	recombinant human IgG1	130-106-489
Jagged2	MHJ2-523	mouse IgG1	130-097-004
KIR2D	NKVSF1	mouse IgG1	130-099-649
KLRG1 (MAFA)	REA261	recombinant human IgG1	130-103-704
LAP (TGF- β 1)	CH6-17E5.1	mouse IgG1	130-096-573
LGR5	DA03-22H2.8	rat IgG2b	130-100-854
LT- β R	REA412	recombinant human IgG1	130-106-749
Melanoma (MCSP)	EP-1	mouse IgG1	130-091-252
HLA-ABC	REA230	recombinant human IgG1	130-101-466
MICA/MICB	6D4	mouse IgG2ak	130-100-882
MSCA-1 (W8B2)	W8B2	mouse IgG1	130-099-650
NKp80	4A4.D10	mouse IgG1	130-100-608
Notch1	REA357	recombinant human IgG1	130-105-926
Notch2	MHN2-25	mouse IgG2ak	130-096-972
Notch3	MHN3-21	mouse IgG1k	130-096-997
O4	O4	mouse IgM	130-099-211
OSCAR	REA494	recombinant human IgG1	130-107-660
Perforin	delta G9	mouse IgG2bk	130-096-569
Plexin-D1	REA542	recombinant human IgG1	130-108-320
PSA-NCAM	2-2B	mouse IgM	130-093-273
PSMA	REA408	recombinant human IgG1	130-106-609
PTK7 (CCK-4)	188B	mouse IgG2a	130-099-660
ROR1	2A2	mouse IgG1k	130-098-320
Siglec-10	5G6	mouse IgG1k	130-103-731
Siglec-5/Siglec-14	REA393	recombinant human IgG1	130-106-068
Siglec-8	7C9	mouse IgG1k	130-098-731
Slan (M-DC8)	DD-1	mouse IgM	130-099-189
SSEA-1	REA321	recombinant human IgG1	130-104-990
SSEA-4	REA101	recombinant human IgG1	130-098-347
SSEA-5	8e11	mouse IgG1k	130-106-719
SUSD2	W5C5	mouse IgG1k	130-106-400
TCR-V α 7.2	REA179	recombinant human IgG1	130-100-178
TCR V β 11	REA559	recombinant human IgG1	130-108-800
TCR V β 14	REA557	recombinant human IgG1	130-108-805
TCR V β 16	REA553	recombinant human IgG1	130-108-767

Antigen	Clone	Isotype	Order no.
TCR V β 23	REA497	recombinant human IgG1	130-107-497
TCR V γ 9	REA470	recombinant human IgG1	130-107-489
TCR-V δ 1	REA173	recombinant human IgG1	130-100-636
TCR-V δ 2	123R3	mouse IgG1	130-099-664
TCR α / β	BW242/412	mouse IgG2b	130-098-859
TCR γ / δ	11F2	mouse IgG1	130-099-393
TdT	19-3	mouse IgG2bk	130-100-749
TIM-1	REA384	recombinant human IgG1	130-106-074
TIM-3	F38-2E2	mouse IgG1k	130-098-934
TLT-2	REA514	recombinant human IgG1	130-107-840
TRA-1-85 (CD147)	REA476	recombinant human IgG1	130-107-158
TSPAN8	REA443	recombinant human IgG1	130-106-855
KLH	IS5-21F5	mouse IgG1	130-098-846
Hapten NP	S43.10	mouse IgG2a	130-098-850
KLH	IS6-11E5.11	mouse IgG2b	130-098-890
KLH	IS5-20C4	mouse IgM	130-099-085
KLH	ES26-14D1.11	rat IgG1k	130-103-034
KLH	ES26-15B7.3	rat IgG2ak	130-103-092
KLH	ES26-5E12.4	rat IgG2bk	130-103-085
KLH	ES26-13D3.4	rat IgMk	130-103-047



Miltenyi Biotec

Miltenyi Biotec GmbH | Friedrich-Ebert-Straße 68 | 51429 Bergisch Gladbach | Germany | Phone +49 2204 8306-0 | Fax +49 2204 85197
 macs@miltenyibiotec.de | www.miltenyibiotec.com

Miltenyi Biotec provides products and services worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS is a registered trademark of Miltenyi Biotec GmbH. Copyright © 2016 Miltenyi Biotec GmbH. All rights reserved.